

[Manuals.plus](#) /

› [Oral-B](#) /

› Oral-B Cross Action Replacement Brush Heads for an Oral-B Electric Toothbrush, Pack of 4, Black Instruction Manual

Oral-B 10069055127183

Oral-B Cross Action Replacement Brush Heads Instruction Manual

Model: 10069055127183

PRODUCT OVERVIEW

The Oral-B CrossAction Black electric toothbrush replacement head is designed to provide a superior clean by removing up to 100% more plaque compared to a regular manual toothbrush. Its round brush head is inspired by dental professional tools, featuring specially engineered bristles for comprehensive tooth coverage. The Max Clean Bristles change color from green to yellow, indicating when it's time to replace the brush head for optimal cleaning performance.



Image: A pack of four Oral-B Cross Action Black replacement brush heads, showcasing the product packaging and the individual brush heads.

SAFETY INFORMATION

Please read and understand the following safety guidelines before using the Oral-B CrossAction Replacement Brush Heads:

- **Choking Hazard:** Keep out of reach of children under 3 years of age. Small parts may pose a choking risk.
- **Allergic Reactions:** This product contains plastic which may cause an allergic reaction in some individuals. Discontinue use if irritation occurs.
- **Sensitive Areas:** Do not use on sensitive or injured oral areas. Consult a dental professional if you have concerns.
- **Usage Instructions:** Always follow the instructions provided with your specific Oral-B electric toothbrush handle.

SETUP AND INSTALLATION

To attach the Oral-B CrossAction Replacement Brush Head to your compatible Oral-B electric toothbrush handle:

1. Ensure your electric toothbrush handle is turned off.
2. Align the replacement brush head with the metal shaft on the toothbrush handle.
3. Push the brush head down firmly until it clicks into place. Do not twist.
4. To remove, simply pull the brush head straight off the handle.

Compatibility: These replacement brush heads are compatible with Oral-B Pro, Smart, and Genius rechargeable toothbrush handles. They are **not** compatible with Oral-B iO Series toothbrushes.



Image: Illustration of the round brush head design, highlighting its ability to improve gum health and provide a thorough clean.

OPERATING INSTRUCTIONS

Once the brush head is securely attached, you can begin brushing:

1. Wet the bristles and apply a small amount of toothpaste.
2. Place the brush head against your teeth at a slight angle towards the gumline.
3. Turn on your Oral-B electric toothbrush.
4. Guide the brush head slowly from tooth to tooth, allowing the rotating and pulsating action to clean each tooth surface. Do not scrub aggressively.

- Many Oral-B electric toothbrushes feature a built-in timer that pulses every 30 seconds to signal when to move to a different quadrant of your mouth, ensuring a complete 2-minute clean.
- After brushing, turn off the toothbrush and rinse the brush head thoroughly under running water.

MAINTENANCE AND REPLACEMENT

For optimal performance and hygiene, regular maintenance of your brush head is essential:

- Daily Cleaning:** After each use, rinse the brush head thoroughly under running water to remove any toothpaste residue and debris.
- Storage:** Store the brush head upright in a clean, dry place to allow it to air dry.
- Replacement Schedule:** Dentists recommend changing your brush head every **3 months**, or sooner if the bristles become splayed or frayed. The Max Clean Bristles on these brush heads will change from green to yellow to visually signal when a replacement is due.

SPECIFICATIONS

Feature	Detail
Brand	Oral-B
Model Number	10069055127183
UPC	069055127186
Color	Black
Size	4 Count (Pack of 1)
Material	Plastic
Dimensions	0.71 x 2.36 x 8.66 inches
Weight	1.44 ounces
Manufacturer	Procter & Gamble - HABA Hub
Country of Origin	Germany

