

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

Manuals.plus /

› [Terralio](#) /

› Terralio Talking Bathroom Scale User Manual

Terralio unknown_SML

Terralio Talking Bathroom Scale User Manual

Model: unknown_SML

INTRODUCTION

This manual provides essential information for the safe and effective use of your Terralio Talking Bathroom Scale. Designed for ease of use, especially for the elderly and visually impaired, this scale offers accurate weight measurements with audible announcements in multiple languages.



Image: Front view of the Terralio Talking Bathroom Scale, showcasing its white glass platform, blue central accent, and clear digital display.

SETUP

Before first use, please follow these steps to set up your scale:

1. **Unpacking:** Carefully remove the scale from its packaging. Inspect for any damage.
2. **Battery Installation:** The scale requires 3 AAA batteries (not included). Locate the battery compartment on the underside of the scale. Insert the batteries, ensuring correct polarity (+/-).
3. **Placement:** Place the scale on a hard, flat, and stable surface. Avoid carpets or uneven flooring, as this can affect accuracy.
4. **Initial Calibration:** After installing batteries and placing the scale, step on it briefly to turn it on, then step off. The scale will display "0.0" and then turn off. This calibrates the scale for its environment.

OPERATING INSTRUCTIONS

Using your Terralion Talking Bathroom Scale is straightforward:

1. **Power On:** Gently tap the scale with your foot or step directly onto the platform. The scale will automatically turn on.
2. **Weighing:** Stand still on the center of the scale platform with bare feet for best accuracy. The display will show your weight, and the scale will announce it audibly.
3. **Language Selection:** To change the spoken language (English, French, German, Spanish), locate the language selection button on the underside of the scale. Press it repeatedly to cycle through the available languages until you hear your desired language.
4. **Volume Adjustment:** The volume can be adjusted (low, high, mute) using the volume control button, also typically located on the underside of the scale.
5. **Unit Conversion:** To switch between kilograms (kg), pounds (lb), and stones (st), use the unit conversion button, usually found near the battery compartment.
6. **Automatic Shut-off:** The scale will automatically turn off after a few seconds of inactivity to conserve battery life.



Image: A close-up view of the scale's LCD screen, clearly showing a weight reading of "56.3".

MAINTENANCE

Proper maintenance ensures the longevity and accuracy of your scale:

- **Cleaning:** Wipe the scale surface with a soft, damp cloth. Do not use abrasive cleaners or immerse the scale in water.
- **Storage:** Store the scale in a dry place at room temperature. Avoid extreme temperatures or high humidity.
- **Battery Replacement:** Replace batteries when the display indicates low power or when the scale fails to turn on. Always replace all batteries at the same time with new ones of the same type.
- **Avoid Overload:** Do not exceed the maximum weight capacity of 160 kg (25 stone).

TROUBLESHOOTING

If you encounter issues with your scale, refer to the following common problems and solutions:

Problem	Possible Cause	Solution
Scale does not turn on.	Batteries are dead or incorrectly installed.	Check battery polarity or replace with new AAA batteries.
Inaccurate readings.	Scale on uneven surface; not calibrated; overloaded.	Place on a hard, flat surface. Step on and off to re-calibrate. Ensure weight is within 160 kg limit.
No audible announcement.	Volume set to mute; language not selected.	Adjust volume using the volume button. Select a language using the language button.
Display shows "Lo".	Low battery power.	Replace all 3 AAA batteries.
Display shows "Err".	Overload or internal error.	Step off the scale. Ensure weight is below 160 kg. If error persists, remove batteries for 1 minute and reinsert.

SPECIFICATIONS

- **Brand:** Terralion
- **Model:** unknown_SML
- **Maximum Weight Capacity:** 160 kg (25 stone)
- **Measurement Increments:** 100 g
- **Measurement Units:** Kilograms (kg), Pounds (lb), Stones (st)
- **Display Type:** LCD with large digits (3.5 cm high)
- **Special Features:** Talking function (4 languages: English, French, German, Spanish), Adjustable volume (low, high, mute), Automatic on/off, Self-adjusting feet.
- **Power Source:** 3 x AAA batteries (not included)
- **Materials:** Glass platform, ABS plastic shell
- **Dimensions:** Approximately 32.5 x 32.5 x 2.5 cm
- **Item Weight:** Approximately 1.5 kg

WARRANTY INFORMATION

Specific warranty details are not provided in the product description. Please refer to the product packaging or contact Terralion customer support for warranty terms and conditions.

CUSTOMER SUPPORT

For further assistance, technical support, or inquiries regarding your Terralion Talking Bathroom Scale, please contact Terralion customer service. Contact information is typically available on the product packaging or the official Terralion website.

Note: As per the provided information, direct contact details are not available. Please refer to the manufacturer's official channels.



© 2024 Terralio. All rights reserved.