

[Manuals.plus](#) /

> [Oral-B](#) /

> Oral-B Genius 9000 Rose Gold Electric Toothbrush User Manual

Oral-B 2724542399054

Oral-B Genius 9000 Rose Gold Electric Toothbrush User Manual

Brand: **Oral-B**

1. PRODUCT OVERVIEW

The Oral-B Genius 9000 in rose gold is an advanced electric toothbrush designed to provide a superior cleaning experience. It incorporates innovative features to help you brush effectively and maintain optimal oral hygiene.

Key features include:

- **Position Detection Technology:** Helps ensure no zone is missed during brushing.
- **Unique Round Head:** Designed to remove up to 100% more plaque compared to a standard manual toothbrush.
- **SmartRing and Pressure Control:** Protects delicate gums by reducing brushing speed and alerting you if you brush too hard.
- **Whitening Capability:** Removes surface stains for a whiter smile from the first day of use.
- **Long-lasting Lithium-Ion Battery:** Provides up to 12 days of brushing on a single charge.
- **Bluetooth Connectivity:** Enables connection to the Oral-B app for real-time feedback and guidance.



Image 1.1: The Oral-B Genius 9000 Rose Gold electric toothbrush, charging stand, brush heads, and smartphone holder.

This package includes the Oral-B Genius 9000 electric rechargeable toothbrush in Rose Gold and four replacement brush heads.

2. SETUP AND FIRST USE

2.1 Charging the Toothbrush

Before first use, fully charge the toothbrush. Place the handle on the charging unit. The battery indicator will show the charging progress. A full charge can take up to 12 hours and provides approximately 12 days of regular brushing.

2.2 Attaching a Brush Head

Firmly push the desired brush head onto the toothbrush handle. Ensure it clicks into place. The Oral-B Genius 9000 is compatible with various Oral-B replacement heads, including CrossAction, 3DWhite, Sensi UltraThin, Sensitive Clean,

Precision Clean, Floss Action, and TriZone.

2.3 Connecting to the Oral-B App (Bluetooth)

1. Download the Oral-B app from your smartphone's app store.
2. Ensure Bluetooth is enabled on your smartphone.
3. Open the Oral-B app and follow the on-screen instructions to pair your Genius 9000 toothbrush.
4. The app provides real-time feedback, brushing guidance, and tracks your brushing habits, utilizing the Position Detection technology.



**NEVER MISS A ZONE WITH
POSITION DETECTION**

Image 2.1: The Position Detection feature helps users cover all areas of their mouth during brushing, as indicated by the on-screen guide.

3. OPERATING THE TOOTHBRUSH

3.1 Brushing Technique

Moisten the brush head and apply a small amount of toothpaste. Place the brush head on your teeth, holding it at a slight

angle towards the gumline. Turn on the toothbrush. Guide the brush head slowly from tooth to tooth, spending a few seconds on each surface. The unique round head oscillates, rotates, and pulsates to break up and remove plaque.



Image 3.1: The Oral-B Genius 9000 toothbrush ready for use in a typical bathroom environment.

3.2 SmartRing and Pressure Control

The SmartRing illuminates in different colors to provide visual feedback. If you apply too much pressure during brushing, the SmartRing will turn red, and the brushing speed will automatically reduce to protect your gums from excessive force. This feature helps prevent over-brushing and promotes healthier gums.



Image 3.2: Illustration of the triple pressure control mechanism on the brush head, designed to protect gums.

3.3 Brushing Modes

The Oral-B Genius 9000 offers multiple brushing modes. The primary mode is "Daily Clean." Refer to the Oral-B app for detailed information on available modes and how to switch between them for a personalized brushing experience, including modes for whitening.



Image 3.3: The Oral-B app interface displaying a whitening progress chart, demonstrating how the toothbrush gently whitens teeth by removing surface stains.

4. MAINTENANCE AND CARE

4.1 Cleaning the Toothbrush

After each use, rinse the brush head thoroughly under running water with the handle switched on. Remove the brush head and clean the handle and brush head separately. Wipe the handle with a damp cloth. Do not immerse the charging unit in water.

4.2 Replacing Brush Heads

Dentists recommend replacing your toothbrush head every 3 months, or sooner if the bristles become splayed or frayed. The indicator bristles on Oral-B brush heads will fade halfway to signal when it's time for a replacement.



REPLACE BRUSH HEADS EVERY 3 MONTH

Image 4.1: Visual guide indicating the need to replace brush heads every three months for optimal performance.

5. TROUBLESHOOTING

If you encounter any issues with your Oral-B Genius 9000 toothbrush, please refer to the following common solutions:

- **Toothbrush not charging:** Ensure the charging unit is properly plugged into a working power outlet and the toothbrush handle is correctly seated on the charger. Check for any debris preventing proper contact.
- **Reduced brushing power:** A low battery can cause reduced power. Fully charge the toothbrush. If the issue persists after charging, the brush head might be worn out and needs replacement.
- **Bluetooth connectivity issues:** Ensure Bluetooth is enabled on your smartphone and the Oral-B app is updated to the latest version. Try restarting both your phone and the toothbrush.
- **SmartRing not responding:** Ensure the pressure sensor is not obstructed. If the issue persists, try resetting the toothbrush by fully discharging and then recharging it.

For more complex issues or persistent problems, please contact Oral-B customer support.

6. SPECIFICATIONS

| Feature | Detail |
|-----------------------|---|
| Product Name | Oral B Genius 9000 Electric Rechargeable Toothbrush - RoseGold |
| Brand | Oral-B |
| Model Number | 2724542399054 |
| ASIN | B07NFL82Z3 |
| Color | Rose Gold |
| Power Source | Battery Powered (Lithium-Ion) |
| Special Features | Brushing Timer, Position Detection, SmartRing, Pressure Control |
| Age Range Description | Adult |
| Product Dimensions | 10.4 x 25.8 x 23.8 cm |
| Item Weight | 815 kg <i>(Note: This weight seems unusually high for a toothbrush, likely a data error in source. Assuming it refers to package weight or a typo for grams.)</i> |
| Manufacturer | Procter & Gamble |
| Toothbrush Head Shape | Round |
| Number of Items | 1 (toothbrush unit) |
| Units | 1 Count |

7. WARRANTY AND SUPPORT

For information regarding product warranty, please refer to the warranty card included with your purchase or visit the official Oral-B website. Oral-B products typically come with a manufacturer's warranty covering defects in materials and workmanship.

If you require technical assistance, have questions about your product, or need to claim warranty service, please contact Oral-B customer support. Contact details can usually be found on the Oral-B website or on the product packaging.

You can also visit the official Oral-B Store for more information:[Oral-B Store on Amazon](#).