

Vivitar MS-6001

Vivitar MS-6001 Hand Held Deep Tissue Massager Instruction Manual

Model: MS-6001

1. INTRODUCTION

Thank you for purchasing the Vivitar MS-6001 Hand Held Deep Tissue Massager. This device is designed to provide targeted relief for muscle aches and pains, promote blood circulation, and aid in muscle recovery. Its long handle and cordless design offer flexibility and ease of use for various body parts, including the back, shoulders, and neck.

Please read this instruction manual thoroughly before using the massager to ensure safe and effective operation. Keep this manual for future reference.



Figure 1: Vivitar MS-6001 Hand Held Deep Tissue Massager (Teal variant shown in packaging). This image displays the product in its retail box, highlighting its long handle design and key features like stress relief, deep tissue massage, and

2. SAFETY INFORMATION

WARNING: To reduce the risk of electric shock, burns, fire, or injury, always follow these important safety precautions.

- This is a battery-powered device. Always remove batteries when the massager is not in use for extended periods.
- Do not use the massager while sleeping or in the shower.
- Keep the device away from heat sources and open flames.
- Use the massager only on the intended body parts (back, muscles, shoulders, neck). Avoid sensitive areas of the body.
- Do not apply excessive pressure during use.
- Consult a healthcare professional before use if you have any medical conditions, are pregnant, or have any concerns about using a massager.
- Keep out of reach of children.

3. SETUP AND BATTERY INSTALLATION

The Vivitar MS-6001 massager requires 4 AA batteries (not included) for operation.

1. Locate the battery compartment on the massager handle.
2. Open the battery compartment cover.
3. Insert 4 new AA batteries, ensuring the correct polarity (+ and -) as indicated inside the compartment.
4. Close the battery compartment cover securely.

Note: Always use fresh batteries for optimal performance. Do not mix old and new batteries, or different types of batteries.

4. OPERATING INSTRUCTIONS

Follow these steps to operate your Vivitar MS-6001 massager:

1. Ensure batteries are correctly installed (refer to Section 3).
2. Hold the massager firmly by its long handle.
3. Press the power button, typically located on the handle, to turn the massager ON.
4. Gently apply the vibrating head of the massager to the desired muscle area. Move the massager slowly over the area, applying light to moderate pressure.
5. Use the long handle to reach difficult areas such as your back and shoulders.
6. To turn the massager OFF, press the power button again.

Recommendation: Start with short massage sessions (5-10 minutes) and gradually increase duration as desired. Avoid using on one spot for too long.

5. MAINTENANCE AND STORAGE

Proper care will extend the life of your massager.

- **Cleaning:** Disconnect the massager from power (remove batteries) before cleaning. Wipe the surface of the massager with a soft, damp cloth. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the massager in a cool, dry place, away from direct sunlight and extreme temperatures. If storing for an extended period, remove the batteries to prevent leakage.
- **Battery Replacement:** Replace all 4 AA batteries when the massager's performance diminishes. Dispose of used batteries responsibly according to local regulations.

6. TROUBLESHOOTING

| Problem | Possible Cause | Solution |
|---------------------------------|----------------------------------------------|-----------------------------------------------------------------------------------------------|
| Massager does not turn on. | Batteries are dead or incorrectly installed. | Check battery polarity. Replace with new AA batteries. |
| Weak vibration. | Low battery power. | Replace all 4 AA batteries with new ones. |
| Unusual noise during operation. | Internal component issue. | Discontinue use and contact customer support if the issue persists after battery replacement. |

If you experience any issues not listed above, please refer to the warranty and support section.

7. SPECIFICATIONS

Model Number: MS-6001

Power Source: 4 x AA Batteries (not included)

Dimensions: 4.7 x 2.2 x 14 inches (11.9 x 5.6 x 35.6 cm)

Weight: 9.91 ounces (281 grams)

Brand: Vivitar

Manufacturer: Sakar International HPC

Special Feature: Travel Size, Cordless, Long Handle

Intended Use: Back, Neck, Shoulders, Muscles

8. WARRANTY AND SUPPORT

For information regarding warranty coverage, technical support, or replacement parts, please refer to the documentation included with your product packaging or visit the official Vivitar website. Please have your model number (MS-6001) and proof of purchase ready when contacting support.

Vivitar Official Website: www.vivitar.com