

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

[manuals.plus](#) /

> [Live Fine](#) /

> [LiveFine Replacement Key User Manual](#)

Live Fine LVFBPDKEY

LiveFine Replacement Key User Manual

Model: LVFBPDKEY

INTRODUCTION

This manual provides essential information for the proper use and maintenance of your LiveFine Replacement Key, model LVFBPDKEY. This key is designed to secure and access specific Live Fine automatic pill dispensers. Please read this manual thoroughly to ensure correct usage and compatibility.

SAFETY INFORMATION

Warning: Choking hazard - small parts. Keep out of reach of children.

PRODUCT OVERVIEW

The LiveFine Replacement Key is a durable accessory crafted from plastic and stainless steel, designed to provide secure access to compatible Live Fine automatic pill dispensers. Its robust construction ensures longevity and reliable performance.



Image: A close-up view of the LiveFine Replacement Key, showcasing its silver finish and design.

- **Seamless Fit:** Ensures secure locking and unlocking of your pill dispenser.
- **Durable Construction:** Made from high-quality plastic and stainless steel for long-lasting use.
- **Essential Accessory:** Provides peace of mind by securing medication within the dispenser.

COMPATIBILITY

It is crucial to verify the compatibility of this replacement key with your specific Live Fine pill dispenser

model before use. This key (Model LVFBPDKEY) is designed for the following Live Fine pill dispenser models:

- LVFPDBTC
- LVFRAPD5
- LVFRAPD5W

Note: This key is NOT compatible with models LVFZWPDC or LVFZWPDF, or other models not listed above. Please check your dispenser's model number to ensure proper fit.



Image: A visual guide illustrating the compatibility of the LVFBPDKEY with specific LiveFine pill dispenser models (LVFPDBTC, LVFRAPD5, LVFRAPD5W) and showing models it is not compatible with.

SETUP AND OPERATION

The replacement key is used to securely lock and unlock the lid of your compatible Live Fine pill dispenser, preventing unauthorized access and ensuring medication integrity.

1. **Verify Compatibility:** Before attempting to use the key, confirm that your pill dispenser model is listed as compatible with the LVFBPDKEY.
2. **Insert Key:** Locate the keyhole on your Live Fine pill dispenser. Gently insert the replacement key into the keyhole.



Image: A hand demonstrating the insertion of the replacement key into the keyhole of a LiveFine pill dispenser.

3. **Locking:** To lock the dispenser, insert the key and turn it clockwise until it clicks into place, securing the lid.
4. **Unlocking:** To unlock the dispenser, insert the key and turn it counter-clockwise. Once unlocked, gently pull the key out and open the dispenser lid.

WIRELESS OPERATION WITH INCLUDED BATTERIES



Included Key Keeps Your Medication Secure

Image: A detailed view of the key insertion point on the LiveFine pill dispenser, showing the mechanism for securing the device.

MAINTENANCE

To ensure the longevity and proper function of your LiveFine Replacement Key:

- Keep the key clean and free from dust or debris. Wipe with a soft, dry cloth if necessary.
- Avoid exposing the key to extreme temperatures or harsh chemicals.
- Always store the key in a safe place to prevent loss or damage. Consider attaching it to a keyring or keeping it in a designated drawer.

TROUBLESHOOTING

If you encounter issues with your LiveFine Replacement Key, consider the following:

- **Key Does Not Fit:** The most common issue is incompatibility. Double-check your pill dispenser's model number against the compatible models listed in the "Compatibility" section of this manual. This key (LVFBPDKEY) is specifically designed for models LVFPDBTC, LVFRAPD5, and LVFRAPD5W. If your model is different, this key will not fit.
- **Key Is Stuck or Difficult to Turn:** Ensure the key is fully inserted into the keyhole. Do not force the key. Check for any obstructions in the keyhole. If the dispenser mechanism is damaged, the key may not operate smoothly.
- **Lost Key:** If your original key is lost, this replacement key serves its purpose. However, always keep

spare keys in a secure, memorable location.

SPECIFICATIONS

Model Number	LVFBPDKEY
Material	Plastic, Stainless Steel
Color	Silver
Dimensions	6.35 x 3.81 x 0.63 cm (2.5 x 1.5 x 0.25 inches)
Weight	5 g (0.18 oz)
Manufacturer	Live Fine

LEGAL DISCLAIMER

Actual product packaging and materials may contain more and different information than what is shown on our website. We recommend that you do not rely solely on the information presented and that you always read labels, warnings, and directions before using or consuming a product.

WARRANTY AND SUPPORT

For information regarding product warranty or customer support, please refer to the documentation included with your original Live Fine pill dispenser or visit the official Live Fine website. Live Fine is dedicated to providing reliable solutions and excellent customer service.