

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

› [MEDca](#) /

› [Pill Crusher - 304 Food Grade Stainless Steel Mortar and Pestle Medicine Grinder Set - Non-Slip Splitter to Easily Crush Medicine Pills Tablets Vitamins to Fine Powder for Adults, Seniors, Dogs, Pets](#)

MEDca MED-215

MEDca Stainless Steel Pill Crusher User Manual

Model: MED-215

INTRODUCTION

The MEDca 304 Food Grade Stainless Steel Mortar and Pestle is designed for efficiently crushing medicine pills, tablets, vitamins, and even kitchen herbs and spices into a fine powder. Its robust construction ensures durability and ease of use for various grinding needs.



Image: The MEDca Stainless Steel Mortar and Pestle set, featuring a bowl-shaped mortar and a cylindrical pestle, both with a brushed metal finish.

PRODUCT FEATURES

- **Material:** Crafted from Type 304 Food Grade Stainless Steel for durability and hygiene.
- **Design:** Features a non-slip base for stability during use and curved surfaces for efficient crushing.
- **Versatility:** Suitable for crushing pills, tablets, vitamins, and various kitchen ingredients like herbs and spices.
- **Ease of Cleaning:** Non-porous stainless steel surface allows for easy cleaning with soap and water.
- **Efficiency:** Designed for maximum crushing efficiency, turning ingredients into fine powder.



**Superior
Crafted
Quality**

- ✔ **Non-Slip**
- ✔ **Easy to Clean**
- ✔ **Easy to Use**
- ✔ **High Quality**



Image: The MEDca Stainless Steel Pill Crusher alongside its product packaging, highlighting features such as Non-Slip, Easy to Clean, Easy to Use, and High Quality.

304 Food Grade Stainless Steel Mortar and Pestle

designed to give
you **maximum
crushing efficiency**
to pulverize
medicine, pills,
kitchen herbs and
spices to the
tiniest bits



Image: The MEDca Mortar and Pestle with text emphasizing its 304 Food Grade Stainless Steel construction and design for maximum crushing efficiency.

- ✔ Corrosion Resistant
- ✔ Keeps Mortar in Place
- ✔ Sturdy & Durable

Slip Resistant

heavy base with an adhesive foam bottom that is slip resistant and **would not scratch your countertops**



Image: An inverted view of the mortar showing its heavy base with an adhesive foam bottom, designed to be slip-resistant and prevent scratching countertops.

SETUP

The MEDca Stainless Steel Mortar and Pestle requires minimal setup. Upon unboxing, follow these steps:

1. **Unpack:** Carefully remove the mortar and pestle from its packaging.
2. **Inspect:** Check both components for any signs of damage during transit.
3. **Clean:** Before first use, wash the mortar and pestle thoroughly with warm water and mild soap. Rinse completely and dry with a clean cloth. Ensure it is fully dry before storage or use.

Easy to Clean

*Just use soap and warm water to wash
and dry with a washcloth*



Image: A hand holding the mortar under running water from a faucet, demonstrating the ease of cleaning the stainless steel surface.

OPERATING INSTRUCTIONS

Using your MEDca Mortar and Pestle for crushing pills or grinding ingredients is straightforward:

1. **Placement:** Place the mortar on a stable, flat, and non-slip surface. The integrated non-slip base will help keep it secure.
2. **Add Ingredients:** Place the pills, tablets, or other ingredients into the mortar. For best results, add a few items at a time rather than overfilling.
3. **Crushing/Grinding:** Hold the pestle firmly. Use a downward pressure combined with a circular or pounding motion to crush and grind the contents against the inner surface of the mortar. The curved design ensures efficient contact.
4. **Consistency:** Continue grinding until the desired consistency (fine powder) is achieved.
5. **Removal:** Once finished, carefully tap the pestle against the rim of the mortar to dislodge any remaining powder. Use a small spoon or brush to remove the crushed material.

Easy to Use

304 food grade stainless steel to give you **the perfect weight to crush pills at ease** while avoiding rusting and deterioration

Ideal for:

- Medicine
- Pills
- Herbs & Spices



Image: A hand holding the pestle over the mortar, illustrating the ease of use for crushing pills or other ingredients.



Maximum Crushing Efficiency

304 food grade stainless steel to give you the perfect weight to crush pills at ease while avoiding rusting and deterioration

Image: A close-up view of the mortar with several white pills and some already crushed powder inside, demonstrating the crushing process with the pestle.



Image: A collage showing various applications of the mortar and pestle, including grinding spices and nuts, preparing vegetable pastes, crushing herbs, and pulverizing medicines.

MAINTENANCE AND CARE

Proper care will ensure the longevity and performance of your MEDca Mortar and Pestle:

- **Cleaning:** After each use, wash both the mortar and pestle with warm water and mild dish soap. Use a soft sponge or cloth.
- **Rinsing:** Rinse thoroughly to remove all soap residue.
- **Drying:** Dry immediately with a clean, soft cloth to prevent water spots and maintain the stainless steel finish. Ensure both parts are completely dry before storage.
- **Storage:** Store in a dry place.
- **Avoid Harsh Cleaners:** Do not use abrasive cleaners, steel wool, or harsh chemicals, as these can scratch or damage the stainless steel surface.

TROUBLESHOOTING

| Problem | Possible Cause | Solution |
|---|--|---|
| Pills not crushing effectively. | Too many pills at once; insufficient pressure; pills are too hard. | Reduce the number of pills; apply more consistent downward pressure; for very hard pills, consider breaking them into smaller pieces first if safe to do so. |
| Mortar slides during use. | Surface is slippery; non-slip base is dirty. | Ensure the surface is clean and dry; clean the non-slip base of the mortar. |
| Residue or odor remains after cleaning. | Incomplete cleaning; strong-smelling ingredients. | Wash thoroughly with warm water and soap, ensuring all surfaces are scrubbed. For persistent odors, a paste of baking soda and water can be used, then rinsed well. |

SPECIFICATIONS

- **Model:** MED-215
- **Material:** 304 Food Grade Stainless Steel
- **Product Dimensions (L x W x H):** 9.4 x 11.9 x 14 Centimetres (approx. 3.7 x 4.7 x 5.5 inches)
- **Item Weight:** 739 g (approx. 1.63 lbs)
- **UPC:** 738246426215
- **Color:** As Shown (Stainless Steel)
- **Finish Type:** Polished

Crafted with curved surfaces to effectively crush ingredients to its smallest bits



Image: A diagram showing the key dimensions of the MEDca Mortar and Pestle, including the height and diameter of the mortar and the length of the pestle.

WARRANTY AND SUPPORT

MEDca stands behind the quality of its products. Your satisfaction is our top priority. For any questions, concerns, or support needs regarding your MEDca Stainless Steel Pill Crusher, please contact our customer service. We are committed to ensuring you have an excellent customer experience.

For further assistance, please refer to the contact information provided with your purchase or visit the official MEDca website.

IMPORTANT INFORMATION

Actual product packaging and materials may contain more and different information than what is shown on our website. We recommend that you do not rely solely on the information presented and that you always read labels, warnings, and directions before using or consuming a product.

