

## Frida Baby Quick-Read Digital Rectal Thermometer

# Frida Baby Quick-Read Rectal Thermometer Instruction Manual

Model: Quick-Read Digital Rectal Thermometer (088)

## 1. INTRODUCTION

The Frida Baby Quick-Read Rectal Thermometer is designed to provide fast and accurate temperature readings for infants and newborns. Its flexible tip and parent-proof stopper ensure safe and comfortable use, while the digital display and memory function offer convenience for monitoring your baby's health.

### Key Features:

- **Precision for Peace of Mind:** Delivers accurate temperature readings in approximately 10 seconds.
- **Memory Recall Function:** Stores up to 10 previous temperature readings for tracking.
- **Gentle, Flexible Tip:** Ensures maximum comfort during use.
- **Water-Resistant Body:** Designed for easy cleaning.
- **Portable Protective Case:** Keeps the thermometer clean and safe.



## All in the Deets

- Memory holds 10 readings
- \*F or \*C
- Keep-clean case
- Water-resistant
- BPA-free

Image: Overview of the Frida Baby Quick-Read Rectal Thermometer's key features, including memory, temperature units, protective case, water resistance, and BPA-free materials.

## 2. IMPORTANT SAFETY INFORMATION

---

Please read all instructions carefully before using the thermometer. Failure to follow these instructions may result in inaccurate readings or injury.

- This thermometer is intended for rectal use only on infants and newborns. Do not use orally, axillary, or for any other purpose.
- Always supervise your child during temperature measurement.
- Do not allow children to play with the thermometer.
- The flexible tip is designed for comfort, but always insert gently. The built-in stopper helps prevent over-insertion.
- Consult a healthcare professional immediately if you have concerns about your child's temperature or health.
- Keep batteries out of reach of children. Swallowing batteries can be fatal.
- Do not attempt to disassemble or repair the thermometer.

## Parent-proof stopper



Image: The Frida Baby Rectal Thermometer highlighting its parent-proof stopper, designed to prevent over-insertion during use.

### 3. PACKAGE CONTENTS

---

Upon opening the package, ensure all items are present and undamaged:

- 1 x Frida Baby Quick-Read Rectal Thermometer
- 1 x Protective Carry Case



Image: The Frida Baby Quick-Read Rectal Thermometer shown alongside its retail packaging, illustrating the product and its included protective case.

## 4. PRODUCT OVERVIEW

---

Familiarize yourself with the components of your Frida Baby Quick-Read Rectal Thermometer:

- **Digital Display:** Shows temperature readings and indicators.
- **Power Button:** Turns the thermometer on and off.
- **Flexible Tip:** The soft, gentle probe for comfortable insertion.
- **Parent-Proof Stopper:** A safety feature to prevent over-insertion.
- **Battery Compartment:** Located at the rear of the device.



Image: A close-up of the Frida Baby Rectal Thermometer held in a hand, displaying a temperature reading of 98.2°F on its digital screen.

## Flexible tip for gentle insertion



Image: A hand gently bending the flexible tip of the Frida Baby Rectal Thermometer, illustrating its design for gentle insertion.

## 5. SETUP AND BATTERY INSTALLATION

---

### Initial Cleaning:

Before first use, and after every subsequent use, clean the thermometer as described in the 'Maintenance and Care' section.

### Battery Installation/Replacement:

The thermometer comes with one unknown battery type included and installed. If the display shows a low battery indicator or does not turn on, replace the battery:

1. Carefully open the battery compartment cover, usually located at the back or side of the thermometer.
2. Remove the old battery.
3. Insert a new battery of the same type, ensuring correct polarity (+/-).
4. Securely close the battery compartment cover.

**Note:** Dispose of used batteries according to local regulations.



## 6. OPERATING INSTRUCTIONS

---

### Taking a Rectal Temperature:

1. Ensure the thermometer is clean. You may apply a small amount of petroleum jelly to the tip for easier insertion.
2. Press the power button to turn on the thermometer. The display will show the last reading briefly, then 'Lo' or '---' indicating it's ready.
3. Position your baby on their back with legs raised, or on their stomach.
4. Gently insert the flexible tip into the rectum, approximately 1/2 to 1 inch (1.5 inches maximum immersion depth), until the parent-proof stopper prevents further insertion. Do not force the thermometer.
5. Hold the thermometer in place until it beeps, indicating the reading is complete (approximately 10 seconds).
6. Carefully remove the thermometer and read the temperature displayed on the screen.
7. Press the power button to turn off the thermometer, or it will automatically shut off after a period of inactivity to conserve battery life.

### Memory Function:

The thermometer automatically stores the last 10 temperature readings. To recall previous readings:

- With the thermometer off, press and hold the power button for a few seconds until 'MEM' appears on the display, followed by the last reading.
- Each subsequent short press of the power button will cycle through the previous readings.

### Fahrenheit/Celsius Switching:

To switch between Fahrenheit (°F) and Celsius (°C):

- With the thermometer off, press and hold the power button for approximately 5 seconds until the unit changes. Release the button.
- The thermometer will now display temperatures in the selected unit.

## Backlit screen for nighttime check-ins



Image: The Frida Baby Rectal Thermometer's backlit screen displaying a temperature of 98.2°F, shown in a dimly lit environment for nighttime visibility.

## 7. MAINTENANCE AND CARE

---

### Cleaning:

The thermometer is water-resistant and designed for easy cleaning. Always clean the thermometer before and after each use.

- Wipe the probe and the entire thermometer body with a soft cloth dampened with rubbing alcohol (70% isopropyl alcohol) or soap and water.
- Rinse the tip under cool running water.
- Do not immerse the entire thermometer in water for extended periods.
- Allow the thermometer to air dry completely before storing.
- Do not use abrasive cleaners, thinners, or benzene for cleaning.

### Storage:

- Store the thermometer in its protective carry case when not in use.
- Keep in a cool, dry place, away from direct sunlight and extreme temperatures.



- Avoid dropping the thermometer or subjecting it to strong shocks.

## 8. TROUBLESHOOTING

Problem	Possible Cause	Solution
Thermometer does not turn on.	Dead battery; Battery incorrectly installed.	Replace battery; Reinstall battery with correct polarity.
Display shows 'Err'.	Measurement error; Temperature outside operating range.	Ensure proper contact and stable environment; Wait for thermometer to return to room temperature.
Inaccurate readings.	Improper insertion; Movement during measurement; Dirty probe.	Ensure correct insertion depth and stability; Clean the probe thoroughly.
Display is dim.	Low battery.	Replace battery.

## 9. SPECIFICATIONS

Feature	Detail
Brand	Frida Baby
Model Name	Quick-Read Digital Rectal Thermometer
Item Model Number	088
Special Feature	Fast Reading System
Display Type	Digital
Response Time	10 seconds
Upper Temperature Rating	104 Degrees Fahrenheit
Age Range (Description)	Infant
Immersion Depth	1.5 Inches
Power Source	Battery Powered (1 Unknown battery included)
Product Care Instructions	Hand Wash Only
Water Resistance	Yes
Memory Function	Stores last 10 readings
Material	Rubber (Outer Material)
Item Weight	2.39 ounces
Item Length	6.75 Inches

## 10. WARRANTY AND SUPPORT

---

For warranty information, product support, or to address any concerns not covered in this manual, please refer to the warranty card included with your product or contact Frida Baby customer service directly. Contact details can typically be found on the manufacturer's official website or product packaging.