

manuals.plus /

- › [zinnor](#) /
- › [zinnor Shockwave Therapy Machine Muscle Massager User Manual](#)

**zinnor B082XY1GR6**

# zinnor Shockwave Therapy Machine User Manual

Model: B082XY1GR6 | Brand: zinnor

## 1. INTRODUCTION

---

This manual provides essential instructions for the safe and effective use of your zinnor Shockwave Therapy Machine. This device is designed primarily for pain relief and to promote blood circulation. It is a non-invasive, convenient, and fast treatment method, offering various professional probes for different applications. The treatment is non-invasive, does not interfere with normal daily activities, and aims for faster and easier healing. It is widely applied to various soft tissue treatments.



## Effectively relieve pain

*Clinically Proven Pain Relief for Minor to Chronic Pain, Muscle Recovery, Joint Pain, Inflammation, Stiffness. Pain relieve and treatment for plantar fasciitis, Achilles tendinitis, shoulder joint tendon inflammation and calcification, groin pain, epicondylitis. The massage effect can effectively reach the fascial layer and effectively reduce the problem of muscle tension.*

Figure 1: zinnor Shockwave Therapy Machine and target areas for pain relief.

## 2. PRODUCT FEATURES

- **Multi-function:** Primarily used for pain relief and blood circulation promotion.
- **Convenient and Fast:** Non-invasive treatment with various professional probes.
- **Safe and Painless:** Non-invasive treatment that does not affect normal life, promoting faster and easier healing.
- **High Effectiveness:** Designed for effective application in various soft tissue treatments.

## Suitable crowd

**The elderly:** low back pain caused by physical decline.

**Sports and fitness:** Inflammation of the fascia such as shoulder joint tendon inflammation and calcification, groin pain, epicondylitis.

**Office clerk:** Back pain caused by sedentary movements, Wrist pain caused by prolonged typing.



Figure 2: The zinnor Shockwave Therapy Machine is suitable for various users seeking pain relief.

### 3. PACKAGE CONTENTS

Upon opening the product case, verify that all components are present and undamaged:

- Main Control Unit
- Control Handle (Applicator)
- Power Cable
- Seven (7) Massage Heads (various shapes and sizes)



Figure 3: The seven interchangeable massage heads.

## 4. SETUP INSTRUCTIONS

Follow these steps to set up your zinnor Shockwave Therapy Machine:

1. **Unpack the Device:** Carefully remove the main control unit, control handle, power cable, and all seven massage heads from the carrying case.
2. **Connect the Control Handle:** Locate the output port on the main control unit (labeled 'OUTPUT' in the video) and securely connect the cable from the control handle to this port.
3. **Connect the Power Cable:** Plug the power cable into the power input port on the back of the main control unit. Then, plug the other end of the power cable into a suitable electrical outlet.
4. **Power On:** Flip the power switch located near the power input port on the back of the main control unit to the 'ON' position. A red indicator light should illuminate.
5. **Attach a Massage Head:** Select the appropriate massage head for your intended treatment. Gently screw the chosen massage head onto the end of the control handle. Ensure it is securely fastened but do not overtighten, as the base is plastic and can be damaged.

Your browser does not support the video tag.

Video 1: Demonstrates unpacking, connecting the control handle and power, and attaching massage heads.

## 5. OPERATING INSTRUCTIONS

---

Once the device is set up, follow these steps for operation:

1. **Select Treatment Program:** On the main control unit's interface, select the desired treatment area or program. The interface features various visual indicators for different applications.
2. **Activate the Device:** With the selected massage head attached, press the green power button located on the control handle to activate the shockwave function. You will hear and feel the device operating.
3. **Apply to Treatment Area:** Gently apply the massage head to the target area on the body. Move the device slowly and steadily over the area.
4. **Adjust Settings (if applicable):** Refer to the main control unit for any adjustable settings such as intensity or frequency, and modify as needed for comfort and effectiveness.
5. **Stop Treatment:** To pause or stop the shockwave, press the green power button on the control handle again.
6. **Change Massage Heads:** To change a massage head, first ensure the device is not active. Unscrew the current head and carefully screw on the new one.



Figure 4: The green power button on the control handle for activation.

## 6. MAINTENANCE

---

Proper maintenance ensures the longevity and optimal performance of your device:

- **Cleaning:** After each use, wipe the control handle and massage heads with a clean, damp cloth. Do not use

harsh chemicals or abrasive cleaners. Ensure the device is unplugged before cleaning.

- **Storage:** Store the device and all accessories in its original carrying case in a cool, dry place, away from direct sunlight and extreme temperatures.
- **Massage Head Care:** When attaching or detaching massage heads, handle them carefully. The base of the massage heads is plastic and can be damaged if mis-threaded or overtightened.

## 7. TROUBLESHOOTING

---

If you encounter any issues with your zinnor Shockwave Therapy Machine, please refer to the following common solutions:

- **Device Not Powering On:** Ensure the power cable is securely connected to both the device and the electrical outlet, and that the power switch on the back of the unit is in the 'ON' position. Check the power outlet with another device to confirm it is functional.
- **No Shockwave Output:** Verify that the control handle cable is firmly connected to the 'OUTPUT' port on the main unit. Ensure a massage head is properly attached to the control handle. Confirm that a treatment program has been selected on the main unit's interface and the green button on the handle has been pressed.
- **Massage Heads Not Attaching Securely:** Ensure you are aligning the threads correctly when attaching massage heads. Do not force them. If a head is damaged due to mis-threading, it may become unusable. Handle with care.
- **Unusual Noises or Malfunctions:** If the device makes unusual noises, emits smoke, or otherwise malfunctions, immediately turn off the power, unplug it, and contact customer support.

## 8. SPECIFICATIONS

---

Feature	Detail
Package Dimensions	16.57 x 13.58 x 4.65 inches
Weight	9.58 Pounds
Manufacturer	zinnor
ASIN	B082XY1GR6
Use for	Whole Body
Power Source	Corded Electric
Brand	zinnor
Color	Blue, Grey, Silver
Special Feature	Rechargeable, Speed Control

## 9. WARRANTY AND SUPPORT

---

zinnor is committed to providing quality products and first-class service. We offer a one-year after-sales service for this product. If you have any questions or require assistance, please contact us directly. We also support a 7-day no-reason refund policy.

If you have not received an electronic version of this manual, please contact us, and we will provide it within 24 hours.