

## Empire Tech B083582KQ8

# Empire Tech LED Eye Device - New Generation User Manual

Model: B083582KQ8

## 1. INTRODUCTION

Thank you for choosing the Empire Tech LED Eye Device - New Generation. This device is designed to provide advanced skincare benefits through LED color therapy, galvanic ion technology, and vibration massage. It aims to support skin regeneration, promote collagen production, and enhance the delivery of nutrients for wrinkle care, whitening, and lifting. Please read this manual thoroughly before using the device to ensure proper operation and maintenance.

## 2. SAFETY INFORMATION

To ensure safe operation and prevent damage to the device or injury to yourself, please observe the following safety precautions:

- Do not use the device on broken, irritated, or infected skin.
- Consult a physician before use if you have any pre-existing medical conditions, especially skin conditions, or if you are pregnant or breastfeeding.
- Keep the device out of reach of children.
- Do not immerse the device in water or any other liquid. Clean with a dry or slightly damp cloth only.
- Do not attempt to disassemble or repair the device yourself. Contact customer support for assistance.
- Use only the charging cable provided with the device.
- Avoid direct eye contact with the LED light.
- Statements regarding this device have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

## 3. PRODUCT OVERVIEW

The Empire Tech LED Eye Device combines multiple technologies for comprehensive skin care:

- **LED Color Therapy:** Utilizes specific light wavelengths to target various skin concerns, aiding in regeneration and collagen production.
- **Galvanic Ion Technology:** Enhances the absorption of skincare products by creating a gentle electrical current that helps nutrients penetrate deeper into the skin.
- **Vibration Massage:** Provides a gentle massage to stimulate blood circulation and relax facial muscles, contributing to a younger and more lifted appearance.



Image: The Empire Tech LED Eye Device, featuring a sleek design with a black body and a gold-colored top and accent. The gold top houses the LED and sensor components.

## 4. SETUP

---

1. **Unboxing:** Carefully remove the device and all accessories from the packaging. Inspect for any visible damage.
2. **Initial Charge:** Before first use, fully charge the device. Connect the charging cable to the device and a suitable power source. The charging indicator (if present) will show charging status.
3. **Cleaning:** Gently wipe the device head with a clean, dry, soft cloth before first use.

## 5. OPERATING INSTRUCTIONS

Follow these steps for optimal use of your Empire Tech LED Eye Device:

1. **Prepare Skin:** Cleanse your face thoroughly and apply your preferred eye cream, serum, or essence to the area you wish to treat. The device works best with a conductive medium.
2. **Activate Device:** The device is designed for convenience and automatically turns on as soon as it makes contact with the skin. There is no manual power button.
3. **Application:** Gently glide the device head over the skin around your eyes, moving in small, circular motions or upward strokes. Focus on areas with fine lines, wrinkles, or where you desire improved absorption of skincare products.
4. **Treatment Duration:** Use for approximately 5-10 minutes per session, or as recommended by your skincare professional.
5. **Deactivation:** The device will automatically turn off when it loses contact with the skin for a prolonged period.
6. **Post-Treatment:** After use, gently wipe the device head clean.

For best results, consistent use is recommended.

## 6. MAINTENANCE AND CARE

Proper care will extend the life of your device:

- **Cleaning:** After each use, wipe the device head with a soft, dry, or slightly damp cloth. Do not use abrasive cleaners, alcohol, or harsh chemicals.
- **Storage:** Store the device in a cool, dry place, away from direct sunlight and extreme temperatures. Keep it in its original packaging or a protective pouch when not in use.
- **Charging:** Recharge the device when the battery is low. Avoid overcharging for extended periods.
- **Avoid Drops:** Handle the device with care to prevent damage from drops or impacts.

## 7. TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on.	Battery is low or depleted. Device not making proper skin contact.	Charge the device fully. Ensure the device head is firmly touching the skin. Apply a conductive skincare product.
Vibration is weak or inconsistent.	Low battery. Insufficient skin contact.	Recharge the device. Ensure full contact with the skin.
Device stops during use.	Loss of skin contact. Battery depleted.	Re-establish firm skin contact. Recharge the device.

If the problem persists after trying these solutions, please contact customer support.

## 8. SPECIFICATIONS

---

- Model:** B083582KQ8
- Brand:** Empire Tech
- Technology:** LED Color Therapy, Galvanic Ion, Vibration Massage
- Power Source:** Rechargeable Battery (details not specified, assume internal)
- Dimensions:** Not specified
- Weight:** Not specified

## 9. WARRANTY INFORMATION

---

Specific warranty details for the Empire Tech LED Eye Device - New Generation are not provided in this manual. Please refer to the product packaging or the retailer's website for information regarding warranty coverage and terms. Keep your purchase receipt as proof of purchase.

## 10. CUSTOMER SUPPORT

---

For technical assistance, troubleshooting beyond this manual, or inquiries about your Empire Tech LED Eye Device, please contact the retailer or manufacturer's customer support channels. Contact information may be found on the product packaging or the official Empire Tech website.

