

Raiuleko R-025

Raiuleko R-025 Combination Body Massager

User Manual

Model: R-025 | Brand: Raiuleko

1. INTRODUCTION

Thank you for choosing the Raiuleko R-025 Combination Body Massager. This device is designed to provide various massage techniques for multi-usage, promoting relaxation and comfort. This manual provides essential information for the safe and effective operation, setup, maintenance, and troubleshooting of your massager. Please read it thoroughly before use and retain it for future reference.

2. SAFETY INFORMATION

No specific safety warnings were provided with this product. However, when using any electrical device, especially one that applies electrical impulses to the body, it is important to observe general safety precautions:

- Do not use if you have a pacemaker, implanted defibrillator, or other electronic medical devices.
- Consult a physician before use if you have any pre-existing medical conditions, including but not limited to heart conditions, epilepsy, or pregnancy.
- Do not apply pads to broken skin, open wounds, or areas with skin irritation.
- Keep out of reach of children.
- Do not use while driving or operating heavy machinery.
- Avoid using near water or in excessively humid environments.
- If you experience any discomfort or adverse reactions, discontinue use immediately.

3. PACKAGE CONTENTS

Upon opening the package, please verify that all components are present:

- 1 x Raiuleko R-025 Main Massager Unit
- 4 x Electrode Pads (self-adhesive)
- 2 x Electrode Connection Cables



Image 1: Raiuleko R-025 massager unit with connected electrode pads and cables.

4. PRODUCT OVERVIEW

4.1 Components

- **Main Unit:** The central control unit with a display screen and buttons.
- **Electrode Pads:** Self-adhesive pads that deliver electrical impulses to the skin.
- **Connection Cables:** Connect the electrode pads to the main unit.

4.2 Controls and Display

Familiarize yourself with the buttons and display on the main unit:

- **ON/OFF Button:** Powers the device on or off.
- **MODE Button:** Cycles through different massage modes (e.g., tapping, kneading, rubbing).
- **TIME Button:** Adjusts the duration of the massage session.
- **STRONG / WEAK Buttons:** Increases or decreases the intensity level of the massage.
- **CH/EN Button:** May switch between channels for different pad outputs or toggle language settings.
- **LCD Display:** Shows current mode, intensity level, and remaining time.

5. SETUP

1. **Install Batteries:** Open the battery compartment on the back of the main unit. Insert three (3) AAA batteries, ensuring correct polarity (+/-). Close the compartment securely.
2. **Connect Electrode Pads:** Plug the connection cables into the designated ports on the top of the main unit. Attach the electrode pads to the other end of the cables.
3. **Prepare Skin:** Ensure the skin where you intend to place the pads is clean, dry, and free of oils or lotions. This helps ensure proper adhesion and conductivity.
4. **Apply Electrode Pads:** Carefully remove the protective film from the electrode pads. Apply the pads firmly to the desired areas of your body. Ensure the pads are not touching each other.

6. OPERATING INSTRUCTIONS

1. **Power On:** Press the **ON/OFF** button to turn on the device. The LCD display will illuminate.
2. **Select Mode:** Press the **MODE** button to cycle through the available massage programs. Observe the display to see the current mode.
3. **Adjust Intensity:** Start with the lowest intensity. Press the **STRONG** button to gradually increase the intensity to a comfortable level. Press the **WEAK** button to decrease intensity. Avoid setting the intensity too high, which may cause discomfort.
4. **Set Time:** Press the **TIME** button to select your desired massage duration. The default time may vary.
5. **Begin Massage:** The device will begin the massage according to your selected mode, intensity, and time.
6. **Power Off:** Once your session is complete, or if you wish to stop, press the **ON/OFF** button to turn off the device.
7. **Remove Pads:** Gently peel the electrode pads from your skin. Reapply the protective film to the pads for storage to maintain their adhesive quality.

7. MAINTENANCE AND CARE

- **Cleaning the Main Unit:** Wipe the main unit with a soft, dry cloth. Do not use abrasive cleaners or immerse the unit in water.
- **Cleaning Electrode Pads:** If the pads lose their stickiness, gently wipe the adhesive surface with a slightly damp cloth (water only). Allow them to air dry completely before reapplying the protective film. Do not use soap or harsh chemicals.
- **Storage:** Store the device and pads in a cool, dry place, away from direct sunlight and extreme temperatures. Always reapply the protective film to the electrode pads after each use to prolong their lifespan.
- **Battery Replacement:** Replace batteries when the display dims or the device's performance decreases. Remove batteries if the device will not be used for an extended period to prevent leakage.

8. TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on	Low or dead batteries; Incorrect battery installation	Replace batteries; Ensure batteries are inserted with correct polarity
No sensation or weak sensation	Pads not properly connected; Pads not adhering well to skin; Intensity too low; Batteries low	Check cable connections; Ensure pads are firmly applied; Increase intensity gradually; Replace batteries
Pads losing stickiness	Pads are dirty; Pads are old/worn out	Gently clean pads with a damp cloth and air dry; Replace with new electrode pads
Skin irritation	Pads applied to sensitive skin; Intensity too high; Allergic reaction	Discontinue use; Reduce intensity; Consult a physician if irritation persists

9. SPECIFICATIONS

Feature	Detail
Product Name	Digital Therapy Machine
Model Number	R-025
Power Source	DC 4.5V (3 x AAA batteries)
Frequency	1-330 Hz
Power Consumption	Less than 180 mW

Feature	Detail
Battery Life (approx.)	45 hours (continuous operation)
Device Dimensions (L x W x H)	153 mm x 70 mm x 18 mm (6.02 x 2.76 x 0.71 inches)
Device Weight	90 g (0.2 lbs)
Total Item Weight	700 g (1.54 lbs)
Material	Mixed Materials
Special Features	Portable
Use For	Body

Name	Digital Therapy Machine
Power	dc4.5v,3aaa(7#)batteries
Frequency	1-330hz
Power consumption	less than 180mw
Life of battery	For 45 hours contionously in normal operation
Size of the device	153mm×70mm×18mm
Weight	90g
Package	1 pcs



The image shows a white, handheld digital therapy massager. At the top, there is a logo for 'Health Herald'. Below the logo is a small LCD screen displaying a human silhouette with two rectangular areas on the lower back, indicating the massage target. Below the screen are several control buttons: an 'ON/OFF' button, a 'MODE' button, and a row of four buttons labeled 'CH/EN', 'TIME', 'STRONG', and 'WEAK'. The bottom of the device has the text 'DIGITAL THERAPY MACHINE' printed on it.

Image 2: Detailed specifications of the Raiuleko R-025 massager.

10. WARRANTY AND SUPPORT

Specific warranty information for the Raiuleko R-025 Combination Body Massager was not provided with this manual. For details regarding warranty coverage, returns, or technical support, please refer to the retailer where the product was purchased or contact the manufacturer directly through their official channels.

