

Manuals.plus /

› Tao Clean /

› TAO Clean Germ Shield UV Toothbrush Sanitizer Instruction Manual

Tao Clean BA-0401-WHT

TAO Clean Germ Shield UV Toothbrush Sanitizer Instruction Manual

INTRODUCTION

The TAO Clean Germ Shield UV Toothbrush Sanitizer is a universal cleaning station designed to kill 99.9% of germs on both manual and electric toothbrush heads using powerful Ultraviolet-C rays. This travel-friendly, battery-operated device ensures your toothbrush is clean and ready for use, promoting better oral hygiene.



Germ Shield

Patented UV-C Sanitizing System • Fast 4 Minute Operation
Accommodates Manual & Electric Brushes • Batteries Included • Travel Friendly



Image: The TAO Clean Germ Shield UV Toothbrush Sanitizer in operation, emitting blue UV-C light to clean a toothbrush, ensuring a superior clean.

KEY FEATURES

- **Universal Compatibility:** Accommodates all manual and electric toothbrushes.
- **Patented Design:** Encloses your toothbrush head for effective UV-C sanitization.
- **Germ Elimination:** Powerful Ultraviolet-C rays kill 99.9% of germs.
- **Battery Operated:** Runs on 2 AA alkaline batteries (included) for convenience.
- **Travel-Friendly:** Compact and portable design for use anywhere.



Image: The TAO Clean Germ Shield, emphasizing its portability and battery-operated nature with 2 AA batteries.

SETUP

1. Locate the battery compartment on the back of the Germ Shield unit.
2. Insert 2 AA alkaline batteries, ensuring correct polarity.
3. Close the battery compartment securely.
4. Place the Germ Shield on a stable, flat surface in your bathroom or desired location.

The Germ Shield runs on 2 AA batteries and is **travel friendly.**



Image: Rear view of the TAO Clean Germ Shield, showing the battery compartment for 2 AA batteries.

OPERATING INSTRUCTIONS

1. After brushing your teeth, rinse your toothbrush thoroughly and shake off any excess water.
2. Open the front door of the Germ Shield UV Sanitizer.
3. Place your toothbrush face down with the bristles facing towards the Germ Shield door. The device is designed to accommodate various brush head sizes.
4. Gently close the door. The sanitizing cycle will activate automatically.
5. A blue light will illuminate around the base of the unit, indicating that sanitization is in progress.
6. The sanitization cycle lasts approximately 4 minutes. After this time, the blue light will automatically shut off, signaling the completion of the cycle.
7. Your toothbrush is now sanitized and ready for its next use.




Image: Step-by-step visual guide on how to insert a toothbrush into the TAO Clean Germ Shield for sanitization.



Image: The blue indicator light of the TAO Clean Germ Shield illuminates during the 4-minute sanitization cycle.

Official Product Videos

Video: Official video demonstrating the ease of use and universal compatibility of the TAO Clean Germ Shield UV Toothbrush Sanitizer with various toothbrush types.

Video: Official product video showcasing the TAO Clean Germ Shield UV Toothbrush Sanitizer, highlighting its design and function.

MAINTENANCE

- **Cleaning:** Regularly clean the interior of the Germ Shield unit to prevent any buildup of toothpaste residue or grime. Use a soft, damp cloth.
- **Battery Replacement:** Replace the 2 AA batteries when the unit no longer activates or the indicator light (if present) signals low power.
- **Drying Toothbrush:** Always shake off excess water from your toothbrush before placing it in the sanitizer to maintain cleanliness and prevent moisture accumulation inside the unit.

TROUBLESHOOTING

- **Unit Not Activating:** If the blue light does not turn on after closing the door, check if the batteries are correctly installed and not depleted. Replace batteries if necessary.
- **Sanitization Cycle Not Completing:** Ensure the door is fully closed and the toothbrush is properly seated within the unit. If the issue persists, try replacing the batteries.

- **Unpleasant Odor/Grime Buildup:** This may occur if toothbrushes are inserted wet. Ensure toothbrushes are shaken dry before placing them in the unit. Clean the interior of the Germ Shield regularly with a damp cloth.

SPECIFICATIONS

- **Product Dimensions:** 4.25 x 4.13 x 4.92 inches
- **Weight:** 11.68 ounces
- **Item Model Number:** BA-0401-WHT
- **Batteries:** 2 AA batteries (included)
- **Manufacturer:** TAO Clean

WARRANTY AND SUPPORT

TAO Clean products come with a **two-year warranty** covering manufacturer-related defects. For support or warranty claims, please contact TAO Clean directly.

- **Email Support:** support@taoclean.com
- **Website:** taoclean.com



Image: Bottom of the TAO Clean Germ Shield packaging, showing contact information and UPC 0859332005722.