

Blu Adult Black-White

Blu Smart Electric Toothbrush User Manual

Model: Adult Black-White

1. INTRODUCTION

This manual provides detailed instructions for the proper use and maintenance of your Blu Smart Electric Toothbrush. Please read this manual thoroughly before using the product to ensure optimal performance and longevity.



Image: The Blu Smart Electric Toothbrush in white, displayed alongside two smartphones showing the interactive live tracking app interface. The app visualizes brushing coverage and provides feedback.

2. WHAT'S IN THE BOX

Upon receiving your Blu Smart Electric Toothbrush package, you will find the following components:

- 1 x Blu Smart Electric Toothbrush Body
- 1 x Wireless Charging Stand
- 2 x Blu Brush Heads



Image: The Blu Smart Electric Toothbrush handle and its wireless charging base. The brush head is detached, showing the connection point.

3. SETUP AND FIRST USE

3.1. Attaching the Brush Head

The brush heads are designed to fit securely onto the top of the toothbrush body. Align the brush head with the metal shaft on the toothbrush handle and push down firmly until it clicks into place.

3.2. Initial Charging

For the first use, place the toothbrush onto the charging stand with the logo side facing outwards. Charge the toothbrush for a full 8 hours to ensure the battery is fully conditioned. After the initial charge, your toothbrush will last up to one month with regular use (brushing twice a day for the recommended two minutes per session).

The charging stand is USB compatible, allowing for convenient charging via a computer or any USB-compatible device, which is ideal for travel.

3.3. Downloading the Blu App

The Blu Smart Electric Toothbrush utilizes an interactive app for enhanced brushing experience and tracking. The app is compatible with both iOS and Android devices.

1. Open your device's App Store (for iOS) or Google Play Store (for Android).
2. Search for "Blu Toothbrush".
3. Locate the app with the Blu logo and download the latest version.

3.4. App Registration and Device Pairing

Upon opening the Blu app for the first time, you will be prompted to sign in or register. You can sign in using your Apple ID, Facebook account, or via email. After successful registration, you will be directed to the home screen. To connect your toothbrush to the app:

1. Navigate to the "Settings" icon located at the lower right-hand corner of the app.
2. Select "Device Management".
3. Click the "Pair Device" button. Ensure your Blu toothbrush is nearby.
4. Any available Blu devices will appear in the device list. Select your toothbrush to pair. If you have multiple brushes (e.g., for family members), ensure you select the correct one and name it accordingly.

Your browser does not support the video tag.

Video: This official video from Blu provides a step-by-step guide on how to set up your Blu toothbrush, including attaching brush heads, initial charging, downloading the app, and pairing your device.

4. OPERATING INSTRUCTIONS

4.1. Brushing Modes

The Blu Smart Electric Toothbrush offers five distinct brushing modes for a tailored experience:

- **Clean:** Standard daily cleaning.
- **White:** For surface stain removal and whitening.
- **Polish:** For polishing and brightening teeth.
- **Soft:** Gentle cleaning for sensitive teeth and gums.
- **Sensitive:** Extra gentle cleaning for highly sensitive areas.

4.2. Using the App for Live Tracking

To begin a brushing session with live tracking:

1. Ensure your smartphone is available and your toothbrush is within reach.
2. From the app's dashboard, press the large "BRUSH" circle button in the center.
3. Select the user who will be brushing (e.g., "You" or a child's named profile).
4. The app will display a visual representation of your teeth and guide you through brushing different zones of your mouth. It will show where you should be brushing along the gum lines and in specific zones.

4.3. Understanding Brushing Feedback

After completing a brushing session, the app will provide a percentage of coverage. Any areas that were missed during brushing will be highlighted in blue on the teeth diagram. You have a 10-minute window to go back and correct these missed areas by re-brushing them.

4.4. Family Connect (Connecting with Kids' Brushes)

The Blu app allows you to monitor the brushing habits of your children. To connect a child's brush:

1. From the dashboard, press the "Connect With Your Kids" button.
2. The app will scan for nearby Blu brushes.
3. Select the child's brush from the list and assign a name to it (e.g., "Thomas"). This allows for individual tracking.

4.5. Brushing History and Statistics

The Blu app stores brushing data for up to 45 days without needing to be connected to the app. You can review past brushing sessions and overall statistics:

- **History:** Click the "History" button at the bottom of the app. You can view a calendar where days with brushing activity are circled in blue. Clicking on a specific date will show the brushing chart for that day, including duration and coverage.
- **Statistics:** The "Statistics" screen provides an overview of your brushing habits, including duration, frequency (e.g., 1.8 times/day), and average coverage percentage.

5. MAINTENANCE

5.1. Brush Head Replacement

It is recommended to replace your Blu brush head every three months. The blue bristles on the brush head will fade to white over time, indicating that it's time for a replacement. The app also features a "Brush Head Monitor" that tracks usage and will remind you when to order new brush heads.

5.2. Cleaning the Toothbrush

Regularly clean your toothbrush handle and brush head with water after each use. Ensure the charging stand is kept dry.

6. SPECIFICATIONS

Feature	Detail
Brand	Blu
Model Number	Adult Black-White
Power Source	Corded Electric (USB compatible)
Brushing Modes	5 (Clean, White, Polish, Soft, Sensitive)
Brush Strokes per Minute	40,000
Timer	Auto 2-minute timer
Waterproof Rating	IPX7
Battery Life	Up to 1 month per full charge
Bristle Firmness	Soft

App Compatibility	iOS and Android
Dimensions	8.7 x 6.3 x 1.5 inches
Weight	1 Pound

7. TROUBLESHOOTING

If you encounter any issues with your Blu Smart Electric Toothbrush or the accompanying app, please refer to the app's "Home Care Instructions" section under settings for common solutions. For further assistance, please contact customer support.

8. IMPORTANT INFORMATION

8.1. Legal Disclaimer

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

9. SUPPORT

For any questions or assistance regarding your Blu Smart Electric Toothbrush, please feel free to email our support team.

Contact Email: Support@Bluetoothbrush.com

You can also visit our website for more information: www.bluetoothbrush.com