

[Manuals.plus](#) /

> [CINCOM](#) /

> CINCOM Hand Massager Instruction Manual

CINCOM CM-026H

CINCOM Hand Massager Instruction Manual

Model: CM-026H

PRODUCT OVERVIEW

The CINCOM Hand Massager is a cordless device designed to provide relief from hand pain, arthritis, and carpal tunnel syndrome through a combination of heat and compression massage. It features multiple massage modes, intensity levels, and heat settings for a customizable and soothing experience.



Image Description: A white CINCOM hand massager with a hand inserted, showing the top control panel with power, intensity, mode, and heat buttons.

SAFETY INFORMATION

Important Safety Precautions:

- Do not use the massager on areas of the hand with open wounds, bruises, or swollen/inflamed skin.
- Consult a healthcare professional before use if you have any pre-existing medical conditions.
- Keep out of reach of children.
- Do not disassemble or attempt to repair the device yourself.
- Use only the provided charging cable.

WHAT'S IN THE BOX

Your CINCOM Hand Massager package includes the following items:



Image Description: An illustration showing the contents of the product box: 1. CINCOM Hand Massager, 2. Product Box, 3. Power Adapter, 4. USB-C Charging Cable, 5. User Manual.

- 1 x CINCOM Hand Massager
- 1 x USB-C Charging Cable
- 1 x Power Adapter
- 1 x User Manual

SETUP AND CHARGING

Initial Charge: Before first use, fully charge the hand massager.

1. Connect the provided USB-C charging cable to the charging port located at the back of the device.
2. Plug the other end of the USB-C cable into the power adapter, and then into a standard wall outlet.
3. The indicator light will show charging status (e.g., blinking while charging, solid when full).
4. A full charge typically takes 4-5 hours and provides 4-5 hours of use.

What's in the Box?

1



3



4



2



5



Image Description: A close-up view of the USB-C charging port located on the back of the hand massager.

OPERATING INSTRUCTIONS

The CINCOM Hand Massager is designed for ease of use with intuitive controls.

Control Panel Overview



✓ Stimulate
Muscles

✓ Ease
Soreness

✓ Increase
Blood Flow

✓ Relieve
Pain

Image Description: A top-down view of the hand massager's control panel, highlighting the Power, Intensity, Mode, and Heat buttons.

- **Power Button (U):** Press and hold for 3 seconds to turn the device ON or OFF.
- **Intensity Button (|||):** Cycles through 3 intensity levels (Low, Medium, High) for air compression.
- **Mode Button (m):** Cycles through 3 massage modes (Soothing, Vitality, Relaxation).
- **Heat Button (Heat):** Toggles between 2 heat settings (Low, High) or turns heat OFF.

Using the Massager

1. **Insert Hand:** Place your hand into the massager, ensuring your fingers are comfortably positioned. The design allows for both left and right hand use, with thumb openings on either side.
2. **Power On:** Press and hold the Power button for 3 seconds.

3. **Select Mode:** Press the Mode button to choose your desired massage program.
4. **Adjust Intensity:** Press the Intensity button to select your preferred compression strength.
5. **Activate Heat:** Press the Heat button to turn on the warming function and cycle through heat levels.
6. **Enjoy:** The massager will automatically shut off after 15 minutes for safety and optimal session duration.

Easy to Control



Intensity Button

3 intensity levels



Mode Button

3 mode programs



Power Button

Press it for 3 seconds



Heat Button

2 heating levels

Image Description: An illustration showing the internal air compression and heat elements targeting different parts of the hand.

Endurance Battery Life



Image Description: An illustration demonstrating how the massager's internal structure optimizes massage for different hand positions and pressure points.

Bigger Size for Your Hands



Non-slip Rubber Feet



Comfy Soft Interior



USB-C Charging Port

Image Description: A visual representation of the two heat levels (Low ~107°F, High ~122°F) and their effect on blood circulation, with a user happily using the device.

MAINTENANCE AND CARE

- **Cleaning:** Wipe the exterior of the massager with a soft, damp cloth. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the massager in a cool, dry place away from direct sunlight and extreme temperatures. Use the provided storage bag for protection.
- **Battery Care:** To prolong battery life, avoid fully discharging the battery frequently. Charge the device regularly, even if not in use for extended periods.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on	Low battery; Power button not pressed correctly	Charge the device fully. Press and hold the Power button for 3 seconds.
Weak compression	Intensity level too low	Increase the intensity level using the Intensity button.
No heat	Heat function is off or set to low	Press the Heat button to cycle through heat levels (Low, High).
Device stops unexpectedly	Automatic shut-off activated (15-minute timer)	This is normal operation. Restart the device if you wish to continue the massage.

SPECIFICATIONS

Feature	Detail
Model Number	CM-026H
Product Dimensions	7.8 x 4.5 x 3.1 inches
Item Weight	1.5 Pounds
Power Source	Battery Powered (Built-in 2500mAh Lithium Battery)
Heat Settings	2 levels (95 - 104 °F / 35 - 40 °C)
Massage Modes	3 selectable modes
Intensity Levels	3 selectable strengths
Automatic Timer	15 minutes
Material	Plastic
Manufacturer	CINCOM

WARRANTY AND SUPPORT

CINCOM Healthcare products are known for their quality and excellent customer service.

- **Warranty Information:** Please refer to the product packaging or the official CINCOM website for detailed warranty terms and conditions.
- **Customer Support:** If you encounter any problems while using the hand massager, or have any questions, please contact CINCOM customer support. Contact details can typically be found in the user manual or on the manufacturer's website.
- **FSA/HSA Eligibility:** This product may be eligible for reimbursement through Flexible Spending Accounts (FSA) or Health Savings Accounts (HSA). Consult your plan administrator for details.