

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

› [Boncare](#) /

› [Boncare 10 Seconds Digital Thermometer for Adults and Kids - Accurate & Safe Rectal Thermometer for Baby with Fever Alarm, Oral Thermometer FSA HSA Eligible \(Dark Blue\)](#)

### Boncare GF-MT501

# Boncare 10 Seconds Digital Thermometer User Manual

MODEL: GF-MT501

Brand: Boncare

## 1. Introduction

---

The Boncare 10 Seconds Digital Thermometer is designed for accurate and fast temperature readings for adults, children, and infants. It is suitable for oral, rectal, or axillary (underarm) use, providing a versatile solution for your family's health monitoring needs. This thermometer features a soft silicone tip and a reliable stainless-steel probe, ensuring both comfort and precision. It is FSA and HSA eligible, making it a convenient choice for personal healthcare.



Image: Boncare 10 Seconds Digital Thermometer in dark blue, shown with its retail packaging.

## 2. Important Safety Information

**Legal Disclaimer:** Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This product is a medical device and should be used according to the instructions provided.

- Before using the health care thermometer, carefully read the product manual.
- To ensure the accuracy of the measurement, avoid touching the probe with your hands or other hot

objects prior to taking the temperature.

- The thermometer is designed for oral, rectal, or axillary use.
- Keep out of reach of small children.



Image: The thermometer is suitable for a wide range of users, including babies, children, adults, and even for veterinary use.

### 3. Product Overview and Key Features

---

The Boncare Digital Thermometer is equipped with several features to ensure ease of use and accurate readings:

- **Soft Silicone Tip:** Provides comfort and safety during use.
- **Stainless-Steel Probe:** Ensures reliable and accurate temperature detection.
- **LCD Screen Display:** Highly legible for easy reading of results.
- **Beeper Alarm:** Notifies you when the reading is complete or if a fever is detected.
- **Memory Display Function:** Recalls the last reading for comparison.
- **Auto-Shut-Off Mechanism:** Conserves battery life.
- **Fahrenheit/Celsius Switch:** Easily switch between °F and °C units.

# High Accuracy

Measurement with  $\pm 0.2^{\circ}\text{F}$



Image: Key components of the Boncare Digital Thermometer, including the ON/OFF button, battery cap, LCD display, memory recall, and temperature sensing probe.

## 4. Setup and Initial Use

Your Boncare Digital Thermometer comes with one LR41 battery pre-installed. To begin using your thermometer:

1. Remove the thermometer from its protective storage case.
2. Press the power button (ON/OFF) to turn on the device. A short beep will sound, and the display will show the last measured temperature for approximately 2 seconds.
3. After the last reading is displayed, the screen will show 'Lo°F' (or 'Lo°C') indicating it is ready for measurement.

## With Storage Case Hygienic and Portable



Image: The thermometer is lightweight and portable, fitting easily into its included storage case for convenience at home or on the go.

### 5. Operating Instructions

The thermometer provides accurate readings in approximately 10 seconds. Ensure the probe is clean before use.

#### Taking a Temperature Reading:

1. Turn on the thermometer by pressing the power button.
2. Place the probe according to the desired measurement method (oral, rectal, or axillary).
3. Hold the thermometer in place until the beeper alarm sounds, indicating the measurement is complete.
4. Read the temperature displayed on the LCD screen.
5. Press the power button again to turn off the thermometer, or it will automatically shut off after approximately 10 minutes of inactivity.

#### Measurement Methods:

- **Oral Use:** Place the probe tip under the tongue. Close your mouth and hold the thermometer in place. Normal oral temperature is typically 96.3°F-99.2°F (35.7°C-37.3°C).

- **Rectal Use:** This method is generally recommended for infants and young children. Gently insert the probe tip into the rectum. Normal rectal temperature is typically 96.3°F-99.2°F (35.7°C-37.3°C).
- **Axillary (Underarm) Use:** Place the probe tip in the armpit, ensuring good skin contact. Hold the arm firmly against the body. Normal underarm temperature is typically 95.4°F-98.6°F (35.2°C-37.0°C).

# Clinical Testing More Convincing

Suitable for **Babies Kids and Adults**




Image: Visual guide demonstrating the oral, rectal, and underarm methods for taking temperature with the Boncare thermometer.

## Switching Temperature Units (°F/°C):

To switch between Fahrenheit (°F) and Celsius (°C), press and hold the power button for approximately 3-5 seconds while the thermometer is off. The display will change units. Release the button to select the desired unit.



Image: Instructions on how to toggle between Fahrenheit and Celsius temperature units on the thermometer's display.

## Fever Alarm:

The thermometer features a fever alarm. If a temperature of 100°F (37.8°C) or higher is detected, the thermometer will emit a series of rapid beeps to alert you.



Image: The thermometer features a fever alarm that alerts the user to elevated temperatures.

### **Official Product Videos:**

No official product videos from the seller are available for embedding at this time.

## **6. Care and Maintenance**

---

### **Cleaning:**

The waterproof tip of the thermometer makes hygienic cleaning simple. You can wipe it with water and clean it with medical alcohol either before or after each use. Do not immerse other parts of the thermometer in water.

# Easy to Clean

## Waterproof Tip



**Note :** The tip part is waterproof.  
(Do not immerse other parts in water)

Image: The waterproof tip allows for easy and hygienic cleaning with water or medical alcohol.

### **Battery Replacement:**

The thermometer requires 1 LR41 battery. When the battery is low, a battery symbol will appear in the lower right corner of the LCD display. To replace the battery, carefully remove the battery cap and insert a new LR41 battery, ensuring correct polarity.

# Low Battery Alarm

Replace Battery when “” appears in the lower right corner of LCD display. “Lo” does not mean the battery is low.



Image: The low battery alarm indicates when the LR41 battery needs to be replaced. The 'Lo' reading does not indicate a low battery.

## Storage:

Store the thermometer in its provided storage case in a cool, dry place, away from direct sunlight and extreme temperatures.

## 7. Troubleshooting

---

- **"Lo" Display:** If the display shows "Lo°F" (or "Lo°C") after turning on, it indicates that the thermometer is ready for measurement and the ambient temperature is below the measurable range. This does not mean the battery is low.
- **Inaccurate Readings:** Ensure the probe is clean and properly positioned for the chosen measurement method. Avoid touching the probe with hands or hot objects before measurement. If issues persist, try replacing the battery.
- **No Display/Not Turning On:** Check if the battery is correctly installed or needs replacement.

## 8. Product Specifications

---

Feature	Detail
Product Dimensions	4.93 x 0.76 x 0.49 inches; 0.37 ounces
Item Model Number	GF-MT501
Batteries	1 LR41 batteries required (included)
Manufacturer	Boncare
Special Feature	Low Battery Indicator, Alarm, Fast Reading System, Auto Off, High Accuracy
Color	Dark Blue
Age Range (Description)	Newborn, infant, toddler, kid, youth, adult
Included Components	Thermometer
Outer Material	Stainless Steel
Display Type	Digital
Product Care Instructions	Wipe Clean

## 9. Warranty and Support

---

Boncare is committed to providing professional and high-quality personal care products. Your Boncare Digital Thermometer comes with a **1-Year Limited Warranty**.

For customer service support or any inquiries, please refer to the contact information provided in your product packaging or visit the official Boncare store online.

An additional [User Guide \(PDF\)](#) is available for download.