

## CINCOM CM028A

# CINCOM Cordless Calf Massager User Manual

Model: CM028A | Brand: CINCOM

---

## INTRODUCTION

---

The CINCOM Cordless Calf Massager (Model CM028A) is designed to provide air compression massage for the calves, feet, and arms. It features multiple modes and intensity levels to help improve circulation, relieve muscle fatigue, and reduce stress. This manual provides essential information for safe and effective use of your device.



*Image: CINCORM Cordless Calf Massager, highlighting its compact and portable design.*

## **SAFETY INFORMATION**

---

Keep out of reach of children. For external use only. Do not use if skin is irritated or injured. Consult a physician before use if you have any pre-existing medical conditions.

## **PACKAGE CONTENTS**

---

- 2 x Massager Wraps (Left and Right)
- 1 x USB-C Charging Cable (Dual USB-C)
- 1 x Wall Adapter
- 1 x Storage Bag
- 1 x User Manual



*Image: All items included in the CINCOM Cordless Calf Massager package.*

## SETUP AND CHARGING

---

Before first use, fully charge both massager wraps.

1. Connect the dual USB-C charging cable to the charging ports on each massager wrap's control panel.
2. Plug the USB-A end of the cable into the provided wall adapter.
3. Plug the wall adapter into a power outlet. The indicator light on the control panel will illuminate to show charging status.
4. Once fully charged (approximately 2.5-3.5 hours for full usage), the indicator light will turn off or change color. Disconnect the charger.

The CINCOM logo is displayed in a stylized, white, italicized font with a registered trademark symbol. The letter 'O' in 'COM' is highlighted with a red-to-white gradient. The background of the entire page is a dark, moody photograph of two women in athletic wear, one of whom is wearing a massage wrap on her arm.

**CINCOM**<sup>®</sup>

***All For The Healthcare***

Discover CINCOM's  
New Release-Advanced  
Relief Meets  
Stylish Comfort

*Image: The USB-C charging port is located on the side of the control panel. Use the provided dual USB-C cable to charge both massager units simultaneously.*


## **WEARING INSTRUCTIONS**

The massager wraps are designed for use on calves, feet, or arms. Ensure correct orientation (Left/Right) as indicated on the wraps.

1. Unzip the massager wrap completely.
2. Place your leg (or arm/foot) inside the wrap.
3. Adjust the Velcro straps to ensure a snug but comfortable fit. Do not overtighten, as the air compression requires some space to inflate.
4. Alternatively, for easier application, open the Velcro straps fully and wrap the unit around your leg, then

secure the Velcro.

## Long Battery Life

 2500mAh\*2 Battery Capacity

 Type-C Fast Charging

 2.5–3.5H Use per Full Charge



Battery  
**2500mAh**

*Image: A person wearing the CINCOS calf massager wraps on both calves, demonstrating proper placement.*



*Image: Step-by-step guide for wearing the massager: 1. Controller outward, wrap with controller facing the outer side of your calf. 2. Gentle Velcro fit, don't fasten too tightly on first use. 3. Power on & Adjust, slightly*

inflate and adjust for comfort. 4. Easy Zip-Off, unzip to take it off effortlessly.

## OPERATING INSTRUCTIONS

---

Each massager wrap has an independent control panel.

1. **Power On/Off:** Press and hold the Power button (circular icon) for 3 seconds to turn the massager on or off.
2. **Select Mode:** Press the 'M' button to cycle through the 3 available massage modes:
  - Sequence Mode: Inflates and deflates chambers sequentially.
  - Circulation Mode: Inflates and deflates chambers in a circulating pattern.
  - Whole Mode: Inflates and deflates all chambers simultaneously.
3. **Adjust Intensity:** Press the Intensity button (bars icon) to cycle through 3 intensity levels (low, medium, high).
4. **Pause Function:** Press the 'PAUSE' button to temporarily stop the massage. Press again to resume.
5. The device has a 20-minute auto shut-off feature for safety and convenience.

*Video: This video demonstrates the application and basic operation of the CINCOM Rechargeable Leg Massager, showing how to wear it and use the control panel.*

## MAINTENANCE AND STORAGE

---

Proper care ensures the longevity of your massager.

- **Cleaning:** Wipe the surface of the massager wraps with a soft, damp cloth. Do not immerse in water or use harsh cleaning agents.
- **Storage:** After use, ensure the massager is turned off and fully charged. Store it in the provided storage bag in a cool, dry place away from direct sunlight.
- Do not store with heavy objects on top to avoid damage to the air bladders.

## TROUBLESHOOTING

---

If you encounter any issues, please refer to the following:

- **Device not turning on:** Ensure the massager is fully charged.
- **Weak compression:** Check if the wraps are securely fastened but not too tight. Adjust intensity level.
- **Uneven compression:** Ensure the wraps are applied evenly and correctly to both legs.
- If problems persist, contact CINCOM customer support.

## SPECIFICATIONS

---

Feature	Detail
Model Number	CM028A
Product Dimensions	14 x 8 x 3 inches
Item Weight	3 Pounds

Feature	Detail
Power Source	Battery Powered
Material	Nylon
Battery Capacity	2500mAh*2 (included)
Usage Time (Full Charge)	2.5-3.5 hours
Auto Shut-off	20 minutes
Modes	3 (Sequence, Circulation, Whole)
Intensity Levels	3 (Low, Medium, High)

## WARRANTY AND SUPPORT

---

CINCOM provides dedicated customer service. For any inquiries or support needs, please refer to the contact information provided in your product packaging or visit the official CINCOM website.

*Legal Disclaimer: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.*