

[Manuals.plus](#) /

› [Signstek](#) /

› Signstek High Frequency Facial Massager with 6 Electrodes Instruction Manual

Signstek A

Signstek High Frequency Facial Massager Instruction Manual

Model: A

INTRODUCTION

Thank you for choosing the Signstek High Frequency Facial Massager. This device is designed to enhance skin appearance through high-frequency electrical currents, promoting cellular metabolism, reducing fine lines, and assisting in acne treatment. Please read this manual thoroughly before use to ensure safe and effective operation.

IMPORTANT SAFETY INFORMATION

WARNING: Improper use can lead to injury. Always follow these safety guidelines:

- Do not use if you have a pacemaker, heart condition, epilepsy, or are pregnant.
- Avoid contact with water. Do not use in wet environments.
- Keep out of reach of children.
- Do not use on broken, irritated, or infected skin.
- Always ensure the electrode is fully inserted into the handle before turning on the device.
- Start with the lowest intensity setting and gradually increase as comfortable.
- Do not leave the device unattended while in operation.
- If you experience any discomfort or adverse reactions, discontinue use immediately.

PACKAGE CONTENTS

Your Signstek High Frequency Facial Massager package includes:

- 1 x High Frequency Handle
- 6 x Glass Electrodes (Mushroom, Tongue, Bent, Comb, Y-shape, Rod)
- 1 x Instruction Manual



Image: The Signstek High Frequency Facial Massager handle with all six included glass electrodes.

PRODUCT FEATURES

The Signstek High Frequency Facial Massager utilizes both Argon and Neon gas-filled electrodes to provide various skin benefits:

- **Argon Electrodes (Violet Light):** Primarily used for acne treatment due to their anti-bacterial properties. Helps to reduce inflammation and promote healing.
- **Neon Electrodes (Orange/Red Light):** Designed for anti-aging benefits, including reducing fine lines, wrinkles, and improving blood circulation for skin rejuvenation.
- **6 Specialized Electrodes:** Each electrode is shaped for specific facial and body areas, ensuring targeted treatment.
- **Adjustable Intensity:** Allows customization of treatment intensity for comfort and effectiveness.
- **Portable Design:** Easy to handle and suitable for home use.



Image: Diagram explaining how high-frequency treatment works, showing instant stimulation, thermal energy, and enriched oxygen for anti-bacterial effects.

SETUP

1. Ensure the device is unplugged and turned off before inserting any electrode.
2. Select the desired glass electrode for your treatment.
3. Gently insert the glass electrode into the high-frequency handle. Push it in firmly until it clicks into place. Do not force it.
4. Plug the device into a standard electrical outlet.

OPERATING INSTRUCTIONS

General Usage Steps:

1. Cleanse your skin thoroughly and apply a thin layer of serum or moisturizer if desired. This helps the electrode glide smoothly and enhances product absorption.
2. Insert the chosen electrode into the handle (refer to Setup).
3. Turn the device on by rotating the intensity dial. Start at the lowest setting.
4. Place your finger on the glass electrode before touching it to your skin. This grounds the current and prevents a sudden shock.
5. Remove your finger and gently place the electrode on the treatment area. You will see a subtle glow and hear a buzzing sound.
6. Move the electrode in circular motions or glide it across the skin for 3-5 minutes per area. For spot treatment, hold the electrode over the blemish for a few seconds.
7. To end the treatment, place your finger back on the electrode, then remove it from your skin. Turn off the device and unplug it.
8. Remove the electrode by gently pulling it out.

Electrode Usage Guide:

The device includes 6 different electrodes, each designed for specific applications:

Power of Argon

anti-bacterial for acne treatment

BENT

Bent tube is for spot area treatment to reduce pimples.



ROD

Rod tube is for body, forehead and other broader areas.



TONGUE

Tongue tube is for sensitive areas, such as dark circles beneath eyes, nose or lips.

Image: Illustration of Neon gas electrodes (Comb, Y-shape, Mushroom) and their primary uses for blood circulation and anti-aging.

- **Mushroom Tube (Neon):** Ideal for broad area treatment on the face, such as cheeks and forehead, to reduce wrinkles and fine lines.
- **Y-Shape Tube (Neon):** Used for toning skin around the neck, arms, and legs.
- **Comb Tube (Neon):** Designed to stimulate blood circulation in the scalp, promoting hair follicle health.



Image: Illustration of Argon gas electrodes (Rod, Bent, Tongue) and their primary uses for anti-bacterial acne treatment.

- **Rod Tube (Argon):** Suitable for broader areas of the body and forehead for general anti-bacterial treatment.
- **Bent Tube (Argon):** Specifically designed for spot treatment of pimples and blemishes.
- **Tongue Tube (Argon):** Perfect for sensitive areas such as under-eye circles, nose, or lips.

Recommended Treatment Frequency:

For optimal results, use the device 2 to 3 times per week. Allow your skin adequate time to respond to the high-frequency stimulation.

MAINTENANCE

- **Cleaning Electrodes:** After each use, wipe the glass electrodes with a soft cloth dampened with rubbing alcohol. Allow them to air dry completely before storing.
- **Cleaning Handle:** Wipe the handle with a dry or slightly damp cloth. Do not immerse the handle in water.
- **Storage:** Store the device and electrodes in a cool, dry place, away from direct sunlight and moisture. Keep them in

their original packaging or a protective case to prevent damage to the glass electrodes.

- **Handling:** Handle glass electrodes with care as they are fragile.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on.	Not plugged in; Power outlet issue; Device malfunction.	Check power connection. Try a different outlet. If still not working, contact customer support.
No light or buzzing from electrode.	Electrode not fully inserted; Intensity too low; Faulty electrode.	Ensure electrode is firmly seated. Increase intensity. Try a different electrode.
Sparkling when touching skin.	Electrode not grounded before contact; Skin too dry.	Always touch the electrode with your finger before applying to skin. Ensure skin is clean and moisturized.
Electrode broke during use/insertion.	Excessive force; Accidental impact.	Glass electrodes are fragile. Handle with care. Contact customer support for replacement options.

SPECIFICATIONS

- **Model Number:** A
- **Brand:** Signstek
- **Power Source:** Corded Electric
- **Number of Electrodes:** 6
- **Package Dimensions:** 32.4 x 21.7 x 5.7 cm
- **Weight:** 810 grams

WARRANTY AND SUPPORT

For any questions regarding your Signstek High Frequency Facial Massager, please do not hesitate to contact our customer service team. Refer to your purchase documentation for specific warranty details.

Customer Service Contact: Please refer to the contact information provided with your purchase or visit the official Signstek website.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question