

Revitive Essential Circulation Booster

Revitive Essential Circulation Booster User Manual

Model: Essential Circulation Booster (5827-RLV-CA)

INTRODUCTION

The Revitive Essential Circulation Booster is a drug-free medical device designed to provide relief from milder occasional symptoms associated with a sedentary or inactive lifestyle. It utilizes professional-strength OxyWave technology to stimulate muscles in the legs and feet, aiming to alleviate tired, aching legs and feet, and help maintain leg health.

IMPORTANT SAFETY INFORMATION

Please read all safety warnings and instructions before using the Revitive Essential Circulation Booster. Failure to do so may result in injury or damage to the device.

- **Do not use if:**
 - You are fitted with an electronic implant, such as a pacemaker or AICD.
 - You are being treated for, or have the symptoms of, a Deep Vein Thrombosis (DVT).
 - You are pregnant.
- Always consult your doctor before using this device if you have any medical conditions or concerns.
- Do not rely solely on the information presented here; always read labels, warnings, and directions provided with the product packaging.

PACKAGE CONTENTS

The Revitive Essential Circulation Booster package typically includes:

- Revitive Essential Circulation Booster unit
- Power adapter
- Instruction Manual

SETUP

1. **Unpack the Device:** Carefully remove the Revitive Essential Circulation Booster and all accessories from the packaging.
2. **Position the Device:** Place the unit on a flat, stable surface in front of a comfortable chair where you can sit and relax.
3. **Connect Power:** Insert the power adapter into the power port on the side of the Revitive unit, then plug the adapter into a wall outlet. Ensure the device is charged if necessary before first use.



Figure 1: The Revitive Essential Circulation Booster main unit, showing the footpads and central control panel.

OPERATING INSTRUCTIONS

Follow these steps for effective use of your Revitive Essential Circulation Booster:

1. **Sit Comfortably:** Sit down in a chair with your feet able to comfortably reach the device.
2. **Place Feet:** Place your bare feet on the footpads of the Revitive unit. Ensure full contact for optimal stimulation.
3. **Power On:** Press the power button on the device to turn it on.
4. **Adjust Intensity:** Use the '+' and '-' buttons on the central control panel to gradually increase the intensity level. Start at a low level and increase until you feel a comfortable, strong muscle contraction in

your calves. The device offers 1-50 intensity levels.

5. **Relax and Enjoy:** Allow the device to work for the recommended session duration (typically 20-30 minutes). You can read, watch TV, or engage in other seated activities during this time.
6. **Post-Treatment:** After each session, allow adequate time for your muscles to recover, similar to after physical exercise.



Figure 2: A user demonstrating the correct posture for using the Revitive Essential Circulation Booster, with bare feet placed on the footpads.

REVITIVE

CIRCULATION BOOSTER

ESSENTIAL

- ▶ Relieves tired, aching legs and feet
- ▶ Helps maintain leg health
- ▶ Actively improves circulation



Figure 3: Front view of the Revitive Essential, highlighting its features including OxyWave technology and 1-50 intensity levels.



Figure 4: Illustrative diagram demonstrating the patented Wide-Pulse Stimulation technology targeting leg muscles.



Figure 5: Diagram illustrating how the device helps maintain leg vein health and relieves tired, aching legs and feet.

MAINTENANCE

To ensure the longevity and proper functioning of your Revitive Essential Circulation Booster, follow these simple maintenance guidelines:

- **Cleaning:** After each use, wipe the footpads and the main unit with a clean, damp cloth. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the device in a cool, dry place away from direct sunlight and extreme temperatures.
- **Inspection:** Periodically inspect the power cord and the device for any signs of damage. If damage is found, discontinue use and contact customer support.

TROUBLESHOOTING

If you encounter issues with your Revitive Essential Circulation Booster, consider the following:

- **No Power:**
 - Ensure the power adapter is securely plugged into both the device and a working wall outlet.
 - Check if the power outlet is functional by plugging in another appliance.
- **No Sensation/Weak Sensation:**

- Ensure your bare feet are making full contact with both footpads.
- Increase the intensity level gradually using the '+' button.
- Ensure your feet are clean and moisturized, as dry skin can reduce conductivity.

- **Device Stops During Use:**

- Check the power connection.
- The device may have an automatic shut-off feature after a set duration; simply restart if needed.

If problems persist, please contact customer support.

SPECIFICATIONS

Model Name	Essential Circulation Booster
Item Model Number	5827-RLV-CA
Product Dimensions	61.5 x 38 x 41.5 cm
Weight	1.74 kg
Power Source	Corded Electric
Material	Acrylonitrile Butadiene Styrene
Colour	White, Green, and Black
Manufacturer	Actegy Limited
Country of Origin	Korea, Republic of

WARRANTY AND SUPPORT

The Revitive Essential Circulation Booster comes with a **2-year manufacturer's warranty**.

For warranty claims, technical support, or further inquiries, please contact Actegy Limited customer service. Refer to the product packaging or the official Revitive website for the most current contact information.

OFFICIAL PRODUCT VIDEOS

No official product videos from the seller were provided in the product data for embedding in this manual.