



[Manuals.plus](#) /

> [MiSMON](#) /

> MiSMON High Frequency Facial Machine (Model OLB-BP-001a) Instruction Manual

MiSMON OLB-BP-001a

MiSMON High Frequency Facial Machine (Model OLB-BP-001a) Instruction Manual

Your guide to effective and safe use of your MiSMON facial device.

1. INTRODUCTION

The MiSMON High Frequency Facial Machine is a versatile skincare device designed for home use. It integrates multiple functions to support various skin care needs, including cleansing, lifting, product absorption, anti-aging, and acne treatment. This manual provides detailed instructions for proper operation and maintenance to ensure optimal results and device longevity.



Image 1: MiSMON High Frequency Facial Machine with its charging base.

2. PRODUCT OVERVIEW

The device features a compact, ergonomic design with an intuitive interface. It combines Radio Frequency (RF), Electrical Muscle Stimulation (EMS) micro-current, acoustic vibration, and 5 LED light therapies (Green, Yellow, Pink, Red, Blue) to address various skin concerns.

Key Components:

- **Probe for EMS/RF:** The head of the device containing the metal contacts for micro-current and radio frequency.
- **LED Light:** Integrated lights on the probe for various light therapy modes.
- **LED Screen:** Displays current mode, level, and battery status.
- **Level Button:** Adjusts the intensity level of the selected mode.
- **Mode Button:** Cycles through the different operational modes.
- **Power Switch:** Turns the device on or off.



Image 2: Labeled diagram of the MiSMON Facial Machine's main components.

3. SETUP

3.1 Charging the Device

1. Remove the facial machine from its packaging.
2. Connect the USB charging cable to the charging base and a suitable power adapter.
3. Place the device onto the charging base. The LED screen will indicate charging status.
4. Fully charge the device before its first use. A full charge typically takes several hours. The battery indicator on the screen will show full bars when charging is complete.

3.2 Initial Preparation

- Ensure your face is clean and free of makeup before use.
- For optimal results and to ensure proper function of RF and EMS, apply a water-based gel, lotion, essence, or facial mask to your skin.



Image 3: The device is USB rechargeable and designed for portability.

4. OPERATING INSTRUCTIONS

The MiSMON Facial Machine offers 5 distinct modes, each combining different technologies for specific skincare benefits. Each mode is recommended for 5-10 minutes per session.

4.1 Powering On/Off

- To power on, press and hold the 'MODE' button for a few seconds until the screen illuminates.
- To power off, press and hold the 'MODE' button again for a few seconds.

4.2 Mode Selection and Usage

Press the 'MODE' button to cycle through the available modes. Use the 'LEVEL' button to adjust intensity (1-5 levels).

Your browser does not support the video tag.

Video 1: Demonstration of the multifunctional facial massager, showcasing various modes and application techniques.

1. Clean Mode (Green LED + Vibration):

Use with cotton pads and makeup remover or cleanser. This mode utilizes green light and acoustic vibration for deep cleansing, helping to remove dirt and impurities from the skin. Recommended duration: 5 minutes.

2. Lift Mode (Yellow LED + RF + EMS Micro-current):

Apply a facial mask. This mode combines yellow light, radio frequency, and EMS micro-current to promote skin lifting and tightening, enhancing the absorption of mask ingredients. Recommended duration: 10 minutes.

3. Lead-in Mode (Pink LED + RF + EMS Micro-current):

Use with toner, cream, lotion, or essence. Pink light, radio frequency, and EMS micro-current work together to improve the absorption of skincare products, aiding in moisturizing and locking in hydration. Recommended duration: 10 minutes.

4. Anti-aging Mode (Red LED + RF + Vibration):

This mode uses red light, radio frequency, and acoustic vibration to stimulate collagen regeneration, targeting anti-aging and anti-wrinkle benefits. Recommended duration: 10 minutes.

5. Remove Acne Mode (Blue LED):

Blue light therapy helps to shrink pores and maintain oil-water balance, assisting in the treatment of acne and breakouts. Recommended duration: 5 minutes.

4.3 Application Technique:

- Gently move the device over your face and neck, avoiding the eye area.
- For best results, move the probe head from the inside to outside, and from the bottom to top of your face.
- Ensure the device maintains full contact with your skin during use.

How to Use



Image 4: Visual guide on the proper usage of the device.

5. MAINTENANCE

5.1 Cleaning

- After each use, clean the device head with a damp cloth or a soft, dry cloth.
- Do not immerse the device in water or use abrasive cleaners.

5.2 Storage

- Store the device in a cool, dry place, away from direct sunlight and extreme temperatures.
- Keep out of reach of children.

6. TROUBLESHOOTING

If you encounter issues with your MiSMON Facial Machine, please refer to the following common solutions:

- **Device does not turn on:** Ensure the device is fully charged. Press and hold the 'MODE' button for a longer duration (approximately 3-5 seconds) to activate it.
- **No sensation or weak effect:** Verify that the device is in full contact with your skin. Ensure you have applied a water-based conductive gel, lotion, or essence, as this is crucial for the RF and EMS functions to work effectively. Increase the intensity level using the 'LEVEL' button.
- **Lights are not visible:** Some LED lights may be less visible in bright environments. Ensure you are in a moderately lit area. Confirm the device is powered on and in the correct mode.
- **Battery drains quickly:** Ensure the device is fully charged before each session. If the issue persists, contact customer support.

7. SPECIFICATIONS

- **Model Number:** OLB-BP-001a
- **Product Dimensions:** 7 x 1.96 x 1.96 inches
- **Weight:** 15.84 ounces

- **Battery:** 1 Lithium Ion battery (included)
- **Material:** Plastic
- **Manufacturer:** MiSMON

8. WARRANTY AND SUPPORT

For any product inquiries, technical support, or warranty information, please contact MiSMON customer service through the retailer where the product was purchased. Please note that statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.