

[Manuals.plus](#) /

› [Philips Sonicare](#) /

› Philips Sonicare One Brush Head Replacement Manual (BH1022/04)

Philips Sonicare BH1022/04

Philips Sonicare One Brush Head Replacement Manual

Model: BH1022/04

1. INTRODUCTION

This manual provides essential information for the Philips Sonicare One Brush Heads (BH1022/04). These replacement brush heads are specifically designed to maintain the performance and effectiveness of your Philips One electric toothbrush, ensuring a consistent and bright smile. For optimal oral hygiene, regular replacement of brush heads is crucial.



Image: Philips One by Sonicare Brush Heads (BH1022/04) in packaging, alongside two individual brush heads.

2. INSTALLATION AND SETUP

2.1. Compatibility

These brush heads are compatible exclusively with Philips One by Sonicare toothbrush handles. Using them with other Sonicare models or non-Philips One toothbrushes may result in an improper fit and reduced performance.

Compatible with
Philips One
by Sonicare
toothbrush
handles only



Image: Illustration showing three Philips One by Sonicare toothbrush handles (coral, white, and navy blue) with brush heads, emphasizing compatibility.

2.2. Attaching the Brush Head

1. Ensure your Philips One toothbrush handle is clean and dry.
2. Align the new brush head with the metal shaft on the toothbrush handle.
3. Gently push the brush head down onto the shaft until it snaps securely into place. Do not twist.
4. Verify that the brush head is firmly attached and there is no significant gap between the head and the handle.

Easily snap-on

to Philips One by
Sonicare toothbrushes



Image: A visual guide demonstrating the simple snap-on mechanism for attaching a Philips One by Sonicare brush head to its handle.

3. OPERATION

The Philips One by Sonicare brush heads feature micro-vibrations and tapered bristles made of soft nylon. These elements work together to gently polish teeth, contributing to a brighter smile. Once installed, the brush head operates in conjunction with your Philips One toothbrush handle's power settings.

Micro vibrations

and tapered bristles
polish teeth
for a brighter
smile



Image: A close-up of a person using a Philips One by Sonicare toothbrush, illustrating the gentle micro-vibrations and tapered bristles in action.

3.1. Product Overview Videos

Your browser does not support the video tag.

Video: An official Philips Sonicare overview of the Philips One Battery Toothbrush, showcasing its features and design.

Your browser does not support the video tag.

Video: An official Philips Sonicare video highlighting the Philips One Rechargeable Toothbrush, demonstrating its use and benefits.

4. MAINTENANCE AND REPLACEMENT

4.1. Daily Care

After each use, rinse the brush head thoroughly under running water to remove any toothpaste residue and debris. Allow it to air dry.

4.2. Replacement Schedule

For optimal performance and hygiene, dentists recommend replacing your brush head every 3 months. Worn bristles are less effective at cleaning and can harbor bacteria.

Dental Professionals recommend replacing your brush head **every three months for best results**



Image: A visual reminder to replace your toothbrush head every three months for best results, featuring a calendar icon.

4.3. Color Options

Philips One brush heads are available in various colors, allowing you to match your toothbrush handle or mix and match for a personalized look.

Brush heads match
colorful handles . . .
**or mix &
match
to create
a unique
design**



Image: An assortment of Philips One by Sonicare brush heads and handles in different colors, demonstrating the ability to match or mix colors.

5. TROUBLESHOOTING

If you encounter issues with your Philips One brush head, consider the following:

- **Improper Fit:** Ensure the brush head is securely snapped onto a Philips One by Sonicare handle. These brush heads are not designed for other Philips Sonicare models, which may lead to a loose fit or reduced vibration.
- **Reduced Effectiveness:** If bristles appear splayed or worn, it's time for a replacement. Old brush heads are less effective at cleaning.
- **No Vibration:** Confirm the brush head is correctly attached and the toothbrush handle is charged and functioning.

6. PRODUCT SPECIFICATIONS

Feature	Detail
Model Number	BH1022/04
Product Dimensions	0.79 x 2.36 x 7.87 inches
Item Weight	0.64 ounces
Material	Nylon bristles
Age Range	Adult
Bristle Type	Tapered, soft nylon
Cleaning Action	Micro-vibrations for gentle polishing
Quantity	Two (2) brush heads per pack

7. WARRANTY AND SUPPORT

For detailed warranty information, product registration, or further technical support regarding your Philips Sonicare One brush heads or toothbrush, please visit the official Philips Sonicare website or contact their customer service directly. Keep your purchase receipt for any warranty claims.

Official Philips Sonicare Store: [Visit Store](#)

