

[Manuals+](#)

[Q & A](#) | [Deep Search](#) | [Upload](#)

REVITIVE 5606-5574AQ-UK

Revitive Medic Knee Circulation Booster User Manual

Model: 5606-5574AQ-UK

INTRODUCTION

The Revitive Medic Knee Circulation Booster is designed to help relieve knee pain and discomfort. Utilizing patented OxyWave and IsoRocker technologies, it aims to boost blood flow and provide targeted relief. This manual provides essential information for the safe and effective use of your device.



Image: The Revitive Medic Knee Circulation Booster, shown with its packaging and electrode pads.

SAFETY INFORMATION

Before using the Revitive Medic Knee Circulation Booster, please read all safety warnings carefully. Do not use this device if you are:

- Fitted with an electronic implanted device such as a heart pacemaker or AICD.
- Being treated for, or have the symptoms of, an existing deep vein thrombosis (DVT).
- Pregnant.

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

SETUP

1. **Unpack:** Carefully remove all components from the packaging.
2. **Battery Installation:** Ensure 2 AA batteries (included) are correctly installed in the device.

3. **Placement:** Place the Revitive unit on the floor in front of your chair.
4. **Power On:** Press the power button on the device or remote control to turn it on.

OPERATING INSTRUCTIONS

Using the Foot Pads:

For general circulation boosting and leg pain relief, place both bare feet on the foot pads of the device. Adjust the intensity using the '+' and '-' buttons on the device or remote control until you feel a comfortable muscle contraction.



Image: A user's feet positioned on the Revitive Medic Knee Circulation Booster during operation.

Using the Electrode Pads for Knee Pain Management:

The Revitive Medic Knee comes with electrode pads for targeted pain relief and muscle strengthening programs for the knee. Connect the electrode pads to the device and apply them to the skin around the knee area as directed in the included quick start guide or as advised by a healthcare professional.

REVITIVE

OxyWave Technology

Patented stimulation for a
clinically proven therapy



Image: Electrode pads from the Revitive Medic Knee applied to a person's knee for targeted therapy.

Adjust the intensity for the electrode pads separately from the foot pads. The device features an intuitive display to show current intensity levels and remaining session time.



Image: A user operating the Revitive Medic Knee Circulation Booster with the included remote control.

MAINTENANCE

- **Cleaning:** Wipe the device with a soft, damp cloth. Do not use abrasive cleaners or immerse the unit in water.
- **Storage:** Store the device in a cool, dry place away from direct sunlight when not in use.
- **Electrode Pads:** Replace electrode pads when they lose their stickiness to ensure optimal conductivity.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on.	Low or dead batteries; incorrect battery installation.	Replace batteries or ensure they are installed correctly.

Problem	Possible Cause	Solution
No sensation during use.	Intensity too low; poor skin contact; dry skin.	Increase intensity gradually; ensure feet/pads are firmly on the device; moisturize skin if dry.
Electrode pads not sticking.	Pads are worn out or dirty.	Clean skin before application; replace electrode pads.

If problems persist, please contact customer support.

SPECIFICATIONS

Feature	Detail
Product Dimensions	14.17 x 3.74 x 15.35 inches
Item Model Number	5606-5574AQ-UK
Batteries	2 AA batteries required (included)
Item Weight	1.76 Kilograms (3.88 Pounds)
Power Source	Battery Powered
Use For	Knee
Massager Form	Massage Belt

WARRANTY & SUPPORT

For information regarding product warranty, please refer to the documentation included with your purchase or visit the official REVITIVE website. For technical support or further assistance, please contact REVITIVE customer service. You can also find more information and support at the [REVITIVE Store on Amazon](#).