

Zeerkeer 9A03088

Zeerkeer Electronic Automatic Pill Dispenser with Timer

Model: 9A03088

1. INTRODUCTION

This instruction manual provides detailed guidance for the proper use and maintenance of your Zeerkeer Electronic Automatic Pill Dispenser. This device is designed to assist with timely medication management through timed reminders and organized storage. Please read this manual thoroughly before operating the device to ensure safe and effective use.

2. PRODUCT FEATURES

- **Timed Reminders:** Features timed vibration and audible reminders for medication.
- **Digital Display:** Clear digital display shows time and timer settings.
- **Durable Construction:** Made from high-quality, food-grade ABS and PP materials, ensuring safety and durability.
- **Partitioned Storage:** Multiple compartments for easy classification and organization of different medications.
- **Moisture and Dust Protection:** Designed to protect medication from moisture, dust, and deterioration.
- **Compact and Portable:** Slim and lightweight design for easy carrying in a pocket or bag.
- **Clock and Timer Function:** Displays time in 12-hour or 24-hour format and includes a countdown timer function.

3. PACKAGE CONTENTS

Upon opening the package, verify that all items are present and undamaged:

- Zeerkeer Electronic Automatic Pill Dispenser (1 unit)

Note: Batteries may not be included and might need to be purchased separately.

4. PRODUCT OVERVIEW

Familiarize yourself with the main components of your pill dispenser:

Prevent pill deterioration

Compact & seal design



① Seamless Compartments



② Strong hinge



③ Compact buckle

Image: This image displays the Zeerkeer Electronic Automatic Pill Dispenser. On the left, the main unit with a digital display showing time and control buttons (HOUR, MODE, MIN). On the right, the pill dispenser is open, revealing four internal compartments, each containing different colored pills. Below, three detailed close-up images highlight key features: (1) "Seamless Compartments" showing the smooth, easy-to-clean interior, (2) "Strong Hinge" illustrating the durable connection between the lid and base, and (3) "Compact Buckle" demonstrating the secure closing mechanism.

The device consists of a main unit with a digital display and control buttons, and a hinged lid that opens to reveal four internal compartments for medication storage. It features a secure buckle and a robust hinge for durability.

5. SETUP

5.1. Battery Installation

The pill dispenser requires batteries for operation. Locate the battery compartment, usually on the back or bottom of the device. Insert the required batteries, ensuring correct polarity (+/-). Close the battery compartment securely.

5.2. Initial Time Setting

After battery installation, the display may show a default time or be blank. To set the current time:

1. Press the **MODE** button to enter time setting mode. The hour digits will likely flash.
2. Use the **HOUR** button to adjust the hour.
3. Use the **MIN** button to adjust the minutes.
4. Press **MODE** again to confirm the time and exit time setting mode.
5. Select 12-hour or 24-hour format if available, typically by pressing **MODE** or **HOUR** while in time display mode.

6. OPERATING INSTRUCTIONS

6.1. Setting Medication Reminders (Alarms)

The device supports up to 5 daily medication reminders. To set an alarm:

1. From the main time display, press the **MODE** button repeatedly until "ALARM 1" (or similar indicator) appears.
2. The hour digits for the alarm will flash. Use the **HOUR** button to set the desired alarm hour.
3. Use the **MIN** button to set the desired alarm minutes.
4. Press **MODE** to save "ALARM 1" and proceed to "ALARM 2" if needed, or exit alarm setting mode.
5. Repeat for up to 5 alarms.
6. To activate or deactivate an alarm, there may be a dedicated switch or a combination of button presses (refer to on-screen indicators).

6.2. Using the Countdown Timer

The device includes a countdown timer function (up to 23 hours 59 seconds).

1. From the main time display, press **MODE** until the countdown timer display appears (e.g., "00:00").
2. Use the **HOUR** and **MIN** buttons to set the desired countdown duration.
3. Press **MODE** or a dedicated START/STOP button to begin the countdown.
4. When the countdown reaches zero, an alarm will sound or vibrate.

6.3. Filling the Compartments

Open the lid of the pill dispenser. Carefully place your medications into the four separate compartments. Ensure the lid closes securely to protect the contents.



Image: This image shows the Zeerker Electronic Automatic Pill Dispenser with its lid open, revealing four distinct compartments. Each compartment is filled with different types and colors of pills, demonstrating its organizational capacity. Above the dispenser, two icons are visible: one representing "Vibration" (a battery with vibration waves) and another representing "Voice" (an alarm clock), indicating the reminder features of the device.

7. MAINTENANCE

7.1. Cleaning

To clean your pill dispenser:

1. Ensure all medication is removed from the compartments.
2. Wipe the exterior and interior surfaces with a soft, damp cloth.
3. For stubborn spots, use a mild soap solution, then wipe clean with a damp cloth and dry thoroughly.
4. Do not immerse the device in water or use harsh chemical cleaners, as this can damage the electronic components.

7.2. Battery Replacement

When the display dims or the alarms become weak, it is time to replace the batteries. Follow the battery installation steps (Section 5.1) to replace old batteries with new ones.

8. TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not power on or display is blank.	No batteries, incorrect battery installation, or dead batteries.	Check battery installation, ensure correct polarity, or replace batteries.
Alarms are not sounding/vibrating.	Alarms not set, alarms deactivated, or low battery.	Verify alarm settings, ensure alarms are activated, or replace batteries.
Time display is incorrect.	Time not set correctly or battery removed/replaced.	Reset the current time following Section 5.2.
Lid does not close securely.	Obstruction in hinge or buckle area, or overfilled compartments.	Check for debris, ensure compartments are not overfilled, and inspect hinge/buckle for damage.

9. SPECIFICATIONS

Model: 9A03088

Brand: Zeerkeer

Material: ABS + PP (Acrylonitrile Butadiene Styrene + Polypropylene)

Color: Blue

Dimensions: 6.5 cm (L) x 6.5 cm (W) x 2 cm (H)

Weight: Approximately 60 grams

Functions: 12/24 hour time display, 5-time 24-hour timing vibration reminder, 23-hour 59-second countdown reminder.

UPC: 699900277805



Image: This image provides a top-down view of the Zeerkeer Electronic Automatic Pill Dispenser, clearly showing its digital display and control buttons. Overlaid on the image are precise dimensions: 6.5 cm (2.55 inches) for both length and width, and 2 cm (0.78 inches) for height, illustrating the compact size of the device.

10. SAFETY INFORMATION

- Keep the device out of reach of children and pets.
- Do not expose the device to extreme temperatures, direct sunlight, or moisture.
- Do not attempt to disassemble or repair the device yourself. Contact customer support if issues arise.
- Dispose of old batteries responsibly according to local regulations.
- This device is intended for organizing and reminding about medication. It is not a substitute for professional medical advice. Always consult with a healthcare professional regarding your medication regimen.

11. WARRANTY AND SUPPORT

For warranty information or technical support, please refer to the contact details provided with your purchase or visit the official Zeerkeer website. Please have your model number (9A03088) and purchase date available when contacting support.