

[manuals.plus](#) /

› [Comfytemp](#) /

› [Comfytemp Weighted Heating Pad for Back & Cramps, Father's Day Gifts, Birthday Gifts for Mom Dad Women Men, 2.2lb Electric Neck Shoulder Heat Pad for Pain Relief, FSA HSA Eligible, 2H Auto-Off, 12x24 Navy Blue 12x24 Inch](#)

Comfytemp K9040-BT

Comfytemp Weighted Heating Pad User Manual

Model: K9040-BT | Brand: Comfytemp

1. PRODUCT OVERVIEW

The Comfytemp Weighted Heating Pad is designed to provide effective pain relief and comfort. Its unique weighted design ensures a snug fit, enhancing heat transfer for various body parts.

- **Weighted Design:** 2.2lb weighted element for a conformable fit and better heat transfer, ideal for back pain, cramps, muscle tension, and soreness.
- **Adjustable Heat Levels:** Features 3 heat levels (Low, Mid, High) for personalized therapy.
- **2-Hour Auto-Off:** Equipped with an automatic shut-off function after 2 hours for energy conservation and safety, preventing overheating.
- **Soft & Machine Washable:** Made with premium flannel for maximum comfort and easy cleaning. The controller is detachable for washing.
- **Moist Heat Therapy:** Can be used with moist heat by lightly spraying water on the surface for deeper penetration.



Image: The Comfytemp Weighted Heating Pad in Navy Blue, showing its plush surface and the attached controller.

2. SAFETY INFORMATION

Please read all safety instructions carefully before use to prevent injury or damage.

- Do not use while lying on the heating pad to avoid overheating and potential burns.
- Do not fold the heating pad during use. Ensure it is laid flat to allow for proper heat dissipation.
- Always unplug the heating pad and allow it to cool down completely before storing.
- Do not use on infants, persons insensitive to heat, or incapacitated persons.
- Avoid using pins or other metallic objects with the heating pad.
- Do not use if the heating pad appears damaged in any way.
- Consult a physician before use if you have any medical conditions.



Do Not Use While Lying On It !



Do Not Fold During Use!



Lay The Heating Pad Cool Down Before Storing

Image: Visual guide illustrating important safety warnings: Do not use while lying on it, Do not fold during use, and Lay the heating pad cool down before storing.

3. SETUP & OPERATION

3.1 Setup

1. Remove the heating pad from its packaging.
2. Connect the detachable controller to the heating pad. Ensure it clicks securely into place.
3. Plug the heating pad into a standard power outlet.

Your browser does not support the video tag.

Video: Official Comfytemp video demonstrating the setup and basic operation of the weighted heating pad, including connecting the controller and power-on sequence.

3.2 Operation

1. Press the power button on the controller to turn on the heating pad.
2. Select your desired heat level (Low, Mid, High) by pressing the temperature button. The indicator lights will

show the current setting.

3. Place the heating pad on the area of your body that requires heat therapy. Ensure it is laid flat and not folded.
4. The heating pad will automatically shut off after 2 hours. To manually turn off, press and hold the power button for 3 seconds.
5. For moist heat therapy, lightly spray water on the surface of the heating pad before use.



Image: Close-up of the heating pad's controller, highlighting the 2-hour auto-off timer, 3 heat levels, and overheating protection features.

4. PRODUCT FEATURES & BENEFITS

The Comfytemp Weighted Heating Pad offers several features designed for optimal comfort and pain relief:

- **Targeted Pain Relief:** Effectively relieves back pain, period cramps, muscle tension, and soreness.
- **Weighted Comfort:** The 2.2lb weight helps the pad conform to your body, providing gentle pressure and enhancing the therapeutic effect, similar to a warm hug.
- **Versatile Application:** Suitable for use on the back, neck, shoulders, abdomen, knees, and other areas

experiencing discomfort.

- **Premium Flannel Material:** Super soft to the touch, providing maximum comfort during use.
- **FSA/HSA Eligible:** Approved for Flexible Spending Account (FSA) and Health Savings Account (HSA) use.

Help Soothe Menstrual Cramps



WEIGHTED

- ✓ XL wrap the abdomen
- ✓ Fit better, more heat
- ✓ Gentle pressure
- ✓ Relax muscles
- ✓ Like a warm hug



Effective for soothing menstrual cramps.

Weighted Makes It All Difference

Better cover your body with additional glass beads



WEIGHTED



Weighted design ensures better body coverage and heat transfer.

Targeted Pain Relief for Muscle Soreness

Inner beads ensure it keeps in place to the body rather than stiffly away



WEIGHTED



Targeted pain relief for muscle soreness in the shoulder area.



Versatile use for back, neck, shoulders, knees, and feet.

5. CARE & MAINTENANCE

Proper care ensures the longevity and effectiveness of your heating pad.

- **Machine Washable:** Detach the controller before washing. The heating pad can be washed by hand or machine.
- **Drying:** Do not tumble dry. Air dry completely before reattaching the controller or storing.
- **Storage:** Ensure the heating pad is unplugged and cool before folding and storing in a dry place.

Safe and Trusted



Image: Diagram showing the internal layers of the heating pad, including premium flannel, glass beads, heating wires, and protective polyester.

6. SPECIFICATIONS

Feature	Detail
Model Number	K9040-BT
Product Dimensions	24 x 12 x 0.39 inches
Item Weight	2.2 Pounds
Material	Flannel
Color	Navy Blue
Heat Levels	3 (Low, Medium, High)
Auto-Off Timer	2 Hours

Manufacturer	Comfytemp
Country of Origin	China

7. TROUBLESHOOTING

If you encounter any issues with your Comfytemp Weighted Heating Pad, please refer to the following common troubleshooting steps:

- **No Heat:** Ensure the power cord is securely plugged into both the heating pad and a working electrical outlet. Check if the power button on the controller is pressed and the indicator light is on.
- **Pad Not Heating Evenly:** Ensure the heating pad is laid flat and not folded or bunched up during use. The weighted beads should be evenly distributed.
- **Automatic Shut-off Before 2 Hours:** Verify that the heating pad is not folded or covered in a way that could trigger overheating protection. Unplug and replug to reset.
- **Controller Feels Warm:** It is normal for the controller to become slightly warm during use. If it becomes excessively hot or emits a burning smell, unplug immediately and discontinue use.

If the problem persists after attempting these steps, please contact Comfytemp customer support for further assistance.

8. WARRANTY & SUPPORT

Comfytemp products are manufactured with quality and care. For any inquiries or support regarding your weighted heating pad, please refer to the contact information provided with your product packaging or visit the official Comfytemp website.

Please retain your proof of purchase for warranty claims. Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

For the most up-to-date information and support, please visit the [Comfytemp Store on Amazon](#).