

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

[manuals.plus](#) /

› [veicomtech](#) /

› [veicomtech Ionic Foot Bath Detox Machine User Manual](#)

## veicomtech Ionic Foot Bath Detox Machine

# veicomtech Ionic Foot Bath Detox Machine User Manual

Model: Ionic Foot Bath Detox Machine

## INTRODUCTION

The veicomtech Ionic Foot Bath Detox Machine is designed to provide a convenient and effective foot spa experience. This system utilizes a mild electric current to electrolyze minerals in water, generating ions that interact with the water to facilitate a detoxifying process. Regular use may contribute to enhancing your immune system, balancing blood sugar and pressure, improving memory and sleep, and providing significant pain relief.

This manual provides comprehensive instructions for the safe and proper setup, operation, and maintenance of your ionic foot bath machine.

## IMPORTANT SAFETY INFORMATION

**Please read all instructions carefully before using the device. Failure to do so may result in injury or damage to the product.**

- Do not immerse the main control unit in water.
- Keep out of reach of children.
- Consult a healthcare professional before use if you have any medical conditions, are pregnant, or have electronic implants (e.g., pacemakers).
- Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.
- Ensure the power supply is disconnected before cleaning or performing maintenance.

## PACKAGE CONTENTS

Your veicomtech Ionic Foot Bath Detox Machine package includes:

- 1 x Main Control Unit
- 1 x Foot Basin
- 2 x Arrays
- 100 x Tub Liners
- 1 x Power Adapter



Image: The complete package contents of the veicomtech Ionic Foot Bath Detox Machine.

## SETUP INSTRUCTIONS

1. Place the foot basin on a flat, stable surface.
2. Line the foot basin with one of the provided tub liners. This helps maintain hygiene and simplifies cleanup.
3. Fill the basin with warm water to a comfortable level, ensuring it covers the array once placed.
4. Connect the power adapter to the DC input port on the main control unit.
5. Connect one or both arrays to the 'ARRAY' ports on the main control unit using the 6.35mm (1/4 inch) adapters. Ensure a secure connection.
6. Place the array(s) into the water-filled basin.
7. Position the main control unit securely on the designated slot of the foot basin.

## 4 STEPS TO ENJOY FOOT SPA



Image: Visual guide illustrating the four steps for setting up your foot spa.

### Setup Video Guide

Your browser does not support the video tag.

Video: This official product video demonstrates the correct setup procedure for the veicomtech Ionic Foot Bath Detox Machine, including connecting the components and preparing the foot bath.

### OPERATING INSTRUCTIONS

1. Once the unit is set up and plugged in, the 6 LEDs on the control unit will flash, indicating standby status.
2. Place your feet into the water.
3. Add a small amount of salt to the water. This enhances the conductivity for the ionic process.
4. Press the On/Off button on the control unit to start a session. The machine is programmed for a 30-minute session.

# ONE TOUCH TO START



Image: The control unit, highlighting the simple one-touch operation to begin a session.

# 30 MINS TIMER

Each light counts 5 mins



Image: Illustration of the 30-minute timer, where each light on the control unit signifies 5 minutes of the session.

## Understanding Water Color Changes

It is important to note that the water in the foot bath will change color during a session. This is a natural result of the electrolysis process, not solely due to toxins being drawn from the body. The device uses a mild electric current to electrolyze minerals present in the water (such as salts), generating oxidized substances (like iron hydroxide, chlorine, etc.). These substances cause the water to change color. Therefore, it is normal for the water to change color even without feet in the basin.

## What is an Ionic Foot Spa?

An ionic foot spa is a machine used for body dextoxifcation through the feet by electrolytic reaction producing ions that attract and neutralize toxins in the body in a tub of salt water.



## How does it work ?

### Water Color Changes

The device utilizes the principle of electrolysis (chemical reaction), a process that breaks down water into oxygen and hydrogen. These gases are attracted by electrodes - metal array parts leading to an oxidation reaction. The final color change is the result of combination of toxins eliminated & oxidation reaction.

So it is normal that the water color may change even without feet in.

Image: Explaining the science behind water color changes during an ionic foot spa session.

# PROCESS OF IONIC DETOXIFICATION



Image: A visual representation of the water changing color during the ionic detoxification process.

## MAINTENANCE AND CARE

- After each use, dispose of the used tub liner and thoroughly clean the foot basin with mild soap and water.
- Clean the array(s) after each session. Rinse them with clean water and gently scrub any residue. Ensure they are completely dry before storage.
- Store the main control unit and arrays in a dry, cool place away from direct sunlight.
- Replace arrays as needed. The lifespan of an array depends on usage and water quality.

## TROUBLESHOOTING

### Device Not Powering On or Functioning

- **Check Power Connection:** Ensure the power adapter is securely plugged into both the wall outlet and the DC input port on the main control unit.

- **Verify Array Connection:** Make sure the array(s) are correctly connected to the 'ARRAY' ports using the proper 6.35mm (1/4 inch) adapter. Using an incorrect adapter (e.g., a 3-pin adapter) will prevent the unit from functioning. Refer to the setup video for correct adapter usage.

## Water Not Changing Color

- **Add Salt:** Ensure you have added a sufficient amount of salt to the water. Salt is crucial for the electrolysis process.
- **Check Array Condition:** Over time, arrays may wear out. If the array is old or damaged, it may not effectively facilitate the ionic process. Replace the array if necessary.
- **Water Mineral Content:** The degree of color change can vary based on the mineral content of your local water.

## SPECIFICATIONS

Feature	Detail
Product Dimensions	17.52 x 13.07 x 7.09 inches
Item Weight	6.19 Pounds
Model Number	Ionic Foot Bath Detox Machine
Manufacturer	veicomtech
Power Source	Corded Electric
Use For	Feet
Color	Blue Panel

## WARRANTY AND SUPPORT

For warranty information and customer support, please refer to the contact details provided with your product packaging or visit the official veicomtech website. Please retain your proof of purchase for warranty claims.