

[Manuals.plus](#) /

> [Braun](#) /

> Braun Oral-B Genius X D706 Electric Toothbrush Replacement Power Handle Instruction Manual

## Braun D706

# Braun Oral-B Genius X D706 Electric Toothbrush Replacement Power Handle Instruction Manual

Model: D706 | Type: 3771

This manual provides essential information for the safe and effective use of your Braun Oral-B Genius X D706 Replacement Power Handle.

## 1. INTRODUCTION

This document serves as an instruction manual for the Braun Oral-B Genius X D706 Replacement Power Handle. It is designed to guide users through the setup, operation, maintenance, and troubleshooting of the device. Please read this manual thoroughly before first use and retain it for future reference.

## 2. PRODUCT OVERVIEW

The Braun Oral-B Genius X D706 is a replacement power handle for an electric toothbrush system. It features 6 brushing modes and is powered by a Li-Ion battery. This handle is designed to be used with existing Oral-B chargers and brush heads.



manual for specific mode descriptions.

### 3. **Brushing Technique:**

- Wet the brush head and apply a small amount of toothpaste.
- Place the brush head on your teeth at a slight angle towards the gumline.
- Guide the brush head slowly from tooth to tooth, allowing it to clean each tooth surface for a few seconds.
- Do not scrub. Let the toothbrush do the work.
- Brush all four quadrants of your mouth evenly.

4. **Smart Features:** If your original Genius X system included a smartphone app or pressure sensor, this replacement handle is designed to integrate with those features.

## 5. MAINTENANCE

---

- **After Each Use:** Rinse the brush head thoroughly under running water with the handle switched on. Remove the brush head and clean both the handle and the brush head separately.
- **Cleaning the Handle:** Wipe the handle with a damp cloth. Do not immerse the handle in water.
- **Brush Head Replacement:** Dentists recommend replacing your brush head every 3 months, or sooner if the bristles are splayed or worn.
- **Battery Care:** To maintain the rechargeable battery's maximum capacity, unplug the charging unit and fully discharge the handle at least once every 6 months. Then, recharge it fully.

## 6. TROUBLESHOOTING

---

Problem	Possible Cause	Solution
Toothbrush does not turn on.	Battery is discharged.	Place the handle on the charger for at least 30 minutes. Ensure the charger is plugged in.
Toothbrush does not charge.	Charger not connected or faulty. Handle not seated correctly.	Check charger connection. Ensure handle is firmly placed on the charger. Try a different outlet.
Reduced brushing power.	Low battery. Worn brush head.	Recharge the handle. Replace the brush head.
Excessive vibration or noise.	Brush head not properly attached.	Remove and reattach the brush head, ensuring it clicks into place.

If the problem persists, please contact customer support.

## 7. SPECIFICATIONS

---

**Model:** D706

**Type:** 3771

**Power Source:** Rechargeable Li-Ion Battery

**Brushing Modes:** 6

**Color:** White

**Manufacturer:** Braun

**Package Dimensions:** 8.27 x 1.18 x 1.18 inches

**Weight:** 8 ounces

## 8. WARRANTY INFORMATION

---

This replacement power handle is typically covered by a manufacturer's warranty. Please refer to the warranty documentation included with your original Oral-B Genius X product or visit the official Braun Oral-B website for detailed warranty terms and conditions. The standard warranty period for Oral-B power handles is generally two years from the date of purchase.

*Note: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This disclaimer is provided as per legal requirements and does not directly pertain to the operation of the toothbrush handle.*

## 9. SUPPORT

---

For further assistance, product information, or to purchase replacement parts and accessories, please visit the official Braun Oral-B website or contact their customer service department.

- **Online Support:** [www.oralb.com](http://www.oralb.com)
- **Manufacturer:** Braun