



Manuals.plus /

› InnovaGoods /

› InnovaGoods Anti-Aging Eye Massager V0103319 User Manual

## InnovaGoods V0103319

# InnovaGoods Anti-Aging Eye Massager V0103319 User Manual

Model: V0103319

## 1. INTRODUCTION

Thank you for choosing the InnovaGoods Anti-Aging Eye Massager. This device is designed to provide a comprehensive eye care experience through phototherapy and vibration, aiming to improve the appearance of the skin around the eyes. Please read this manual carefully before use to ensure proper operation and maintenance.

## 2. SAFETY INFORMATION

- Do not use if you have any eye conditions, skin sensitivities, or are pregnant. Consult a doctor if unsure.
- Keep out of reach of children.
- Do not immerse the device in water or other liquids.
- Do not use if the device is damaged.
- Use only the provided charging cable.
- Discontinue use if you experience any discomfort or adverse reactions.

## 3. PRODUCT OVERVIEW

The InnovaGoods Eye Massager combines gentle vibration with phototherapy to target the delicate skin around the eyes. It features a compact, portable design for convenient use.



Image: InnovaGoods Eye Massager with its protective cap and two USB charging cables.

- **Massage Head:** Applies phototherapy and vibration.
- **Power Button:** Turns the device on/off and cycles through modes.
- **Charging Port:** For connecting the USB charging cable.
- **Protective Cap:** Covers the massage head when not in use.

## 4. SETUP

### 4.1. Charging the Device

1. Ensure the device is turned off before charging.
2. Connect the small end of the provided USB cable to the charging port on the massager.

3. Connect the larger end of the USB cable to a USB power source (e.g., computer, wall adapter).
4. The indicator light will illuminate during charging and change color or turn off when fully charged.



Image: The eye massager connected to a laptop via its USB charging cable, indicating the charging process.

## 5. OPERATING INSTRUCTIONS

### 5.1. Preparing for Use

- Ensure your face and the area around your eyes are clean and dry.
- For enhanced results, apply a small amount of eye cream or serum to the area you wish to treat.



Image: A woman has applied small dots of eye cream around her eye, preparing for the massager's use.

## 5.2. Powering On and Mode Selection

1. Press and hold the power button for a few seconds to turn on the device. It will typically start in a default mode.
2. Press the power button repeatedly to cycle through the available modes:
  - **Red Light Phototherapy + Vibration:** Designed for anti-aging and collagen stimulation.
  - **Blue/Purple Light Phototherapy + Vibration:** May help with soothing and improving skin tone.
  - **Vibration Only:** Gentle massage without light therapy.



Image: A woman is using the eye massager with red light therapy around her eye area.



Image: A close-up view of the eye massager emitting a purple light during use on the skin.

### **5.3. Application**

- Gently place the massage head on the skin around your eyes.
- Move the device in slow, uniform, circular motions or gentle strokes around the area to be treated.
- Avoid direct contact with the eyeball.
- Recommended usage time is 5-10 minutes per session, 1-2 times daily.



Image: A woman gently massaging the skin under her eye with the device, after applying cream.

#### **5.4. Powering Off**

Press and hold the power button for a few seconds to turn off the device after use.

### **6. MAINTENANCE**

---

#### **6.1. Cleaning**

- After each use, gently wipe the massage head with a soft, damp cloth.
- Do not use abrasive cleaners or immerse the device in water.
- Ensure the device is completely dry before storing.

#### **6.2. Storage**

- Store the device in a cool, dry place, away from direct sunlight and extreme temperatures.
- Keep the protective cap on the massage head when not in use.



Image: The eye massager standing upright on a wooden surface next to a folded towel and decorative flowers, with its protective cap removed.

## 7. TROUBLESHOOTING

| Problem                            | Possible Cause                                  | Solution   |
|------------------------------------|---|--|
| Device does not turn on.           | Low battery.<br>Device not properly charged.    | Charge the device fully. Ensure the charging cable is securely connected.                      |
| Vibration is weak or inconsistent. | Low battery.<br>Device malfunction.             | Charge the device fully. If the problem persists, contact customer support.                    |
| Light therapy not working.         | Incorrect mode selected.<br>Device malfunction. | Cycle through modes using the power button. If the problem persists, contact customer support. |

## 8. SPECIFICATIONS

---

- **Model Number:** V0103319
- **Brand:** InnovaGoods
- **Dimensions:** 2.3 x 2.3 x 11.8 cm
- **Weight:** 90 g
- **Material:** Plastic
- **Features:** Anti-aging, Anti-dark circles treatment, Illuminator, Cordless
- **Power Source:** Rechargeable (USB)

## 9. WARRANTY AND SUPPORT

---

This product comes with a standard manufacturer's warranty. For detailed warranty information, please refer to the warranty card included with your purchase or visit the official InnovaGoods website.

For technical support, troubleshooting assistance, or any questions regarding your InnovaGoods Anti-Aging Eye Massager, please contact InnovaGoods customer service through their official channels.

**Website:** [www.innovagoods.com](http://www.innovagoods.com) *(Example link, please verify official website)*