

Zeerkeer 3-in-1 Acupuncture Pen

Zeerkeer 3-in-1 Electronic Acupuncture Massage Pen User Manual

Model: 3-in-1 Acupuncture Pen (159ce221-ad30-49c4-9c73-714024c5d22f)

1. INTRODUCTION

The Zeerkeer 3-in-1 Electronic Acupuncture Massage Pen is designed to provide a needle-free acupuncture experience for pain relief and meridian energy stimulation. This device combines modern technology with traditional principles to offer a convenient and effective solution for personal wellness. Please read this manual thoroughly before use to ensure proper operation and safety.

2. SAFETY INFORMATION AND IMPORTANT WARNINGS

- **Individual Sensitivity:** Everyone experiences stimulation differently. Always begin at intensity level 1 and gradually increase to a comfortable level. Higher intensity does not necessarily equate to better results.
- **Usage Frequency:** Use the device 1-2 times per day, for 10-15 minutes per session. A typical course of treatment is 10 days.
- **Electrical Sensation:** The acupuncture pen operates on the principle of single-phase electricity. Therefore, sensations like "Electrical Stimulation" or "Acupuncture Sensation" are normal and indicate the device is functioning as intended. Do not be alarmed by these sensations.
- **Static Electricity:** If your body's static electricity is high, you may experience a stronger shock sensation. Before using the product, ensure your hands are dry and touch a grounded object to discharge any static electricity.
- **Initial Charge:** Fully charge the device before its first use.
- **Disclaimer:** Statements regarding dietary supplements and health devices have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. Consult a healthcare professional for medical advice.

3. PACKAGE CONTENTS

Please verify that all items are present in your package:

- 1 x Zeerkeer Electronic Acupuncture Pen

- 3 x Interchangeable Massage Heads (Node, Spheroidal, Scraping)
- 1 x USB Charging Cable

4. PRODUCT FEATURES

- **Versatile Massage Heads:** Includes three metal massage heads for various applications: Node point for partial body pains, Spheroidal for quick pain alleviation, and Scraping for acupuncture and blood circulation.
- **Adjustable Intensity:** Features 9 adjustable intensity levels to suit individual comfort and needs.
- **Convenient Charging:** USB rechargeable with an LED display for intensity indication.
- **Energy Saving:** Equipped with an auto turn-off feature to prevent overheating and conserve battery life.
- **Pain Relief:** Designed to reduce chronic pain, improve mood, promote relaxation, enhance skin condition, and support blood circulation.
- **Wide Application:** Suitable for relieving pain in areas such as the shoulder, knee, nerves, ankle, and muscles without the use of needles.
- **User-Friendly:** Features auto-detection of meridian energy points, adjustable sensitivity, intensity, and frequency controls for personalized treatment.

5. SETUP

1. **Charge the Device:** Before first use, connect the USB charging cable to the pen and a suitable USB power source. The LED display will indicate charging status. A full charge typically takes 2-3 hours.
2. **Select Massage Head:** Choose one of the three interchangeable massage heads (Node, Spheroidal, or Scraping) and securely attach it to the pen.
3. **Prepare for Use:** Ensure your hands are dry. If you are prone to static electricity, touch a grounded object before handling the pen to minimize potential shock sensations.

9 gear free adjustment

For more personal and efficient use



2-3 gear is recommended when first using Comfortable button



3s no-operation at 0 gear for power off

Image: The Zeerkeer Acupuncture Pen shown with its dimensions and the included USB charging cable.

6. OPERATING INSTRUCTIONS

1. **Power On:** Press the On/Off button to turn on the device. The LED display will light up, showing the current intensity level.
2. **Adjust Intensity:** Use the "+" and "-" buttons to adjust the intensity level from 1 to 9. Always start at level 1 and gradually increase to a comfortable setting.
3. **Apply to Skin:** Gently place the chosen massage head onto the desired area of your skin. The pen is designed to detect meridian energy points automatically.
4. **Laser Function (if applicable):** To activate the laser function, press and hold the "+" button for three seconds.
5. **Massage Technique:** Move the pen slowly over the skin, focusing on areas of discomfort or known acupuncture points.
6. **Session Duration:** Use for 10-15 minutes per session, 1-2 times daily.
7. **Power Off:** Press the On/Off button to turn off the device. The pen also features an auto turn-off if no

operation is detected at 0 gear for 3 seconds.



Image: A detailed diagram labeling the strength gear display, charging port, pulse indicator, On/Off button, intensity controls, and replaceable massage tip.



Image: A hand holding the acupuncture pen, demonstrating the activation of the laser function by long-pressing the '+' button for three seconds, with a red laser dot visible on the palm.

2 massage tips

Different tip aims at different body part massage

Oblate
Contact with large body area,
softer when used



Ball type
Stronger when used, good results



Image: Illustration of the two main massage tips: the Oblate tip for large body areas and softer use, and the Ball type for stronger, more focused results.



Image: A diagram illustrating various acupuncture points on the human body, such as Shaofu Point, Neiguan, Waiguan, Knee acupoint, Life gate, Dazhui, Zhiyang, and Zusanli, indicating areas for pen application.

7. MAINTENANCE

- **Cleaning:** After each use, wipe the massage heads and the body of the pen with a clean, dry cloth. Do not use abrasive cleaners or immerse the device in water. The device is not water resistant.
- **Storage:** Store the pen in a cool, dry place, away from direct sunlight and out of reach of children.
- **Battery Care:** To prolong battery life, avoid fully discharging the battery frequently. Recharge the device when the battery indicator is low.

8. TROUBLESHOOTING

- **No Power:** Ensure the device is fully charged. If the device does not turn on after charging, contact customer support.
- **Strong Shock Sensation:** This can be due to high static electricity in the body. Before use, ensure your hands are dry and touch a grounded object. Also, ensure the massage head is in full contact with the skin.
- **Discomfort at Higher Intensities:** If using the device above medium power causes discomfort or a sensation of discharge into your arm, reduce the intensity level. It is recommended to stay within comfortable levels.
- **Device Not Working:** Check if the massage head is securely attached. Ensure the device is charged. If issues persist, refer to the warranty and support section.

9. SPECIFICATIONS

Material	ABS / Plastic
Power	0.6W

Pulse Frequency	6HZ
Pulse Width	160-180US
Pulse Peak	0-500 V
Working Current (with laser)	Less than 150MA
Battery Capacity	450 mAh (1 Lithium Polymer battery included)
Charging Time	2-3 hours
Continuous Use Time	3-4 hours
Item Dimensions (L x W x H)	1" x 1" x 8" (approx. 3.4cm x 3.4cm x 17.9cm)
Power Source	Battery Powered
Special Features	Portable, Rechargeable, Cordless
Water Resistance Level	Not Water Resistant
Manufacturer	Zeerkeer

10. WARRANTY AND SUPPORT

This product does not have explicit warranty information provided in the product details. For any warranty claims, technical support, or operational questions, please contact the seller or manufacturer directly through your purchase platform. Keep your purchase receipt or order number handy for faster service.