

[Manuals.plus](#) /

- › [Amazon Basic Care](#) /
- › Amazon Basic Care Temple Touch Digital Thermometer Instruction Manual

Amazon Basic Care Temple Touch Digital Thermometer

Amazon Basic Care Temple Touch Digital Thermometer Instruction Manual

PRODUCT OVERVIEW

The Amazon Basic Care Temple Touch Digital Thermometer is designed for quick and accurate temperature readings. It offers a non-invasive method suitable for all ages, from newborns to adults, providing results in under 10 seconds. The device supports dual measurement options in both Fahrenheit and Celsius.

Key Features

- **FAST RESULTS:** Get temperature readings in under 10 seconds with this instant thermometer.
- **NON-INVASIVE MEASUREMENT:** Use this temporal thermometer comfortably on the temple, minimal disruption whether awake or asleep, perfect for babies to adults and even sensitive users.
- **UNIVERSAL READING:** This thermometer offers versatile temperature readings in both Fahrenheit and Celsius for global compatibility.
- **CONVENIENCE ASSURED:** This digital thermometer for kids and adults arrives ready to use with battery included, saving time and ensuring immediate functionality, ideal for busy parents and caregivers.
- **FOR ACCURATE RESULTS:** Hold to temple, wait for beep & read result.



amazon

basic+care

Non-invasive

Meets accuracy
standards for
clinical
thermometers

Temple Touch

Digital
Thermometer

Results in under
10 seconds

- Temple measurement -
awake or asleep
- Reads in °F & °C



1 Digital Thermometer

Image: Amazon Basic Care Temple Touch Digital Thermometer packaging and the thermometer itself.

GETTING STARTED

Battery Installation

The thermometer requires 2 AAA batteries (included). Ensure correct polarity when inserting the batteries into the compartment located at the base of the device.



Image: The thermometer with its battery compartment open, illustrating where to insert the two AAA batteries.

Initial Setup

Before first use, allow the thermometer to adjust to room temperature for at least 30 minutes. The ambient temperature should be between 10°C-40°C (50°F-104°F) for optimal performance. It is not recommended for outdoor use in extreme cold or hot conditions.

HOW TO USE YOUR THERMOMETER

Taking a Temperature Reading

1. Press the power button once to turn on the device.
2. Place the probe tip gently on the temple area.

3. Press and release the 'Start' button. Wait for a long beep, then read the result displayed on the screen.
4. Wait at least 1 minute between consecutive measurements for accuracy.

R.A.T.E.™ Technology

Rapid Accurate Temperature Establishment—utilizing a unique sensor, R.A.T.E.™ rapidly samples the heat flowing from the blood vessels to the skin surface and accurately translates it to the core body temperature.

Image: A diagram illustrating the correct placement of the thermometer on the temple for accurate readings, and an incorrect placement example.

Changing Temperature Units (°F/°C)

To switch between Fahrenheit (°F) and Celsius (°C), ensure the device is off. Press and hold the power button for 5 seconds until the unit changes.

Switching Measurement Modes (Forehead/Object)

The thermometer supports both Forehead and Object modes. To switch between modes, press the mode button (often indicated by a person/house icon) while the device is on. The display will show the active mode.

Your browser does not support the video tag.

Video: A quick demonstration of the Amazon Basic Care Temple Touch Digital Thermometer in use, highlighting its ease and speed.

Your browser does not support the video tag.

Video: See the thermometer provide a temperature reading in just 6 seconds, demonstrating its rapid measurement

capability.

CARE AND MAINTENANCE

Cleaning

The thermometer should be cleaned by hand washing only. Do not immerse the device in water.

Storage

Store the thermometer in a cool, dry place. Avoid extreme temperatures and direct sunlight.

TROUBLESHOOTING

If you encounter any issues with your thermometer, please refer to the detailed user manual included in the packaging or contact customer support.

PRODUCT SPECIFICATIONS

- **Product Dimensions:** 5.55 x 2.3 x 1.45 inches; 1.76 ounces
- **Item model number:** Temple Touch Digital Thermometer
- **Date First Available:** January 31, 2023
- **Manufacturer:** Amazon.com Services LLC.
- **ASIN:** B09MFKK719
- **Brand:** Amazon Basic Care
- **Special Feature:** Fast Reading System
- **Color:** White
- **Age Range (Description):** All Ages
- **Included Components:** Thermometer
- **Outer Material:** Plastic
- **Specification Met:** UL
- **Display Type:** Digital
- **Connectivity Technology:** Digital
- **Product Care Instructions:** Hand Wash Only



Image: The Amazon Basic Care Temple Touch Digital Thermometer displayed next to a ruler, indicating its physical dimensions in both centimeters and inches.

WARRANTY AND CUSTOMER SUPPORT

This product includes warranty and service information within the packaging. For questions or comments, please contact Amazon Basic Care customer support at 1-877-485-0385.