

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

[manuals.plus](#) /

> [AUSSNICE](#) /

> [AUSSNICE Cordless Electric Pill Crusher Grinder Pulverizer - Instruction Manual](#)

AUSSNICE 10019A

AUSSNICE Cordless Electric Pill Crusher Grinder Pulverizer

MODEL: 10019A

[Introduction](#)

[Safety Information](#)

[Product Components](#)

[Setup](#)

[Operation](#)

[Maintenance](#)

[Warranty &](#)

[Support](#)

[Troubleshooting](#)

[Specifications](#)

1. Introduction

The AUSSNICE Cordless Electric Pill Crusher Grinder Pulverizer is designed to effortlessly crush various medications, vitamins, and supplements into a fine powder. This device is particularly beneficial for individuals who experience difficulty swallowing pills, including the elderly, children, and pets. Its cordless design and simple one-button operation ensure convenience and ease of use.



Image 1: Overview of the AUSSNICE Cordless Electric Pill Crusher highlighting its features.

2. Safety Information

- Read all instructions before using the appliance.
- Keep the appliance out of reach of children.
- Ensure the grinding cup is securely attached and the lid is properly closed before operation to prevent splashing.
- Do not operate the appliance for more than 20 seconds continuously to prevent overheating. The device has an automatic overheat protection feature.
- Avoid contact with moving parts. Always unplug or switch off the device before cleaning or maintenance.
- Do not immerse the main body of the appliance in water or other liquids. The grinding cup is removable and washable.
- Use only the provided USB charging cable.
- This appliance is intended for crushing pills and tablets only. Do not use it for other purposes.

3. Product Components

Familiarize yourself with the parts of your pill crusher:

- **Main Body:** Contains the motor and battery.
- **Removable Grinding Cup:** Stainless steel cup with blades for crushing pills.
- **Sealing Lid:** Inner lid to prevent powder leakage.
- **Transparent Top Lid:** Outer lid that activates the grinding mechanism and allows observation.
- **USB Charging Cable:** For recharging the internal battery.
- **Dual-use Brush and Spoon:** For cleaning and scooping powder.

Transparent Top Lid

Visible grinding process



Removable Grinding Cup

5.7oz large capacity, dishwasher safe



Sealing Lid

No more powder messy

Safety Lock

Easy and safe locking



Main Body

Image 2: Exploded view illustrating the main components of the pill crusher.

4. Setup

4.1 Initial Charging

Before first use, fully charge the pill crusher. Connect the provided USB charging cable to the charging port on the main body and plug the other end into a compatible USB power source. The battery indicator lights will show the charging progress.



Image 3: Battery indicator on the main body.

4.2 Assembly

1. Ensure the main body is clean and dry.
2. Place the removable grinding cup onto the main body. Twist clockwise to lock it securely into place.
3. Place the pills into the grinding cup.
4. Cover the grinding cup with the sealing lid, then place the transparent top lid over it.

5. Operating Instructions

Follow these steps to crush your pills:

1. **Prepare Pills:** Place the desired amount of pills into the removable grinding cup. The large capacity can hold over 100 pills at once.
2. **Assemble:** Securely attach the grinding cup to the main body by twisting clockwise. Place the sealing lid, then the transparent top lid.
3. **Activate Grinding:** Gently press down on the transparent top lid. The grinder will start operating.
4. **Control Fineness:** Continue pressing the lid for varying durations to achieve your desired powder consistency (coarse, medium, or fine). The device will automatically stop after 20 seconds of continuous operation to prevent overheating. Release and press again if further grinding is needed.
5. **Retrieve Powder:** Once grinding is complete, lift your hand from the lid. Carefully remove the transparent top lid, then the sealing lid, and finally detach the grinding cup from the main body by twisting counter-clockwise. Pour the fine powder into a suitable container.



Image 4: Step-by-step guide for operating the pill crusher.



Image 5: Example of coarse versus fine powder consistency.

5.1 Video Demonstration

Your browser does not support the video tag.

Video 1: Official demonstration of the AUSSNICE Cordless Electric Pill Crusher Grinder in use, showing the grinding process and results.

6. Maintenance

6.1 Cleaning

The removable grinding cup can be rinsed directly under running water for easy cleaning. Use the provided dual-use brush to clean any residual powder from the blades and inside the cup. Ensure all parts are completely dry before reassembly or storage.



Image 6: Cleaning the removable grinding cup under running water.

Do not wash the main body of the appliance. Wipe it with a damp cloth if necessary.

7. Troubleshooting

Problem	Possible Cause	Solution
Device does not turn on.	Low battery; Grinding cup or lid not properly secured.	Charge the device fully. Ensure the grinding cup is twisted clockwise to lock and the top lid is pressed down firmly.
Pills are not finely ground.	Insufficient grinding time; Too many pills in the cup.	Press the lid for a longer duration (up to 20 seconds per cycle). Reduce the quantity of pills per grinding cycle.
Device stops during operation.	Overheat protection activated.	Allow the device to cool down for a few minutes before resuming operation. Avoid continuous grinding for more than 20 seconds.
Powder leaks during grinding.	Lids not properly sealed.	Ensure both the sealing lid and transparent top lid are correctly placed and pressed down firmly before starting.

8. Specifications

- **Model Number:** 10019A
- **Product Dimensions:** 3.35 x 3.34 x 6.89 inches
- **Item Weight:** 1.53 Pounds
- **Battery:** 1 Lithium Ion battery (included), 2000mAh
- **Motor Power:** 200W
- **Blade Material:** 304 Stainless Steel
- **Grinding Capacity:** 5.7 oz (approx. 170ml), holds over 100 pills
- **Operation:** One-button, lid-activated switch
- **Safety Features:** Overheating protection, double sealing cover

9. Warranty & Support

For warranty information or technical support, please refer to the product packaging or contact AUSSNICE customer service directly. Keep your purchase receipt for warranty claims.