Unold 86010

UNOLD 86010 Wärm Electric Hot Water Bottle User Manual

Model: 86010 | Brand: Unold

1. INTRODUCTION

This manual provides comprehensive instructions for the safe and effective use of your UNOLD 86010 Wärm Electric Hot Water Bottle. Please read all instructions carefully before use and retain this manual for future reference. This electric hot water bottle is designed to provide soothing warmth for various body areas such as the stomach, back, and neck, offering a safe and convenient alternative to traditional hot water bottles.

2. IMPORTANT SAFETY INSTRUCTIONS

To ensure safe operation and prevent injury or damage, always follow these safety precautions:

- **No Water Filling Required:** This electric hot water bottle is pre-filled and sealed. Do not attempt to open or fill it with water, as this can be dangerous and will void the warranty.
- **Automatic Shut-off:** The device features an automatic shut-off function that activates when the temperature reaches approximately 65°C. This prevents overheating.
- **Overheating Protection:** An integrated overheating protection mechanism further enhances safety by preventing excessive temperatures.
- **Electrosmog-Free After Charging:** Once fully charged and disconnected from the power supply, the hot water bottle emits no electrosmog.
- **Power Cord and Charger:** Only use the provided charger and power cord. Inspect the cord and charger for any damage before each use. Do not use if damaged.
- **Keep Dry:** Do not immerse the hot water bottle or charger in water or any other liquid. Keep away from wet environments.
- **Supervision:** Do not leave the device unattended while charging. Keep out of reach of children and individuals with reduced physical, sensory, or mental capabilities unless supervised.
- **Intended Use:** Use the product only for its intended purpose of providing warmth. Do not use for medical diagnosis or treatment without consulting a healthcare professional.
- **Storage:** Store the device in a cool, dry place when not in use.

3. PRODUCT OVERVIEW

The UNOLD 86010 Wärm Electric Hot Water Bottle is a wireless heating device designed for comfort and ease of use. It comes with a soft fabric cover and can be used with an optional belly/back bandage for secure placement.

Image 3.1: The UNOLD 86010 Wärm Electric Hot Water Bottle in beige. This image shows the main product unit, a soft, rectangular heating pad.

Components:

- Electric Hot Water Bottle Unit
- Charging Cable with Power Adapter
- Fabric Bag / Cover
- Belly / Back Bandage (optional accessory for secure wearing)

Image 3.2: This image illustrates the UNOLD 86010 Wärm Electric Hot Water Bottle being worn with the belly/back bandage, alongside the separate charging unit.

4. SETUP AND CHARGING

Before first use, and whenever the device needs reheating, follow these steps:

1. **Prepare for Charging:** Ensure the hot water bottle is placed on a flat, non-flammable surface.
2. **Connect Charger:** Locate the charging port on the hot water bottle. Connect the charging cable to the port, ensuring a secure fit.

Image 4.1: Close-up view of the charging port located on the electric hot water bottle.

3. **Plug In:** Plug the power adapter into a standard electrical outlet. The charging indicator light on the charger will illuminate, indicating that charging has begun.

Image 4.2: The charging unit for the UNOLD 86010 Wärm Electric Hot Water Bottle.

4. **Charging Time:** The device typically takes approximately 10 minutes to heat up to its operating temperature of about 65°C. The indicator light will turn off or change color once charging is complete and the automatic shut-off activates.
5. **Disconnect:** Once fully charged, disconnect the power adapter from the outlet first, then disconnect the charging cable from the hot water bottle. The device is now ready for use.

5. OPERATING INSTRUCTIONS

After charging, the UNOLD 86010 Wärm Electric Hot Water Bottle is ready to provide warmth.

1. **Placement:** Place the heated hot water bottle on the desired body area, such as the stomach, back, or neck.

Image 5.1: A person comfortably using the electric hot water bottle on their stomach for soothing warmth.

2. **Using the Bandage:** If using the optional belly/back bandage, insert the hot water bottle into the bandage pocket and secure it around your body for hands-free use.

Image 5.2: This image shows a person using the electric hot water bottle on their neck, highlighting key features such as automatic shut-off, overheating protection, 10-minute heating time to 65°C, no electrosmog, and several hours of warmth depending on ambient temperature.

3. **Duration of Warmth:** The hot water bottle will maintain its warmth for several hours, with the exact duration depending on the ambient temperature and how it is used.
4. **Reheating:** When the warmth diminishes, repeat the charging process as described in Section 4.

6. MAINTENANCE AND STORAGE

Proper care ensures the longevity and safe operation of your electric hot water bottle.

- **Cleaning:** Disconnect the device from power before cleaning. Wipe the surface of the hot water bottle and charger with a soft, damp cloth. Do not use abrasive cleaners or immerse the device in water.
- **Fabric Cover:** The fabric bag/cover can be gently hand-washed or machine-washed on a delicate cycle with mild detergent. Ensure it is completely dry before reattaching to the hot water bottle.
- **Storage:** When not in use, store the hot water bottle and charger in a cool, dry place, away from direct sunlight and extreme temperatures. Ensure the device is completely cool before storing.

7. TROUBLESHOOTING

If you encounter issues with your UNOLD 86010 Wärme Electric Hot Water Bottle, refer to the following common problems and solutions:

Problem	Possible Cause	Solution
---------	----------------	----------

Problem	Possible Cause	Solution
Device does not heat up.	• Charger not properly connected. • Power outlet not working. • Charger or cable is faulty.	• Ensure charger is securely connected to both the hot water bottle and the power outlet. • Test the power outlet with another device. • Inspect charger and cable for damage. If damaged, discontinue use and contact support.
Device heats up but quickly cools down.	• Insufficient charging time. • Ambient temperature is very low.	• Ensure the device charges for the full 10 minutes until the indicator light signals completion. • In colder environments, the heat retention may be shorter. Consider using the fabric cover or bandage for better insulation.
Charging indicator light does not turn on.	• Charger not connected. • Power outlet issue. • Faulty charger.	• Verify all connections. • Check power outlet. • Contact customer support if the charger appears faulty.

If the problem persists after attempting these solutions, please contact Unold customer support for further assistance.

8. SPECIFICATIONS

Model Number:	UNOLD-86010
Brand:	Unold
Product Dimensions:	10.63 x 7.28 x 1.97 inches (27 x 18.5 x 5 cm)
Item Weight:	3.09 Pounds (1.4 kg)
Material:	Nylon (outer cover)
Maximum Temperature:	Approx. 65°C
Heating Time:	Approx. 10 minutes
Safety Features:	Automatic Shut-off, Overheating Protection

9. WARRANTY AND SUPPORT

For warranty information, technical support, or service inquiries, please contact Unold customer service. Details

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