

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

› [Braun](#) /

› [Braun Silk-expert Pro 3 PL3133 IPL Hair Removal Device User Manual](#)

Braun PL3133

Braun Silk-expert Pro 3 PL3133 IPL Hair Removal Device

User Manual

1. INTRODUCTION

The Braun Silk-expert Pro 3 PL3133 is an Intense Pulsed Light (IPL) hair removal device designed for safe and effective hair reduction at home. It utilizes advanced IPL technology to provide long-lasting smooth skin, with visible results typically seen within 4 weeks of consistent use. The device features Skin pro (SensoAdapt™) technology, which automatically adapts to your skin tone to ensure optimal power and safety.

This manual provides essential information on how to use and maintain your device for best results and safety.

2. SAFETY INFORMATION

Please read all safety warnings and instructions carefully before using the device. Failure to follow these instructions may result in injury or damage to the device.

- **Do not use the device with excessive pressure.** Apply gentle contact to the skin.
- **Keep the device away from water.** Do not use it near bathtubs, showers, basins, or other vessels containing water.
- **Do not use the device while it is connected to a power source if it is designed for battery operation.** (Note: This device is primarily corded for operation. Refer to power instructions for proper use.)
- **Consult a doctor before use** if you have any pre-existing skin conditions, medical conditions, or are taking medications that may affect skin sensitivity.
- **Do not use on irritated, sunburned, or broken skin.**
- **Avoid direct eye exposure to the IPL flashes.** The device is designed with a skin contact sensor to prevent accidental flashes when not in full contact with the skin.
- **Keep out of reach of children.**

3. PACKAGE CONTENTS

Ensure all items are present in the package:

- 1x Braun Silk-expert Pro 3 IPL Device
- 1x Standard Head
- 1x Precision Head
- 1x Venus Extra Smooth Razor
- 1x Storage Pouch
- 1x Power Supply (corded)



Image: Contents of the Braun Silk-expert Pro 3 PL3133 package, including the IPL device, two treatment heads, a Venus razor, and a storage pouch.

4. SETUP

1. **Read the User Manual:** Before first use, thoroughly read this entire user manual to understand the device's functions, safety precautions, and usage guidelines.
2. **Prepare the Device:** Connect the power supply cord to the device and then plug it into a suitable electrical outlet. The device is designed for corded operation.
3. **Select the Treatment Head:** Choose the appropriate head for the area you intend to treat. The standard head is suitable for larger areas like legs and arms, while the precision head is ideal for smaller, sensitive areas such as the face and bikini line.

5. OPERATING INSTRUCTIONS

Follow these steps for effective and safe IPL hair removal:

1. **Prepare Your Skin:** Ensure the skin area to be treated is clean, dry, and free of any lotions, creams, or deodorants. Shave the area thoroughly before each treatment session.
2. **Power On:** Press the power button to turn on the device. The indicator lights will illuminate.
3. **Skin Tone Adaptation:** The Braun Silk-expert Pro 3 features Skin pro (SensoAdapt™) technology, which continuously reads your skin tone and automatically adjusts the light intensity for optimal safety and efficacy. You do not need to manually adjust the intensity.
4. **Apply to Skin:** Place the treatment window of the device firmly against your skin. Ensure full contact to activate the skin contact sensors. The device will not flash if not in full contact.
5. **Initiate Flash:** Press the flash button. A light flash will be emitted.
6. **Treatment Modes:**
 - **Gliding Mode:** For larger areas, hold down the flash button and slowly glide the device over your skin. The device will emit continuous flashes.
 - **Stamp Mode:** For smaller or sensitive areas, place the device on the skin, press and release the flash button, then lift and move to the next area.
 - **Gentle Mode:** For sensitive areas or first-time users, activate the gentle mode (if available on your model) to reduce light intensity for a more comfortable treatment.
7. **Complete Treatment:** Continue treating the desired areas. The device can deliver up to 100 flashes per minute for fast treatment. A full body treatment can be completed in approximately 13 minutes.
8. **Power Off:** After completing your session, press and hold the power button to turn off the device. Unplug it from the electrical outlet.

Güvenli ve neredeyse acısız.

Skin-Protection LITE sensör güvenli ve rahat bakım sağlar.



Cilt tonlarınızı okur.

Her flaşın gücünü otomatik olarak uyarlar.

3 güç seviyesi.

Image: The Skin-Protection LITE sensor automatically reads skin tones and adjusts flash power for safe and comfortable treatment across 3 power levels.

Etkili ve Güçlü.

⌚ 13 dk Sadece 13 dakikada tam vücut* bakımı.



Sadece haftada bir bakım



Neredeyse acısız.



Özel bölgelerde bile.



Erkekler için uygundur.



* *Bacaklar (Alt ve üst) Koltuk Altı, Bikini Bölgesi, Üst dudak

Image: The device offers effective and powerful treatment, completing a full body session in approximately 13 minutes with weekly use. It is almost painless, suitable for sensitive areas, and can be used by men.

Bir kez alın,
istediđiniz sıklıkta
bakım yapın.



- ✓ Tam vücut bakımı, bikini **bölgelerinde bile.**
- ✓ Aylık masraf yok, ek yenileme gerekmez
- ✓ Zaman kazanın, evinizde bakım yapın.
- 🇩🇪 Uzun süre dayanıklı. Almanya'da **tasarlanmıştır.**

Image: The Braun Silk-expert Pro 3 provides full body care, including bikini areas, with no monthly expenses. It saves time by allowing home treatments and is designed in Germany for durability.

Uzun süre kalıcı pürüzsüz cilt.
Klinik olarak **kanıtlanmış etkinlik.**



Image: Clinically proven effectiveness of the Braun Silk-expert Pro 3, showing significant hair reduction on both armpits and legs after consistent use.

Silk·expert
Pro

3

1 yıla kadar
pürüzsüz
cilt.*

3 kullanımda gözle
görülür sonuçlar.

*Önerilen program uygulandığında. Bireysel sonuçlar farklılık gösterebilir.



Image: Achieve up to 1 year of smooth skin with the Braun Silk-expert Pro 3, with visible results in as few as 3 uses.

6. MAINTENANCE

Proper maintenance ensures the longevity and optimal performance of your device:

- **Cleaning:** After each use, ensure the device is turned off and unplugged. Use a soft, dry cloth to wipe the treatment window and the device body. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the device and its accessories in the provided pouch in a cool, dry place, away from direct sunlight and moisture.
- **Inspection:** Regularly inspect the device and power cord for any signs of damage. If any damage is found, discontinue use and contact customer service.

7. TROUBLESHOOTING

If you encounter issues with your Braun Silk-expert Pro 3, refer to the following common troubleshooting tips:

- **Device not turning on:** Ensure the power cord is securely connected to both the device and a working electrical outlet.
- **Device not flashing:** Check that the treatment window is in full and firm contact with your skin. The skin contact sensors must be fully engaged for the device to flash. Ensure the device is not in an error state (refer to indicator lights).
- **Skin irritation:** If you experience redness or irritation, ensure you are using the device on clean, shaved skin. Consider using the gentle mode or reducing the frequency of treatments. If irritation persists, discontinue use and consult a healthcare professional.
- **No visible results:** IPL treatment requires consistent use over several weeks. Ensure you are following the recommended treatment schedule. Results vary by individual and hair type.

For further assistance, please refer to the full user manual or contact Braun customer support.

8. SPECIFICATIONS

Feature	Specification
Model Number	PL3133 (Article No. 411925)
Brand	Braun
Dimensions (L x W x H)	10 x 10 x 15 cm
Weight	1.11 kg (1113.4 Grams)
Color	White
Power Source	Corded Electric (Mains Powered)
Special Features	Portable, SensoAdapt™ Skin Tone Sensor, Gentle Mode
Recommended Use	Face, Body
Included Components	Precision Head, Storage Pouch, Venus Razor

9. WARRANTY AND SUPPORT

For detailed warranty information, please refer to the warranty card included with your product or visit the official Braun website. Braun offers customer support for any questions or issues you may encounter with your device.

Online Support: Visit the official Braun website for FAQs, troubleshooting guides, and product registration.

Contact Customer Service: Refer to the contact information provided in your product packaging or on the Braun website for direct assistance.