



Manuals.plus /

- › XSOUL /
- › XSOUL Dermaflash Pro IPL Hair Removal Device User Manual

## XSOUL SN02

# XSOUL Dermaflash Pro IPL Hair Removal Device

## USER MANUAL

### 1. Introduction

The XSOUL Dermaflash Pro is an advanced Intense Pulsed Light (IPL) hair removal device designed for at-home use. It offers a convenient and effective method for long-term hair reduction on various body areas, including armpits, back, legs, arms, face, and bikini line. This device utilizes IPL technology to target hair follicles, leading to a significant reduction in hair growth over time. It is a corded device, ensuring consistent power during use.



Image: The XSOUL Dermaflash Pro IPL hair removal device, shown with its accompanying green carrying bag and product packaging.

## 2. Important Safety Information

Please read all safety instructions carefully before using the device. Failure to follow these instructions may result in injury or damage to the device.

- The XSOUL IPL hair removal device is a plug-in device and is not wireless. Ensure it is properly connected to a power source during operation.
- Do not use the device if your skin has been artificially or naturally tanned, or if it is burned from sun overexposure.
- Do not use the device during pregnancy or menstruation.
- Do not use the device on freckles or on tattoos.
- Start treatment at Light Intensity Level 1 and gradually increase as needed, based on your skin sensitivity and

comfort.

### Skin Tone & Body Hair Type Chart

The effectiveness and safety of IPL treatment depend on your skin tone and hair color. Refer to the chart below to determine if this device is suitable for you.

		<b>SKIN TONE &amp; BODY HAIR TYPE CHART</b>					
Skin Tone \ Hair Color							
White/Grey		✗	✗	✗	✗	✗	✗
Light Blonde		✗	✗	✗	✗	✗	✗
Red		✗	✗	✗	✗	✗	✗
Dark Blonde Light Brown		✓	✓	✓	✓	✓	✗
Brown		✓	✓	✓	✓	✓	✗
Dark Brown		✓	✓	✓	✓	✓	✗
Black		✓	✓	✓	✓	✓	✗

✓ = It's match!   ✗ = Try a different hair removal method

## Safe, Gentle, Effective

**NOTE:** Different skin sensitivity and exposure settings may lead to different user experience.

Start from Light Intensity Level 1 and work your way up as needed.

**DO NOT** use if your skin has been artificially or naturally tanned or is burned from sun overexposure.

**DO NOT** use during pregnancy or menstruation.

**DO NOT** use on freckles or on tattoos.

Image: A chart illustrating compatible skin tones and hair colors for IPL hair removal. Green checkmarks indicate suitability, while red crosses indicate incompatibility.

## IPL Compatibility Guide

Hair Color	Skin Tone Compatibility
White/Grey, Light Blonde, Red	Not suitable for IPL treatment.
Dark Blonde, Light Brown, Brown, Dark Brown, Black	Generally suitable for IPL treatment, depending on specific skin tone. Refer to the chart for details.

Note: IPL is most effective on darker hair colors and lighter skin tones. Consult the provided chart for precise compatibility.

### 3. Product Features

- **Advanced IPL Technology:** Features 6 energy levels and a maximum output of 15.3 Joules (J) for effective hair reduction.
- **High Flash Capacity:** Boasts over 999,999 flashes, providing years of reliable use.
- **Dual Operating Modes:** Includes both automatic and manual modes for comprehensive hair removal across different body areas.
- **Enhanced Design:** The 2023 Enhanced Version features a deeper outlet to protect the skin during treatment.
- **Painless Application:** Designed for a comfortable and painless hair removal experience.
- **Dermatologist Approved:** The device has been dermatologically tested and approved for safe use.

# STRONGER ENERGY FASTER HAIR REMOVAL



hair reduction  
in 3 weeks



15.25J higher  
energy output

MEDICAL GRADE  
SKIN CARE



Image: Illustration highlighting the device's strong energy output (15.25J) and potential for hair reduction within 3 weeks.

  
permanent  
hair removal

  
eyes  
protection

  
years  
of usage

  
safe &  
painless

  
6 energy  
levels

**AUTO**  
  
manual &  
auto modes

# 6 LEVELS 2 MODES



Image: Close-up of the device's display, indicating over 999,999 flashes available, 6 energy levels, and the presence of manual and auto modes.

# OLD VERSION

**HIGH HEAT  
SHALLOW**

# NEW VERSION

**DEEPER OUTLET  
PROTECT SKIN**



Image: Visual comparison between an older and newer version of the IPL device, emphasizing the deeper light outlet in the new version for enhanced skin protection.

# DERMATOLOGIST APPROVED



Image: The IPL device in use on a person's leg, accompanied by text confirming it is "Dermatologist Approved."

## 4. Setup

1. **Prepare the Skin:** Before treatment, ensure the skin area is clean, dry, and free of any lotions, creams, or deodorants. Shave the hair in the treatment area to ensure the light energy is absorbed by the hair follicle, not the surface hair.
2. **Connect Power:** Plug the device into a suitable electrical outlet. The device is corded and requires continuous power during operation.
3. **Wear Eye Protection:** Always wear the provided protective eyewear (or suitable alternatives) to shield your

eyes from the intense light flashes.

4. **Inspect Device:** Before each use, inspect the device for any signs of damage. Do not use if damaged.



Image: The complete set of items included with the XSOUL Dermaflash Pro, showing the device, power adapter, user manual, disposable razors, and protective glasses.

## 5. Operation

1. **Power On:** Press and hold the power button to turn on the device.
2. **Select Energy Level:** Use the appropriate buttons to select your desired energy level (1-6). Begin with Level 1 and increase gradually based on comfort and skin reaction.
3. **Choose Mode:** The device offers both automatic and manual modes.
  - **Manual Mode:** Ideal for small or sensitive areas like the face, armpits, or bikini line. Press the flash button for each pulse.
  - **Automatic Mode:** Suitable for larger areas like legs or back. Hold the flash button, and the device will automatically emit flashes when in full contact with the skin.
4. **Apply to Skin:** Place the device's treatment window firmly against the skin. Ensure full contact for the device to flash.
5. **Perform Treatment:** For manual mode, press the flash button. For automatic mode, glide the device slowly over the skin, ensuring continuous contact. Overlap treatment areas slightly to avoid missing spots.
6. **Treatment Frequency:** For optimal results, it is generally recommended to use the device 2 times per week for the initial few weeks. Hair reduction may be noticeable within 3 weeks of consistent use.
7. **Post-Treatment Care:** After treatment, you may apply a soothing lotion or gel to the treated area if desired. Avoid sun exposure on treated areas immediately after use.



Image: A composite image demonstrating the versatility of the IPL device for use on various body parts, including armpits, arms, legs, and the face.

## 6. Maintenance

- **Cleaning:** After each use, wipe the treatment window with a clean, dry, lint-free cloth. Do not use abrasive cleaners or solvents. Ensure the device is unplugged before cleaning.
- **Storage:** Store the device in its original packaging or the provided carrying bag in a cool, dry place, away from direct sunlight and moisture. Keep out of reach of children.

- **Cord Care:** Do not wrap the power cord tightly around the device. Store the cord loosely to prevent damage.

## 7. Troubleshooting

Problem	Possible Cause	Solution
Device does not turn on.	Not properly plugged in; power outlet issue; device malfunction.	Ensure the power cord is securely plugged into both the device and a working electrical outlet. Try a different outlet. If the issue persists, contact customer support.
Device flashes intermittently or not at all.	Treatment window not in full contact with skin; device in manual mode; flash count depleted.	Ensure the treatment window is pressed firmly and completely against the skin. If in manual mode, press the flash button. Check the remaining flash count on the display.
Skin irritation after use.	Energy level too high; sensitive skin; improper skin preparation.	Reduce the energy level for subsequent treatments. Ensure skin is clean and dry before use. Apply a soothing lotion after treatment. If irritation persists, discontinue use and consult a dermatologist.

## 8. Specifications

- **Model Name:** SN02
- **Product Dimensions:** 6.5 x 4.9 x 3.07 inches
- **Item Weight:** 9.12 ounces
- **Power Source:** Corded Electric
- **Special Feature:** Portable
- **Manufacturer:** XSOUL
- **ASIN:** B09WH69WSR

## 9. Warranty and Support

XSOUL stands behind the quality of its products. Your Dermaflash Pro IPL Hair Removal Device comes with comprehensive after-sales support.

- **Warranty:** The product includes a 36-month (3-year) warranty from the date of purchase.
- **Money-Back Guarantee:** A 100-day money-back guarantee is offered, ensuring your satisfaction.

For any questions, concerns, or to inquire about returns or exchanges, please contact our professional customer service team:

**Email:** [help@xsoulofficial.com](mailto:help@xsoulofficial.com)

