

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

> [MLAY](#) /

> [MLAY IPL Hair Removal Device User Manual - Model 1](#)

## MLAY 1

# MLAY IPL Hair Removal Device User Manual

Model: 1

<a href="#">Introduction</a>	<a href="#">Safety Information</a>	<a href="#">Package Contents</a>	<a href="#">Device Overview</a>	<a href="#">Setup</a>
<a href="#">Operation</a>	<a href="#">Maintenance</a>	<a href="#">Troubleshooting</a>	<a href="#">Specifications</a>	<a href="#">Support</a>

## 1. INTRODUCTION

Welcome to the MLAY IPL Hair Removal Device user manual. This guide provides essential information for the safe and effective use of your device.

The MLAY IPL Hair Removal Device utilizes advanced Intense Pulsed Light (IPL) technology to target hair follicles, aiming for permanent hair reduction. It is designed for home use, offering a convenient solution for hair removal on various body areas.

## 2. IMPORTANT SAFETY INFORMATION

Read all instructions before use. Failure to follow safety guidelines may result in injury or device damage.

### 2.1 General Safety Precautions:

- Always wear the provided safety glasses during operation to protect your eyes from intense pulsed light.
- Do not use on tattooed areas, moles, or birthmarks.
- Avoid direct exposure to the eyes.
- Keep out of reach of children.
- Consult a physician before use if you have any skin conditions, medical implants, or are pregnant/nursing.
- Do not use on recently tanned skin or if you have a history of skin sensitivity to light.

### 2.2 Skin and Hair Type Compatibility:

The MLAY IPL device is most effective on lighter skin tones with darker hair. It is not suitable for very dark skin tones or very light hair colors (white, gray, red, golden) as the IPL technology targets melanin in the hair follicle.



Image description: A table illustrating skin tone and hair color compatibility with IPL hair removal. It shows that the device is applicable for white, ivory tint, beige, and light brown skin with black, dark brown, and brown hair. It is not applicable for dark brown/black skin or white, gray, red, or golden hair.

### 3. PACKAGE CONTENTS

---

Verify that all items are present in your package:

- 1 x MLAY IPL Hair Removal Device (Main Body)
- 1 x AC Adapter
- 1 x Safety Glasses
- 1 x Razor
- 1 x User Manual

### 4. DEVICE OVERVIEW

---

Familiarize yourself with the components of your MLAY IPL Hair Removal Device.



Image description: Front view of the MLAY IPL Hair Removal Device, showing the power button, intensity level indicators, and the main flash button.



Image description: Rear view of the MLAY IPL Hair Removal Device, showing ventilation holes.



Image description: The MLAY IPL Hair Removal Device connected to its power adapter, with the light emission window visible at the bottom.

#### 4.1 Key Components:

- **Power Interface:** For connecting the AC adapter.
- **Gear Display:** Indicates the selected power level.
- **Power Button:** Long press to turn on/off; short press to switch power levels.
- **Lighting Button:** Long press for automatic light output (Auto Mode); short press for manual light output (Manual Mode).
- **Out of the Mouth of Light (Flash Window):** Area where the IPL light is emitted.

## Single flash or double flash depending on muscle

Choose sparkle with your skin



Image description: A diagram highlighting the power interface, gear display, power button, lighting button, and flash window on the MLAY IPL device.

### 4.2 Internal Structure:

The device features a high-quality quartz lamp tube and filter glass for stable energy output and skin protection.



Image description: An exploded view diagram showing the internal components of the IPL device, specifically the quartz burner and the filter glass.

## 5. SETUP

1. **Prepare the Skin:** Shave the area to be treated two hours before use. Ensure the skin is clean, dry, and free of any lotions, creams, or deodorants.
2. **Connect Power:** Plug the AC adapter into the device and then into a power outlet.
3. **Wear Safety Glasses:** Always wear the provided safety glasses before turning on the device.

## 6. OPERATING INSTRUCTIONS

### 6.1 Power On/Off and Intensity Levels:

1. **Power On:** Long press the Power Button to turn on the device. The gear display lights will illuminate.
2. **Adjust Intensity:** Short press the Power Button to cycle through the 5 available power levels. It is recommended to start at Power Level 1 for your first treatment and gradually increase as your skin tolerates it.
3. **Power Off:** Long press the Power Button again to turn off the device.

## Suitable skin tone for use

Skin and hair color comparison table:

Skin color	Hair color	Black	Dark brown	Brown	Golden	Red	White	Gray
White		✓	✓	✓	✓	✗	✗	✗
Ivory tint		✓	✓	✓	✓	✗	✗	✗
Beige		✓	✓	✓	✓	✗	✗	✗
Light brown		✓	✓	✓	✓	✗	✗	✗
Brown		✓	✓	✓	✓	✗	✗	✗
Dark brown / black		✗	✗	✗	✗	✗	✗	✗

✓ Applicable ✗ Not applicable

Image description: A diagram illustrating three energy levels for the IPL device, labeled 'Low level care delicate skin,' 'High-grade effect is remarkable,' and an intermediate level. It also shows icons for various body parts where the device can be used.

### 6.2 Treatment Modes (Auto/Manual):

The device offers two modes for light emission:

- **Manual Mode:** Short press the Lighting Button to emit a single flash. This mode is suitable for smaller, precise areas like the upper lip, chin, or bikini line.
- **Auto Mode:** Long press the Lighting Button to activate continuous flashing. This mode is ideal for larger areas such as legs, arms, back, and chest.

### 6.3 Treatment Procedure:

1. Place the flash window firmly against the skin in the desired treatment area. Ensure full contact.
2. Select your preferred mode (Manual or Auto) and activate the flash.
3. Move the device slowly and methodically across the skin, ensuring each section is treated without overlapping or missing spots. Treat the entire area 3-4 times per session.
4. After completing the session, apply a plant-based moisturizer to the treated area.



Image description: A woman using the MLAY IPL Hair Removal Device on her leg, demonstrating the application process.

## 6.4 Treatment Schedule:

Consistent use is key for optimal results. A typical course of treatment involves 7-9 sessions.

- **Weeks 1-3:** One treatment per week.
- **Weeks 4-9:** One treatment every 2-3 weeks.
- **Maintenance:** After the initial course, perform a treatment every two months as needed for hair regrowth.

Individual results may vary. Some users may experience redness after treatment, which is normal and usually subsides quickly. Apply sunblock if treated areas will be exposed to the sun.

## 7. MAINTENANCE AND CARE

- **Cleaning:** After each use, ensure the device is powered off and unplugged. Gently wipe the flash window with a soft, dry cloth. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the device in a cool, dry place, away from direct sunlight and moisture.

## 8. TROUBLESHOOTING

If you encounter issues with your MLAY IPL Hair Removal Device, refer to the following common problems and solutions:

Problem	Possible Cause	Solution
Device does not turn on.	Not properly plugged in; power outlet issue.	Ensure the AC adapter is securely connected to both the device and a working power outlet.
No flash emitted.	Flash window not in full contact with skin; device in standby mode.	Ensure the flash window is pressed firmly against the skin. Short press the Lighting Button for manual flash or long press for Auto Mode.
Skin irritation or discomfort.	Intensity level too high; unsuitable skin type.	Reduce the intensity level. Review the skin and hair type compatibility chart. Discontinue use if irritation persists.

For further assistance, please contact customer support.

## 9. SPECIFICATIONS

- **Model Number:** 1
- **Product Dimensions:** 6.25 x 2.8 x 1.63 inches
- **Weight:** 1.94 Pounds
- **Manufacturer:** Shenzhen Mareal Technology Co., Ltd.
- **Technology:** IPL (Intense Pulsed Light)
- **Power Levels:** 5
- **Flash Capacity:** Unlimited Flashes
- **Modes:** Manual, Auto

## 10. WARRANTY AND SUPPORT

---

For any product-related inquiries, technical support, or warranty claims, please contact MLAY customer service. Refer to the contact information provided in your product packaging or visit the official MLAY website.

We aim to provide a satisfactory solution as soon as possible for any concerns you may have.

---