

## Oral-B 427063

# Oral-B Vitality Pro Rechargeable Electric Toothbrush User Manual

Model: 427063

## 1. INTRODUCTION

---

The Oral-B Vitality Pro electric toothbrush is designed to provide a superior clean compared to a traditional manual toothbrush. It features professional cleaning technology, offering effective plaque removal and gum protection. This manual provides essential instructions for the safe and effective use of your new electric toothbrush.

## 2. SAFETY INFORMATION

---

Please read all instructions before using this appliance. Keep this manual for future reference.

- This product is for external oral use only.
- Avoid contact with eyes.
- Keep away from open flames or direct heat sources.
- Do not immerse the charger in water or other liquids.
- Do not modify or repair the product. Contact authorized service for repairs.
- Supervise children when using the toothbrush.

## 3. SETUP

---

### 3.1 Unpacking

Carefully remove all components from the packaging. Ensure you have the toothbrush handle, one brush head, and the charger.

### 3.2 Charging the Toothbrush

Before first use, charge the toothbrush for at least 16 hours. For subsequent charges, a full charge typically takes 12-16 hours.

1. Plug the charger into an electrical outlet.

2. Place the toothbrush handle on the charging unit. The charging indicator light (if present) will illuminate.
3. Once fully charged, the toothbrush is ready for use.

*Note: The product packaging may contain different information than shown on our site. Please always read the label, warnings, and instructions provided on the product before using or consuming it.*

## 4. OPERATING INSTRUCTIONS

---

### 4.1 Attaching the Brush Head

Firmly push the brush head onto the toothbrush handle until it clicks into place.

### 4.2 Brushing Your Teeth

1. Wet the brush head and apply a small amount of toothpaste.
2. Place the brush head on your teeth, holding it at a slight angle towards the gumline.
3. Press the power button to turn on the toothbrush.
4. The Vitality Pro features 3 brushing modes:
  - **Daily Clean:** For everyday comprehensive cleaning.
  - **Sensitive:** Gentle cleaning for sensitive areas.
  - **Super Sensitive:** Extra gentle cleaning for very sensitive gums.

Select your desired mode by pressing the power button repeatedly.

5. Move the brush head slowly from tooth to tooth, allowing it to clean each surface. Do not scrub.
6. Brush your teeth for two minutes, utilizing the built-in timer. The timer signals every 30 seconds to indicate when to move to a different quadrant of your mouth, ensuring an even clean.
7. After brushing, press the power button to turn off the toothbrush.



Image: The Oral-B Vitality Pro electric toothbrush highlighting its three brushing modes (Daily Clean, Sensitive, Super Sensitive), the 4x30 second professional timer, and its long-lasting battery feature.

## 5. MAINTENANCE AND CLEANING

### 5.1 After Each Use

1. Rinse the brush head thoroughly under running water with the handle switched on.
2. Turn off the handle and remove the brush head.
3. Rinse both the brush head and the handle separately under running water.
4. Wipe the handle dry with a cloth.

### 5.2 Replacing the Brush Head

Dentists recommend replacing your brush head every 3 months, or sooner if the bristles become splayed or worn. The indicator bristles on the brush head will fade from green to yellow, signaling when it's time for a replacement.

### 5.3 Charger Cleaning

Unplug the charger before cleaning. Wipe the charger with a damp cloth. Do not immerse the charger in water.



Image: Important tips for responsible use and environmental care, including turning off the tap during brushing, unplugging the charger when not in use, the product's recyclable cardboard packaging, and its European origin.

## 6. TROUBLESHOOTING

Problem	Possible Cause	Solution
Toothbrush does not turn on.	Battery is discharged.	Place the toothbrush on the charger for at least 16 hours.
Toothbrush does not charge.	Charger not plugged in correctly or faulty outlet.	Ensure the charger is securely plugged into a working electrical outlet. Try a different outlet.
Reduced brushing performance.	Worn brush head or low battery.	Replace the brush head. Fully charge the toothbrush.
Bristles are splayed.	Brush head is old or too much pressure is applied.	Replace the brush head. Apply less pressure during brushing.

## 7. SPECIFICATIONS

---

**Model Number:** 427063

**Product Name:** Oral-B POC BAD Vitality Pro D103 Black HBOX 1CT

**Dimensions:** 11.3 x 8.4 x 22.4 cm

**Weight:** 246.65 grams

**Power Source:** Battery powered (Rechargeable)

**Battery Type:** 1 AAA battery (included) - *Note: This information appears to be a discrepancy in the source data, as this is a rechargeable electric toothbrush that comes with a charger. The device contains an internal rechargeable battery.*

**Special Features:** Pressure Sensor, Brushing Timer

**Age Range:** Adults

**Material:** Plastic

**Color:** Vitality Pro Protect X Clean Black

**Number of Items:** 1 (toothbrush handle, 1 brush head)

**Brush Head Shape:** Round

**Manufacturer:** Procter & Gamble

**Country of Origin:** Hungary

**What's in the Box:** Charger

## 8. WARRANTY AND SUPPORT

---

Information regarding product warranty and customer support was not available in the provided product data. Please refer to the product packaging or the official Oral-B website for detailed warranty terms and contact information for support.

For further assistance, you may visit the official Oral-B store: [Oral-B Store](#)