

[Manuals.plus](#) /

› [Oral-B](#) /

› Oral-B Vitality Electric Toothbrush User Manual

Oral-B NEW Model

Oral-B Vitality Electric Toothbrush User Manual

Model: NEW Model | Brand: Oral-B

INTRODUCTION

The Oral-B Vitality Electric Toothbrush is designed to provide a superior clean compared to a regular manual toothbrush. Its advanced 2D cleaning action oscillates and rotates to effectively remove plaque. Featuring a dentist-inspired round brush head and an in-handle timer, it helps ensure a thorough and consistent brushing experience.

SAFETY INFORMATION

Do not use this product on children under the age of 3.

Replace the brush head every 3 months or sooner if the brush head becomes worn.

WHAT'S INCLUDED

- One (1) Oral-B Vitality Electric Toothbrush
- Two (2) 3D White Replacement Brush Heads
- One (1) Charger



Image: The Oral-B Vitality Electric Toothbrush package, including the blue and white handle, charging base, and two brush heads.

SETUP

1. Remove the toothbrush from its packaging.
2. Charge the toothbrush before initial use. Place the toothbrush handle onto the charging base. The toothbrush can be kept on the charger when not in use to ensure it is always ready.
3. Attach a brush head by pushing it onto the handle until it clicks into place.

OPERATING INSTRUCTIONS

1. Wet the toothbrush head and apply toothpaste.
2. Place the toothbrush head in your mouth, against your teeth.
3. Press the power button to start brushing.
4. Brush for two minutes, following the 30-second quadrant timer. The timer will briefly pause or stutter every 30 seconds to signal you to move to the next quadrant of your mouth, ensuring a complete clean.
5. Rinse the toothbrush and your mouth after brushing.



Image: A person demonstrating the use of the Oral-B Vitality Electric Toothbrush.

Personalized Brush Modes



**DAILY
CLEAN**



SENSITIVE

Image: Illustration of the two personalized brush modes: Daily Clean and Sensitive.

Dentist-Inspired Round Brush Head



Image: Visual representation of the built-in 2-minute timer feature.

MAINTENANCE

Replacing the Brush Head

For optimal performance and hygiene, replace your Oral-B brush head every 3 months, or sooner if the bristles become splayed or worn. To replace, simply pull the old brush head off the handle and push a new one on until it clicks.

Replace Brush Head Every 3 Months



Image: Reminder to replace the brush head every 3 months for effective cleaning.

Cleaning

After each use, rinse the brush head thoroughly under running water with the handle switched on. Remove the brush head and clean the handle and brush head separately. Wipe the handle with a damp cloth. Do not immerse the charger in water.

TROUBLESHOOTING

Toothbrush Stuttering or Slowing Down

If your toothbrush appears to stutter or slow down during brushing, this is typically not a flaw. The Oral-B Vitality Electric Toothbrush features an in-handle timer that provides a brief pause or stutter every 30 seconds to indicate when to move to the next quadrant of your mouth. This ensures you brush evenly for the dentist-recommended 2 minutes.

Charging Issues

If the toothbrush is not charging, ensure the charging base is securely plugged into a working power outlet. Verify that the toothbrush handle is correctly seated on the charging base. If issues persist, contact customer support.

SPECIFICATIONS

Feature	Detail
Product Dimensions	4.33 x 2.13 x 8.27 inches; 10.93 ounces
Item Model Number	NEW Model
Manufacturer	Procter and Gamble
Brand	Oral-B
Age Range (Description)	Adult
Special Feature	Brushing Timer
Power Source	Battery Powered
Item Firmness Description	Soft

WARRANTY & SUPPORT

The Oral-B Vitality Electric Toothbrush is offered with a 30-day Money Back Guarantee. Please refer to the official Oral-B Website for full terms & conditions regarding warranty and support.



Image: Visual representation of the 30-day Money Back Guarantee.

For further assistance, visit the [Oral-B Store on Amazon](#) or contact Oral-B customer service.

© 2024 Oral-B. All rights reserved.