

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

› [Oral-B](#) /

› [Oral-B Precision Clean Electric Toothbrush Heads with CleanMaximiser Technology - User Manual](#)

Oral-B 4210201435389

Oral-B Precision Clean Electric Toothbrush Heads with CleanMaximiser Technology - User Manual

Brand: Oral-B | Model: 4210201435389

INTRODUCTION

This user manual provides comprehensive instructions for the Oral-B Precision Clean Electric Toothbrush Heads. Learn about the CleanMaximiser technology, proper installation, usage, and maintenance to achieve effective plaque removal and maintain optimal oral hygiene. These brush heads are compatible with most Oral-B electric toothbrushes, excluding Pulsonic and iO models.

Key Features

- **CleanMaximiser Technology:** Green bristles turn yellow to signal optimal time for replacement, ensuring better cleaning performance.
- **Superior Plaque Removal:** Designed to remove up to 100% more bacterial plaque in hard-to-reach places compared to a traditional manual toothbrush.
- **Precision Cleaning:** Features three different lengths of bristles for precise cleaning, wrapping each tooth individually.
- **Wide Compatibility:** Compatible with all Oral-B electric toothbrush handles, except for Pulsonic and iO models.
- **Convenient Packaging:** Supplied in recyclable mailbox packaging for easy delivery.

SETUP AND INSTALLATION

Follow these steps to properly attach your Oral-B Precision Clean brush head to your electric toothbrush handle:

1. Gently pull off the existing brush head from your Oral-B electric toothbrush handle.
2. Align the new Oral-B Precision Clean brush head with the metal shaft on the toothbrush handle.
3. Push the brush head down firmly until it clicks into place. Ensure it is securely fastened before use.



Image showing the Oral-B Precision Clean Electric Toothbrush Heads packaged for convenient delivery, designed to fit through a standard mailbox.

OPERATING INSTRUCTIONS

For optimal cleaning performance, follow these guidelines when using your Oral-B Precision Clean brush head:

1. Wet the brush head and apply a small amount of toothpaste.
2. Place the brush head on your teeth, ensuring it covers one tooth at a time.
3. Activate your Oral-B electric toothbrush.
4. Guide the brush head slowly from tooth to tooth, allowing the rotating and oscillating action to clean each tooth surface. Do not scrub aggressively.
5. Brush all surfaces of your teeth: outer, inner, and chewing surfaces.
6. Brush for the dentist-recommended two minutes, typically indicated by your toothbrush's built-in timer.



A close-up view of the Oral-B Precision Clean brush head actively cleaning a tooth surface, illustrating its effectiveness in removing plaque.



Diagram highlighting the Oral-B Precision Clean brush head's design, featuring three distinct bristle lengths for comprehensive cleaning and reaching difficult areas.

MAINTENANCE AND REPLACEMENT

Proper maintenance ensures the longevity and effectiveness of your brush heads.

- **CleanMaximiser Technology:** The green bristles on your Oral-B Precision Clean brush head are designed to fade and turn yellow over time. This color change is a visual indicator that the brush head is wearing out and signals the optimal time for replacement to maintain superior cleaning performance.
- **Replacement Frequency:** For best results and hygiene, replace your brush head every 3 months. You may need to replace it sooner if the bristles become frayed or spread apart.
- **Cleaning the Brush Head:** After each use, rinse the brush head thoroughly under running water to remove any toothpaste residue and debris. Allow it to air dry.

COMPATIBILITY

Oral-B Precision Clean brush heads are designed for broad compatibility with Oral-B electric toothbrush handles.

- These brush heads are compatible with **all Oral-B electric toothbrush handles**, except for Oral-B Pulsonic and Oral-B iO models.

- Always ensure your specific Oral-B handle is listed as compatible before purchasing replacement heads.



An illustration demonstrating the compatibility of Oral-B Precision Clean brush heads with various Oral-B electric toothbrush handles, explicitly noting exclusion of Pulsonic and iO models.

SPECIFICATIONS

Feature	Detail
Brand Name	Oral-B
Model Name	Oral-B Precision Clean
Item Model Number	4210201435389
Number of Pieces	16
Unit Count	16.0 Count
Material Type	Cardboard (packaging)
Bristle Material	Nylon

Feature	Detail
Item Firmness Description	Medium
Color	White
Age Range Description	Adult
Manufacturer	Procter & Gamble
Package Dimensions	10.51 x 8.15 x 0.87 inches
Item Weight	111.2 Grams (3.92 ounces)

SAFETY INFORMATION AND WARNINGS

Please read and adhere to the following safety guidelines:

- For adult use only.
- Keep out of reach of children.
- Do not use if you have an allergy to nylon, as the bristles are made of nylon.
- Always follow the instructions provided with the product and your Oral-B electric toothbrush handle.
- Consult your dentist or physician if you experience any discomfort or adverse reactions during use.