

LifeBasis LifeBasis Electric Toothbrush

LifeBasis Electric Toothbrush User Manual

Model: LifeBasis Electric Toothbrush

1. INTRODUCTION

This manual provides essential information for the safe and effective use of your LifeBasis Electric Toothbrush. Please read these instructions thoroughly before first use and retain for future reference.

2. SAFETY INFORMATION

- Brush gently to prevent damage to teeth and gums.
- Do not rinse the main body of the toothbrush with water or immerse it.
- Ensure the USB interface remains dry at all times.

3. PRODUCT FEATURES

3.1 High-Frequency Vibration



The toothbrush operates at 2300 vibrations per minute for effective cleaning.

3.2 Interchangeable Brush Heads

The kit includes four multifunctional brush heads designed for various cleaning needs:

- **Silicone Brush:** For sturdy dental plaque.
- **Flat Brush:** For food residue on the tooth surface.
- **Sharp Brush:** Designed to clean tooth tartar between the teeth.
- **Silicone Cup Brush:** For polishing teeth and gum care.

Calculus Scraper



Silicone Brush
sturdy dental plaque



Flat Brush
food residue on tooth surface



Sharp Brush
tooth tartar between the teeth



Silicone Cup Brush
polish teeth and care gum



Overview of the four included brush head types.

Sharp Brush



The sharp brush head is designed to reach between teeth.

3.3 Cleaning Modes

The toothbrush features three distinct cleaning modes, selectable via a single button:

- **Clean Mode:** Standard cleaning.
- **Medium Mode:** Moderate intensity cleaning.
- **Sensitive Mode:** Gentle cleaning for sensitive teeth and gums.

Three tooth brushing mode

One Button

Clean mode

Medium mode

Sensitive mode



Select from three cleaning modes using the single control button.

3.4 Power Supply Options

The device can be powered in two ways:

- Via the included universal micro USB cable and an adapter plug (for home use, keep USB plug in when in use).
- Using 2 AA batteries (not included).

2 POWER SUPPLY WAYS

1. universal micro USB cable (included) during home use
(Please keep the USB plug in when use)
2. or by 2 AA batteries (excluded)



The toothbrush supports both USB and AA battery power.

4. SETUP

1. **Install Brush Head:** Securely attach the desired brush head onto the machine.
2. **Powering the Device:**
 - For USB power: Connect the provided USB cable to the toothbrush and an adapter plug. Ensure the USB interface is dry before connecting.
 - For battery power: Insert 2 AA batteries into the battery compartment.

5. OPERATION

1. **Apply Toothpaste:** Apply a small amount of toothpaste, polishing paste, or tooth powder onto the brush head.
2. **Brushing Technique:** Touch the cleaning head to your teeth. Move it slowly back and forth, staying on each tooth for a

few seconds. Do not press hard to avoid damage to teeth or gums.

3. **Selecting Mode:** Use the single button on the handle to cycle through the Clean, Medium, and Sensitive modes.

5.1 Instructional Video

Your browser does not support the video tag.

Watch this video for a visual guide on how to use the LifeBasis Electric Toothbrush.

6. MAINTENANCE

1. **After Each Use:** Remove the brush head for cleaning and drying. Place it back into the storage box once dry.
2. **Water Resistance:** The product has an IPX3 waterproof rating, meaning it is protected against spraying water. However, the bottom USB interface must not come into contact with water. Do not rinse the entire toothbrush with water after brushing.
3. **Brush Head Replacement:** Replace the brush head regularly for optimal hygiene and performance.

7. TROUBLESHOOTING

- **Device Not Turning On:** Ensure the USB cable is properly connected and the adapter plug is receiving power, or check if the AA batteries are correctly inserted and have sufficient charge.
- **Reduced Performance:** Replace the brush head if it appears worn. Ensure the power source is stable.
- **Difficulty Changing Modes:** Press the control button firmly. If issues persist, ensure the device is adequately powered.

8. SPECIFICATIONS

Feature	Detail
Model Name	LifeBasis Electric Toothbrush
Product Dimensions	1.3 x 1.3 x 8 inches
Weight	5.93 ounces
Power Source	Battery Powered (USB cable/adapter plug or 2 AA batteries)
Number of Cleaning Modes	3
Operation Mode	Automatic
Bristle Material	Silicone
Included Components	4 Multifunctional Brush Heads, USB cable, adapter plug
Water Resistance	IPX3

9. WARRANTY AND SUPPORT

For warranty information and customer support, please refer to the documentation included with your purchase or visit the official LifeBasis store online. Continuous use for at least 3 months may be required to observe obvious results.

Legal Disclaimer: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question