



Manuals.plus /

› PHILIPS /

› PHILIPS Sonicare Diamond Clean Rechargeable Toothbrush User Manual

PHILIPS Sonicare Diamond Clean

PHILIPS Sonicare Diamond Clean Rechargeable Toothbrush User Manual

Model: Sonicare Diamond Clean | Brand: PHILIPS

INTRODUCTION

This manual provides essential information for the safe and effective operation, maintenance, and troubleshooting of your PHILIPS Sonicare Diamond Clean Rechargeable Toothbrush. Please read it thoroughly before first use and keep it for future reference.

PACKAGE CONTENTS

Your PHILIPS Sonicare Diamond Clean Rechargeable Toothbrush package includes the following items:

- Two (2) Sonicare Handles (One White, One Black)
- Two (2) Wireless Charging Glasses
- Two (2) Premium Plaque Brush Heads
- Two (2) USB Charging Travel Cases



Image: Contents of the PHILIPS Sonicare Diamond Clean 2-Pack, showing the black and white toothbrush handles, charging glasses, and travel cases.

KEY FEATURES

- **Automatic Operation:** The toothbrush operates in automatic mode for convenience and efficiency.
- **5 Different Settings:** The toothbrush offers 5 different cleaning modes for tailored oral care.
- **Diamond Shaped Bristles:** The toothbrush features diamond shaped bristles for effective cleaning.
- **Rechargeable Battery:** The toothbrush is powered by a rechargeable battery for long-lasting use.
- **Complete Oral Care:** The toothbrush provides complete oral care with its various features and settings.

SETUP GUIDE

1. **Initial Charge:** Before first use, fully charge your Sonicare toothbrush for at least 24 hours. Place the handle into the wireless charging glass, ensuring it is upright and stable. The charging indicator light will illuminate to confirm charging.
2. **Attach Brush Head:** Align the brush head with the metal shaft on the toothbrush handle and push down firmly until it clicks into place.
3. **First Use:** Moisten the bristles and apply a small amount of toothpaste. Place the brush head against your teeth at a slight angle towards the gumline.

OPERATING INSTRUCTIONS

Your Sonicare Diamond Clean toothbrush offers 5 distinct cleaning modes to suit your oral care needs:

- **Clean:** Standard mode for superior everyday cleaning.

- **White:** Removes surface stains for a brighter smile.
- **Sensitive:** Gentle cleaning for sensitive teeth and gums.
- **Gum Care:** Gently massages gums for improved gum health.
- **Deep Clean:** Provides an invigorating deep clean.

To operate:

1. Press the power button once to turn on the toothbrush. It will start in the default "Clean" mode.
2. To change modes, press the power button repeatedly within 2 seconds of turning it on. The mode indicator light will show the selected mode.
3. To turn off the toothbrush, press and hold the power button for 2 seconds, or press the power button once if the toothbrush has been running for more than 2 seconds in a single mode.
4. The toothbrush features a 2-minute Smartimer to ensure recommended brushing time and a QuadPacer that signals when to move to the next quadrant of your mouth.

MAINTENANCE

- **After Each Use:** Rinse the brush head and bristles thoroughly under running water. Wipe the handle with a damp cloth.
- **Weekly Cleaning:** Remove the brush head and rinse the metal shaft area with warm water. Clean the charging glass and travel case with a damp cloth.
- **Brush Head Replacement:** For optimal performance and hygiene, replace your brush head every 3 months, or sooner if the bristles are splayed or worn.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Toothbrush does not turn on.	Battery is depleted.	Place the toothbrush on the charger for at least 24 hours.
Weak vibration or reduced power.	Low battery or brush head is worn.	Recharge the toothbrush. Replace the brush head if it's old.
Charging indicator not lighting up.	Charger not connected or faulty.	Ensure the charging glass/travel case is properly plugged in. Try a different outlet.

SPECIFICATIONS

Feature	Detail
Brand	PHILIPS
Power Source	Battery Powered (Rechargeable)
Age Range (Description)	Adult
Item Firmness Description	Medium
Package Dimensions	15.31 x 14.21 x 3.82 inches
Weight	4.67 Pounds
Manufacturer	Philips
ASIN	B0BF8PLVNP

SAFETY INFORMATION

Important: Keep out of reach of children. Do not use for purposes other than oral hygiene.

LEGAL DISCLAIMER

This product is intended for personal oral hygiene use. Any statements regarding health benefits are for informational purposes only and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

WARRANTY AND SUPPORT

For warranty information, product registration, or technical support, please visit the official PHILIPS website or contact their customer service. Retain your proof of purchase for warranty claims.

You can find more information and support at: [PHILIPS Store on Amazon](#)