Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

[manuals.plus](#) /

› [COMFIER](#) /

› [COMFIER 2-in-1 Shiatsu Foot Massager&Foot Rest,FSA HSA Eligible Foot and Calf Massager with APP Control,Heated Feet Massager for Foot Pain Relief,Gifts for Men Women Off-white](#)

COMFIER CF-5001S-APP

COMFIER 2-in-1 Shiatsu Foot & Calf Massager User Manual

Model: CF-5001S-APP

IMPORTANT SAFETY INFORMATION

Please read all instructions carefully before using this appliance. Failure to follow safety instructions may result in electric shock, fire, or serious injury.

- Consult a physician before use if you have a medical condition, are pregnant, or have reduced sensation.
- Do not use the massager while sleeping or driving.
- Stop use immediately if you experience discomfort or pain.
- Keep the massager away from water and ensure it is properly grounded for electrical safety.
- Do not use the massager on bare skin or with oils or lotions.
- Do not let children play with the massager.

PRODUCT OVERVIEW

The COMFIER 2-in-1 Shiatsu Foot & Calf Massager is a versatile device designed to provide comprehensive massage therapy for your feet and calves, while also functioning as a comfortable ottoman.

Figure 1: COMFIER 2-in-1 Shiatsu Foot & Calf Massager in its primary configuration, showing the foot and calf massage areas and APP control interface.

Key Features:

- **2-in-1 Functionality:** Easily converts between a foot and calf massager and a luxury ottoman foot rest.
- **Multiple Massage Functions:** Offers Shiatsu, compression, rolling massage, and heating functions for comprehensive relief.
- **Calf Massage:** Provides relaxing compression and vibration massage to alleviate leg tension and sore calf muscles.
- **Customizable Settings:** Features 2 modes and 4 compression pressure intensities for personalized massage experiences.
- **Intelligent APP Control:** Conveniently control settings via a mobile application, in addition to direct button control.
- **Adjustable Support Bar:** Allows for a comfortable 20-degree inclined massage angle.
- **Easy to Clean:** Zippered foot sleeves are detachable for hygienic cleaning.
- **Wide Foot Size Compatibility:** Accommodates most foot sizes up to men's size 13.

SETUP

1. **Unboxing:** Carefully remove the massager from its packaging. Ensure all components, including the power adapter and user manual, are present.
2. **Placement:** Place the massager on a flat, stable surface near a power outlet.
3. **Power Connection:** Plug the power adapter into the massager's input port and then into a wall outlet.
4. **Adjusting the Support Bar:** If using as a massager, pull up the support bar located at the bottom of the unit to achieve a 20-degree inclined angle for optimal comfort.

Figure 2: The adjustable support bar allows for a comfortable massage angle.

5. **APP Installation (Optional):** For intelligent control, download the "Comfier" APP from your mobile device's app store. Follow the in-app instructions to register and connect your device via Bluetooth.

OPERATING INSTRUCTIONS

The massager can be operated using the control panel buttons or via the mobile APP.

Control Panel Functions:

Figure 3: Control panel with various massage settings.

- **Power Button:** Press to turn the massager ON/OFF.
- **Foot Shiatsu Button (Paw Icon):** Activates/deactivates the Shiatsu and rolling massage for feet.
- **Leg Mode Button (Boot Icon):** Activates/deactivates compression and vibration massage for calves.
- **Heat Button (Wave Icon):** Toggles heat function ON/OFF for both feet and calves.
- **Intensity Button (Bars Icon):** Cycles through 4 levels of compression intensity (both foot and calf).

Using the Massager:

1. Sit comfortably in a chair with the massager positioned in front of you.

2. Insert your feet into the foot chambers. For calf massage, ensure your calves are positioned correctly within the upper sections.
3. Press the Power Button to turn on the massager.
4. Select your desired massage mode (Foot Shiatsu, Leg Mode) and adjust the intensity and heat settings using the control panel buttons or the APP.
5. The massager has an automatic shut-off feature after 15 minutes. You can adjust this timer via the APP (from 1 to 30 minutes).

MAINTENANCE

Regular maintenance ensures the longevity and hygienic use of your massager.

- **Cleaning:** The zippered foot sleeves can be easily detached for cleaning. Hand wash or machine wash on a gentle cycle with mild detergent. Air dry completely before reattaching.
- **Storage:** When not in use, fold the massager back into its ottoman form and store it in a cool, dry place away from direct sunlight.
- **General Care:** Wipe the exterior of the massager with a soft, damp cloth. Do not use abrasive cleaners or solvents.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Massager does not turn on.	Not properly plugged in or power outlet issue.	Ensure power adapter is securely connected and try a different outlet.
Weak or no massage intensity.	Intensity setting is too low or feet/calves not positioned correctly.	Increase intensity level. Reposition feet/calves to ensure full contact with massage nodes.
Heat function not working.	Heat function is not activated.	Press the Heat Button to activate. Allow a few minutes for heat to build up.
APP not connecting to massager.	Bluetooth is off or massager is out of range.	Ensure Bluetooth is enabled on your device and the massager is within range. Restart both devices if necessary.

SPECIFICATIONS

- **Model Number:** CF-5001S-APP

- **Product Dimensions:** 19 x 16.9 x 18 inches
- **Item Weight:** 22.6 Pounds
- **Material:** Leather
- **Power Source:** Corded Electric
- **Max Foot Size:** Men's US Size 13
- **Automatic Shut-off:** 15 minutes (adjustable via APP)

WARRANTY & SUPPORT

COMFIER products come with a manufacturer's warranty. Please refer to the warranty card included in your packaging for specific terms and conditions.

For any product inquiries, technical support, or warranty claims, please contact COMFIER Customer Service:

- **Website:** [Visit the COMFIER Store on Amazon](#)
- **Email:** Refer to user manual for specific email address.
- **Phone:** Refer to user manual for specific phone number.

If the electric foot massager does not meet your expectations for any reason, you may return it for a refund within 30 days of purchase.

© 2023 COMFIER. All rights reserved.

Related Documents - CF-5001S-APP

 <p>COMFIER FUZZY MASSAGE BOOT (F) FOOT MASSAGER (H) MASSAGEUR DE PIEDS (F) MASSAGGIATORE PIEDI (I) MASSADOR DE PIES (E)</p> <p>CE User Manual CF-5913</p>	<p>Comfier CF-5913 Shiatsu Foot Massager User Manual</p> <p>User manual for the Comfier CF-5913 Shiatsu Foot Massager Machine. Learn how to use this kneading foot massager with heat for pain relief, plantar fasciitis, and general relaxation. Includes operation, safety, maintenance, and warranty details.</p>
 <p>COMFIER Leg and Foot Massager</p> <p>CE User Manual CF-F25H</p>	<p>Comfier CF-F25H Leg and Foot Massager User Manual</p> <p>Comprehensive user manual for the Comfier CF-F25H Leg and Foot Massager, detailing safety instructions, product features, setup and operation, maintenance, troubleshooting, warranty information, and FCC compliance.</p>

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