

Tilcare 1

Tilcare Portable Urinal for Men with Spill-Proof Glow-in-the-Dark Lid

Model: 1

[Introduction](#)

[Features](#)

[Setup](#)

[Operation](#)

[Maintenance](#)

[Troubleshooting](#)

[Specifications](#)

[Warranty & Support](#)

1. INTRODUCTION

The Tilcare Portable Urinal for Men is designed for convenient and hygienic use in various situations, including for bedridden individuals, during travel, or in outdoor settings. This manual provides instructions for proper use, care, and maintenance to ensure optimal performance and longevity of your urinal.

2. KEY FEATURES

2.1 Large Capacity

The urinal features a generous 3000 ml (102 fl oz) capacity, reducing the need for frequent emptying. This makes it suitable for extended use, such as overnight or during long journeys.

102 fl oz (3000ml) size

Built for long nights and extended trips



Figure 2.1: Urinal with 3000 ml capacity, suitable for extended use.

2.2 Leak-Proof Screw Lid

Equipped with a secure screw-on lid, the urinal is designed to prevent spills and leaks, even when tilted or during transport. This feature helps maintain hygiene and cleanliness.

100% LEAK-PROOF SCREW LID

No spills, no mess— peace of mind assured



Figure 2.2: Secure screw lid for spill prevention.

2.3 Glow-in-the-Dark Lid

The lid features a glow-in-the-dark property, enhancing visibility in low-light conditions. This allows for easier location and use in dark environments such as bedrooms, cars, or campsites.

GLOW IN THE DARK LID

Easy to find in those hard to see places



Figure 2.3: Glow-in-the-dark lid for nighttime visibility.

2.4 Portable and Versatile Design

Its lightweight and travel-friendly design makes it suitable for a wide range of users and scenarios, including bedridden patients, truckers, campers, and hospital use. The wide opening facilitates ease of access.

VERSATILE, TRAVEL-FRIENDLY & PORTABLE

Care that goes wherever you go



Figure 2.4: Versatile use in various environments.

2.5 Heavy-Duty Construction

Constructed from thick, durable plastic, the urinal is built to withstand regular use. Its anti-tip-over design enhances stability, minimizing the risk of spills on uneven surfaces.

HEAVY-DUTY CONSTRUCTION

Reliable and Long-Lasting



Generic Product

Thin & Fragile Material
Flimsy Design
Prone to Quick Wear
Leaks Under High Pressure
Breaks Easily with Accidental Drops



Tilcare Urinal

Thick & Durable Plastic
Durable Design
Resistant to Wear and Tear
Handles High Pressure Easily
Withstands Accidental Drops

Figure 2.5: Durable construction for reliability.

2.6 Wide and Comfortable Opening

The urinal features a wide opening designed for ease of use and comfort, ensuring a hassle-free experience.

WIDE & COMFORTABLE OPENING

Designed for ease—use with no hassle



Figure 2.6: Wide opening for comfortable use.

2.7 Easy Pour Handle

An integrated handle allows for effortless and clean disposal of contents, helping to maintain hygiene.

EASY POUR HANDLE

Effortless disposal keeping your hands clean



Figure 2.7: Handle for easy and hygienic disposal.

3. SETUP

1. **Unpack:** Carefully remove the urinal from its packaging.
2. **Inspect:** Check the urinal for any signs of damage or defects. Ensure the screw lid is present and functions correctly.
3. **Initial Cleaning:** Before first use, wash the urinal thoroughly with warm water and mild soap. Rinse completely and allow to air dry.
4. **Secure Lid:** Ensure the screw lid is securely fastened when not in use to prevent spills.

4. OPERATING INSTRUCTIONS

1. **Preparation:** Unscrew the lid and place it aside in a clean area.
2. **Positioning:** Position the urinal comfortably for use. The wide opening is designed for ease of access.

3. **Use:** Use the urinal as needed. The 3000 ml capacity allows for multiple uses before emptying.
4. **Securing After Use:** After use, carefully replace the screw lid and turn clockwise until it is tightly sealed. This prevents leaks and odors.
5. **Disposal:** When ready to empty, hold the urinal by the integrated handle and unscrew the lid. Pour the contents into a toilet or designated disposal area.

5. MAINTENANCE AND CLEANING

Regular cleaning is essential for hygiene and to prolong the life of your urinal.

- **Daily Cleaning:** After each use, or at least once daily, empty the urinal and rinse it thoroughly with warm water.
- **Deep Cleaning:** For a more thorough cleaning, use warm water and a mild disinfectant soap. Use a brush if necessary to reach all internal surfaces.
- **Rinsing:** Rinse the urinal completely with clean water to remove all soap residue.
- **Drying:** Allow the urinal to air dry completely before re-securing the lid. This helps prevent bacterial growth and odors.
- **Storage:** Store the clean and dry urinal in a cool, dry place away from direct sunlight.

6. TROUBLESHOOTING

Problem	Possible Cause	Solution
Leakage from lid	Lid not screwed on tightly; lid or urinal threads damaged.	Ensure the lid is screwed on firmly until hand-tight. Inspect lid and urinal threads for damage; replace if necessary.
Difficulty cleaning thoroughly	Residue buildup in corners or hard-to-reach areas.	Use a long-handled brush specifically designed for cleaning bottles or containers. Ensure thorough rinsing after cleaning.
Glow-in-the-dark lid not glowing	Lid has not been exposed to sufficient light to charge.	Expose the lid to bright light (natural or artificial) for at least 15-30 minutes to allow it to charge before use in darkness.

7. SPECIFICATIONS

Attribute	Detail
Brand	Tilcare
Model Number	1
Capacity	3000 ml (102 fl oz)
Item Weight	9.9 ounces
Package Dimensions	11.85 x 6.89 x 4.37 inches
Color	White
Material	Durable, thick plastic
Special Features	Spill-proof screw lid, Glow-in-the-dark lid, Anti-tip stability, Easy pour handle

8. WARRANTY AND SUPPORT

For any questions regarding your Tilcare Portable Urinal, including warranty information or technical support, please contact Tilcare customer service. Refer to the product packaging or the official Tilcare website for contact details.

© 2026 Tilcare. All rights reserved.