

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

> [Braun](#) /

> [Braun Oral-B iO Ultimate Clean Replacement Brush Head iORBCB-4EL User Manual](#)

Braun iORBCB-4EL

Braun Oral-B iO Ultimate Clean Replacement Brush Head User Manual

Model: iORBCB-4EL | Brand: Braun

INTRODUCTION

This manual provides essential instructions for the proper use and maintenance of your Braun Oral-B iO Ultimate Clean Replacement Brush Head (iORBCB-4EL). This brush head is designed for use with Oral-B iO Series electric toothbrushes, offering superior plaque removal for a comprehensive clean.

KEY FEATURES


- **Ultimate Clean Design:** Features two types of bristle bundles with different lengths to thoroughly remove plaque.
- **Angled Bristles:** Oral-B's unique 16-degree angled bristles are engineered to reach and remove plaque from hard-to-reach areas.
- **iO Series Compatibility:** Specifically designed for Oral-B iO Series electric toothbrushes, ensuring optimal performance.
- **Round Rotation Technology:** Utilizes Oral-B's unique 'round rotation' brush action, inspired by professional dental cleaning tools, to transfer energy efficiently from the iO body to the brush tip.

SETUP AND INSTALLATION

Follow these steps to correctly attach the replacement brush head to your Oral-B iO electric toothbrush handle:

1. Align the replacement brush head with the top of the toothbrush handle.
2. Gently push the brush head down onto the handle until you hear a click, indicating it is securely in place.

16度に傾斜したブラシが 届きにくい歯間まで歯垢を除去



Dentists recommend replacing the brush head every 3 months to maintain optimal cleaning performance. This image illustrates the process of attaching a new brush head to the Oral-B iO handle.

OPERATING INSTRUCTIONS

Once the brush head is attached, you can begin brushing:

1. Apply a small amount of toothpaste to the brush head.
2. Place the brush head against your teeth at a slight angle towards the gumline.
3. Turn on the toothbrush and guide it slowly from tooth to tooth, allowing the brush to do the work. Do not scrub aggressively.
4. Brush for the dentist-recommended two minutes, ensuring all areas of your mouth are cleaned.

99.7%UPの 歯垢除去力で 歯垢の蓄積まで防ぐ



*手磨きとの比較。P&G調べ

The Oral-B iO brush head is designed to remove plaque efficiently, preventing accumulation. This image demonstrates the brush head's action on tooth surfaces.

MAINTENANCE AND STORAGE

Proper maintenance ensures the longevity and effectiveness of your brush head:

- **Replacement Frequency:** Replace the brush head every 3 months, or sooner if the bristles become frayed or worn. This restores 100% of the dental plaque removal power.
- **Cleaning:** After each use, rinse the brush head thoroughly under running water.
- **Storage:** Store the brush head in a dry place to maintain its effectiveness and hygiene.

長さの違う2つの毛束で すみずみまで磨き上げる



Replacing your brush head restores 100% of the plaque removal power, ensuring effective cleaning. This image visually compares a used brush head with a new one.


COMPATIBILITY

The Oral-B iO Ultimate Clean Replacement Brush Head (iORBCB-4EL) is compatible with all Oral-B iO Series electric toothbrushes, including:

- Oral-B iO3
- Oral-B iO4
- Oral-B iO5
- Oral-B iO6
- Oral-B iO7
- Oral-B iO8
- Oral-B iO9
- Oral-B iO10

替えブラシの交換で歯垢除去力100%復活

100%
洗浄力



Oral-B iO offers various brush head types to cater to different dental needs, including plaque removal, gum gentleness, and whitening. This chart provides a comparison of different brush head functionalities.

SPECIFICATIONS

Feature	Detail
Model Name	iORBCB-4EL
Brand	Braun
Colour	Black
Size	Standard
Material	Nylon
Age Range (Description)	Adult
Bristle Material	Nylon
Item Firmness Description	Medium

Number of Pieces	4 Count
Compatible Devices	Oral-B iO Series Electric Toothbrushes

IMPORTANT SAFETY INFORMATION

Warning:

To ensure optimal performance and safety, always use genuine Oral-B replacement brush heads. Product failure or damage resulting from the use of compatible products or parallel import goods from overseas is not covered by the manufacturer's warranty. Using non-genuine parts may also affect the effectiveness of your toothbrush and could potentially cause damage.

WARRANTY AND SUPPORT

For information regarding product warranty, returns policy, or customer support, please refer to the documentation provided with your original Oral-B iO electric toothbrush or visit the official Braun website. Keep your proof of purchase for any warranty claims.