

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

› [RENPHO](#) /

› [RENPHO R-HAP61 Weighted Electric Heating Pad Instruction Manual](#)

## RENPHO R-HAP61

# RENPHO R-HAP61 Weighted Electric Heating Pad

## INSTRUCTION MANUAL

### INTRODUCTION

---

Thank you for choosing the RENPHO R-HAP61 Weighted Electric Heating Pad. This instruction manual provides essential information for the safe and effective use of your new heating pad. Designed for targeted relief, this large 24x33 inch heating pad covers your back, neck, and shoulders, offering soothing heat therapy. Please read this manual thoroughly before first use and retain it for future reference.

### IMPORTANT SAFETY INFORMATION

---

To ensure safe operation and prevent injury or damage to the product, please adhere to the following guidelines:

- This product contains built-in heating wires. Do not bunch, sit on, lay on, or squeeze the pad. Avoid sharp creases, which may cause short circuits and burns.
- Do not use the heating pad directly on bare skin. Always wear a piece of clothing underneath.
- Do not use while sleeping.
- The applicable voltage for this product is 120V and the rated power is 135W.
- The heating pad is ETL Certified and features built-in overheat protection, which automatically shuts off the device if it becomes too hot.
- Keep the controller dry and away from water.

### PRODUCT FEATURES

---

- **Large Coverage Area:** Measures 24x33 inches, designed to cover the full back, neck, and shoulders simultaneously.
- **Weighted Design:** Features weighted edges and a weighted neck collar to ensure the pad stays securely in place and provides gentle pressure for enhanced comfort.
- **Soft Material:** Made from super soft flannel for maximum comfort against the skin.

- **Adjustable Heat Settings:** Offers 6 heat settings to customize your heat therapy experience.
- **Timer Settings:** Includes 3 timing options (30, 60, and 90 minutes) with automatic shut-off for safety and convenience.
- **Digital Controller:** Easy-to-read LCD controller for precise temperature and timer adjustments.
- **Long Power Cord:** A 7-foot-long power cord provides flexibility and allows for movement during use.
- **Machine Washable:** The pad is machine-washable after removing the controller, simplifying cleaning.

## SETUP

1. **Unpack:** Carefully remove the heating pad and controller from its packaging.
2. **Connect Controller:** Ensure the controller cable is securely plugged into the heating pad's connector.
3. **Position the Pad:** Place the heating pad over your neck and shoulders, allowing the main body of the pad to cover your back. Use the integrated ties to secure the pad around your waist for a snug fit.
4. **Plug In:** Insert the power plug into a standard 120V electrical outlet.



Figure 1: RENPHO Weighted Electric Heating Pad, Gray Green. This image shows the heating pad laid flat, highlighting

its large size and shape designed to cover the back, neck, and shoulders. The integrated ties are visible at the waist area.

## OPERATING INSTRUCTIONS

---

The digital controller allows for easy adjustment of heat and timer settings.

1. **Power On/Off:** Press the power button (**U**) to turn the heating pad on or off. The LCD screen will illuminate.
2. **Adjust Heat Settings:** Use the '+' and '-' buttons to cycle through the 6 available heat settings, ranging from 100°F to 150°F. The selected temperature will be displayed on the screen.
3. **Set Timer:** Press the timer button (**Clock icon**) to select a duration of 30, 60, or 90 minutes. The timer will count down, and the pad will automatically shut off once the time expires. If no timer is set, the pad will auto-shutoff after 30 minutes for safety.
4. **Enjoy Heat Therapy:** Allow the pad to warm up. It heats up quickly, typically within 5 seconds.

# Large Size & Full Body Relief



Warms up instantly | Pain relief effectively



Figure 2: Close-up of the RENPHO heating pad's digital controller, showing the LCD display with temperature and timer settings, along with power, temperature adjustment (+/-), and timer buttons.

## Relieve Pain & Enjoy Warmth



Figure 3: An illustration demonstrating the 6 heat settings (from 100°F to 150°F) and 3 timer settings (30, 60, 90 minutes) available on the RENPHO heating pad's controller.

## MAINTENANCE

---

### Cleaning Instructions

The RENPHO heating pad is designed for easy cleaning:

1. **Disconnect:** Always disconnect the controller from the heating pad and unplug the power cord from the outlet before cleaning.
2. **Machine Wash:** The heating pad is machine-washable. Use a gentle cycle with cold water.
3. **Drying:** Do not tumble dry or iron at high temperatures. Air dry the pad completely before reconnecting the controller and using it again.
4. **Frequency:** Since the heating pad has built-in heating wires, avoid frequent machine washing to prolong its lifespan.

## Easy To Clean

The heating pad can hand & machine wash  
please make sure disconnect the plug before washing



Figure 4: A woman is shown placing the disconnected RENPHO heating pad into a washing machine, illustrating the ease of cleaning.

### TROUBLESHOOTING

---

- **Display shows 'F2':** This indicates that the cable of the hand controller is disconnected from the heating pad. To resolve this, simply plug the cable securely between the heating pad and the hand controller.

### SPECIFICATIONS

---

Feature	Detail
Model Number	R-HAP61
Dimensions	24" x 33"
Material Type	Polyester (Super Soft Flannel)
Color	Gray Green
Heat Settings	6 (100°F - 150°F)
Timer Settings	3 (30, 60, 90 minutes) with auto-off
Rated Voltage	AC 120V
Rated Power	135W
Item Weight	2.9 Pounds
Manufacturer	RENPHO
UPC	840191403569

## WARRANTY AND SUPPORT

---

For warranty information, product support, or any questions regarding your RENPHO R-HAP61 Weighted Electric Heating Pad, please refer to the contact information provided in your product packaging or visit the official RENPHO website. Keep your purchase receipt for warranty claims.