



Manuals.plus /

- › Oral-B /
- › Oral-B iO3 Electric Toothbrush User Manual

Oral-B iO G3.1A6.0

Oral-B iO3 Electric Toothbrush User Manual

Model: iO G3.1A6.0

INTRODUCTION

The Oral-B iO Series 3 Electric Toothbrush combines advanced technology with a dentist-inspired round brush head and micro-vibrating bristles to deliver a professional purifying clean. This rechargeable toothbrush is clinically proven to remove 400% more plaque along the gumline compared to a regular manual toothbrush. It features a built-in 2-minute timer with 30-second interval alerts for thorough brushing, and a pressure sensor to indicate optimal cleaning pressure.



Figure 1: Oral-B iO3 Electric Toothbrush with its packaging and charging base.

KEY FEATURES

- **Advanced Plaque Removal:** Clinically proven to remove 400% more plaque along the gumline compared to a regular manual toothbrush.
- **Smart Pressure Sensor:** Indicates optimal cleaning pressure and warns against excessive pressure to protect gums.
- **2-Minute Timer with Quadrant Pacer:** Ensures a dentist-recommended brushing time with 30-second intervals for each quadrant of the mouth.
- **Three Cleaning Modes:** Offers Daily Clean, Sensitive, and Whitening modes for personalized brushing.
- **Round Brush Head Design:** Reaches areas that rectangular manual brushes may miss, mimicking the natural shape of teeth for thorough cleaning.
- **Long-Lasting Lithium-Ion Battery:** Provides extended use between charges.

Remove 100% more plaque in 1 week*

with Oral-B's most advanced technology



Figure 2: Visual comparison illustrating enhanced plaque removal with the Oral-B iO electric toothbrush.

Protect your gums



Figure 3: The smart pressure sensor illuminates green for optimal pressure, ensuring effective and safe cleaning.

Smart features to maximize your clean



Figure 4: Smart features of the Oral-B iO toothbrush, including the 2-minute timer, multiple cleaning modes, and brush head replacement alert.

SETUP

1. **Unpack:** Carefully remove all components from the packaging: the Oral-B iO3 toothbrush handle, one Ultimate Clean brush head, and the charger.
2. **Charge:** Place the toothbrush handle onto the charging base. Ensure the charger is plugged into a suitable power outlet. Allow the toothbrush to fully charge before first use. The battery indicator on the handle will show charging status.
3. **Attach Brush Head:** Align the brush head with the shaft of the toothbrush handle and push it down until it clicks into place.

OPERATING INSTRUCTIONS

1. **Wet Brush Head:** Lightly wet the brush head.
2. **Apply Toothpaste:** Apply a small amount of toothpaste to the brush head.
3. **Power On:** Press the power button on the handle to turn on the toothbrush.
4. **Select Mode:** Use the mode selection button (if available on your model) to cycle through the three cleaning modes:

- **Daily Clean:** For everyday thorough cleaning.
 - **Sensitive:** For gentle yet effective cleaning of sensitive areas.
 - **Whitening:** Designed to polish and whiten teeth.
5. **Brush Teeth:** Guide the brush head slowly from tooth to tooth, allowing the brush to do the work. Do not scrub. The round brush head is designed to clean each tooth individually.
 6. **Pressure Sensor:** The light ring around the brush head will illuminate to provide real-time pressure guidance.
 - **Green Light:** Indicates optimal brushing pressure.
 - **Red Light:** Indicates too much pressure. Reduce pressure to avoid harming gums and teeth.
 7. **Timer Guidance:** The built-in timer will pulse every 30 seconds to indicate when to move to the next quadrant of your mouth. After 2 minutes, the toothbrush will signal completion of the recommended brushing time.
 8. **Power Off:** Press the power button to turn off the toothbrush after brushing.

MAINTENANCE

- **Clean After Use:** After each use, rinse the brush head thoroughly under running water with the toothbrush still switched on. Remove the brush head and rinse the handle and the brush head separately. Wipe the handle dry.
- **Brush Head Replacement:** Replace your brush head every 3 months, or sooner if the brush head becomes worn or the indicator bristles fade. The toothbrush may also provide an alert when it's time to replace the brush head.
- **Charging Base Care:** Clean the charging base regularly with a damp cloth. Ensure it is unplugged before cleaning.
- **Storage:** Store the toothbrush and brush heads in a clean, dry place. The included upright case can be used to hold extra brush heads.

TROUBLESHOOTING

If you experience any issues with your Oral-B iO3 Electric Toothbrush, please consider the following common solutions:

Problem	Possible Cause	Solution
Toothbrush not turning on.	Low battery or not charged.	Place the toothbrush on the charger and ensure it is plugged in. Allow sufficient time for charging.
Brush head not vibrating effectively.	Brush head not properly attached or worn out.	Ensure the brush head is securely clicked into place. Replace the brush head if it is old or worn.
Pressure sensor light constantly red.	Applying too much pressure.	Reduce the pressure applied during brushing. The light should turn green when optimal pressure is achieved.

SPECIFICATIONS

- **Model Number:** iO G3.1A6.0
- **Power Source:** Battery Powered (1 Lithium Ion battery included)
- **Item Firmness Description:** Soft
- **Special Feature:** Pressure Sensor
- **Package Dimensions:** 10.16 x 7.09 x 4.09 inches
- **Item Weight:** 15.2 ounces
- **UPC:** 069055138052

Choose the right brush for you



	iO7	iO5	iO3
Gum Pressure Sensor	3 intensities	3 intensities	3 intensities
Cleaning modes	● ● ● ● ●	● ● ● ● ●	● ● ● ○ ○
Brushing feedback	Real-time (through A.I.)	Real-time (through A.I.)	
Display	✓ Black & white	Led-mode indicator	
3h Fast Charging	✓		
Travel Case	✓		

Figure 5: Feature comparison across different Oral-B iO series models.

WARRANTY AND SUPPORT

This Oral-B iO3 Electric Toothbrush is manufactured by Procter & Gamble. For specific warranty information, product registration, or technical support, please refer to the documentation included with your product or visit the official Oral-B website. Do not use on children under age 3.

Oral-B is the #1 brand used by dentists worldwide, based on regular surveys of a representative worldwide sample of dentists carried out for P&G.



Most used brand
by Dentists worldwide



Figure 6: Oral-B is recognized as the most used brand by dentists worldwide.