

[Manuals.plus](#) /

> [CINCOM](#) /

> CINCOM Hand Massager with Heat and Compression User Manual (Model CM-026H)

CINCOM CM-026H

CINCOM Hand Massager with Heat and Compression User Manual

Model: CM-026H

1. IMPORTANT SAFETY INFORMATION

Warning: Always keep the massager away from direct contact with skin to avoid burns. Do not use if you have a pacemaker or other electronic medical devices. Consult your doctor before use if you have any medical conditions. Keep the massager away from water and do not charge it while in use. Always unplug and turn off the massager when not in use.

2. PRODUCT OVERVIEW

The CINCOM Hand Massager is designed to provide comprehensive hand massage therapy through deep kneading and rhythmic air compression. It targets acupuncture points across fingers, palms, and wrists to support relaxation and circulation.

Smart Electric Hand Massager

Comprehensively care for your hands with multiple massage techniques



Palm Pain



Finger Pain



Wrist Pain



Hand Pain

Conceptual image illustrating the massage points and air compression zones within the hand massager.

- **All-Around Hand Massage Therapy:** Delivers deep kneading with 4 massage heads and rhythmic air compression.
- **Soothing Heat:** Features two heating levels (95–104°F) to help reduce hand fatigue, stiffness, and coldness.
- **Adjustable Modes:** Customize your session with 3 massage modes, 3 intensity levels, and 2 heat options.
- **Cordless and Rechargeable:** Built-in 2500mAh lithium battery for portable use.

3. PACKAGE CONTENTS

Verify that all items are present in the package:

What's in the Box?

1



3



4



2



5



The CINCOM Hand Massager package includes the massager unit, a power adapter, a USB Type-C charging cable, and the user manual.

1. CINCOM Hand Massager Unit
2. USB Type-C Charging Cable
3. Power Adapter
4. User Manual

4. SETUP

Before first use, fully charge the hand massager.

1. Connect the USB Type-C charging cable to the massager's charging port.
2. Plug the other end of the cable into the provided power adapter, and then into a wall outlet.
3. The battery indicator will show charging status. A full charge typically takes 1 hour 45 minutes.
4. The built-in 2500mAh lithium battery lasts 4–5 hours per full charge.

Battery Smart for Longer



The massager features a USB Type-C charging port and a long-lasting battery for cordless operation.

5. OPERATING INSTRUCTIONS

Follow these steps for an optimal massage experience:

Easy to Control



Intensity Button

3 intensity levels



Mode Button

3 mode programs



Power Button

Press it for 3 seconds



Heat Button

2 heating levels

The control panel features intuitive buttons for power, intensity, mode, and heat settings.

1. **Power On:** Press and hold the Power Button (U) for 3 seconds to turn on the device.
2. **Insert Hand:** Place your hand into the massager, ensuring your palm rests comfortably on the internal nodes and your fingers are aligned. The device is designed for both left and right hands.
3. **Select Mode:** Press the Mode Button (m) to cycle through 3 massage modes: Soothing, Vitality, and Relaxation.

4. **Adjust Intensity:** Press the Intensity Button (|||) to choose from 3 intensity levels: Low, Middle, and Highest.
5. **Choose Heat Level:** Press the Heat Button (Heat) to select between 2 heating levels (Low or High) or turn the heat off.
6. **Enjoy Massage:** The massager will operate for 15 minutes before automatically shutting off.
7. **Power Off:** Press and hold the Power Button (U) for 3 seconds to turn off the device manually.

Official CINCOM video demonstrating the hand massager's features and operation.

Demonstration of the Cincom Cordless Hand Massager, highlighting its use in a desk setting.

6. MAINTENANCE

- **Cleaning:** Wipe the exterior of the massager with a soft, damp cloth. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the massager in a cool, dry place away from direct sunlight when not in use.
- **Battery Care:** To prolong battery life, avoid fully discharging the device frequently. Charge it regularly, even if not used for extended periods.

7. TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on	Low battery or power button not pressed correctly.	Charge the device fully. Press and hold the power button for 3 seconds.
No heat or weak heat	Heat function not activated or set to low.	Press the Heat button to cycle through heat levels (Low, High, Off). Ensure the heat function is active.
Weak or inconsistent massage	Low intensity setting or improper hand placement.	Increase the intensity level. Adjust hand position to ensure proper contact with massage nodes.

8. SPECIFICATIONS

Feature	Detail
Product Dimensions	4.3 x 12.5 x 3.1 inches
Item Model Number	CM-026H-BU

Feature	Detail
Batteries	1 AAA batteries required (for remote, if applicable, or internal). Internal: 2500mAh Lithium.
Item Weight	2.7 Pounds
Power Source	Battery Powered
Special Features	Auto Shut-Off, Cordless, Heat, Lightweight, Quiet
Massage Therapy Technique	Compression
Heat Levels	2 (95–104°F)
Massage Modes	3
Intensity Levels	3
Auto Shut-Off Timer	15 minutes
Material Type	Plastic

9. WARRANTY AND SUPPORT

CINCOM offers reliable wellness tools and responsive customer support. For warranty information or technical assistance, please refer to the contact details provided in your product packaging or visit the official CINCOM website.