

## Fxtiaa LS-001

# Fxtiaa 7-in-1 LED Light Therapy Facial Wand User Manual

Models: LS-001 / LS-002

[Introduction](#)[Safety Information](#)[Product Components](#)[Setup](#)[Operation](#)[Maintenance](#)[Warranty &](#)[Troubleshooting](#)[Specifications](#)[Support](#)

## 1. INTRODUCTION

Thank you for choosing the Fxtiaa 7-in-1 LED Light Therapy Facial Wand. This device utilizes natural light waves transmitted by LEDs to deliver energy to the skin, promoting various skin benefits. Designed for home use, it offers a professional-grade skincare experience. Please read this manual thoroughly before use to ensure proper operation and safety.

## 2. SAFETY INFORMATION

- This device is intended for cosmetic use only. Any medical conditions should be discussed with a healthcare professional.
- Do not use on broken, irritated, or infected skin.
- Avoid direct eye contact with the LED lights.
- Keep out of reach of children.
- Do not immerse the device in water or use it in excessively wet environments.
- If you experience any discomfort or adverse reactions, discontinue use immediately and consult a physician.
- Use only the provided charging cable.

## 3. PRODUCT COMPONENTS

The Fxtiaa 7-in-1 LED Light Therapy Facial Wand consists of the following main parts:

- **Massage Head with LED Lights:** This rotating head emits the various colored LED lights and provides warming facial massage.
- **On/Off Button:** Powers the device on and off.
- **Light Color Adjust Button:** Cycles through the 7 different LED light colors.

- **Charge Indicator:** Displays charging status.
- **Type-C Charging Port:** For connecting the charging cable.



Figure 1: Diagram of the Fxtiaa 7-in-1 LED Light Therapy Facial Wand, showing the massage head, LED lights, control buttons, charge indicator, and Type-C charging port.

## 4. SETUP

### 4.1. Charging the Device

1. Connect the provided Type-C charging cable to the charging port on the device.
2. Plug the other end of the cable into a suitable USB power adapter (not included).
3. During charging, the green light indicator will flash.
4. Once fully charged, the red light indicator will remain steadily on.

Ensure the device is fully charged before its first use for optimal performance.

## 5. OPERATION

## 5.1. General Usage Guidelines

- Recommended usage: 3-5 times per week.
- Session duration: 5-10 minutes per session.
- The device features a warming facial massage function that works in conjunction with the light therapy. This helps dilate blood capillaries, allowing light therapy to penetrate deeper into the skin and enhance product absorption.
- The massage head can rotate 90 degrees for ergonomic use on different facial contours.

## 5.2. Step-by-Step Usage

1. **Prepare Skin:** Cleanse your face thoroughly and apply your preferred skincare products (e.g., serum, moisturizer).
2. **Power On:** Press and hold the On/Off button to turn on the device.
3. **Select Light Color:** Press the Light Color Adjust button to cycle through the available 7 light colors. Choose the color suitable for your specific skin concern.
4. **Apply to Skin:** Gently press the massage head against your skin. Move the device slowly across your face, neck, and other desired areas. For specific points, hold the device in place for approximately 10 seconds, repeating 2-3 times.
5. **Session Duration:** Continue for 5-10 minutes.
6. **Power Off:** Press and hold the On/Off button to turn off the device after your session.

## Meet different skin types

- DRY SKIN
- MIXE SKIN
- OILY SKIN
- NEUTRAL SKIN

## HOW TO USE ?

1



Clean face skin  
Apply skin care products

2



Turn on the facial wand  
choose light of appropriate color  
use for 10-15 minutes

3



Gently press the massage head  
Stay at each point for about 10 seconds  
and repeat for 2-3 times

4



According to the above method  
massage the cheeks  
chin and other parts will make your skin better

Figure 2: Illustrated guide for using the Fxtiaa 7-in-1 LED Light Therapy Facial Wand.

### 5.3. LED Light Functions

The device offers 7 distinct LED light colors, each targeting different skin concerns:

| Light Color  | Wavelength (nm) | Primary Function                                                                      |
|--------------|-----------------|---------------------------------------------------------------------------------------|
| Red Light    | 630nm           | Increases blood circulation, stimulates collagen production, promotes skin whitening. |
| Green Light  | 520nm           | Improves pigmentation, reduces wrinkles and fine lines, anti-aging benefits.          |
| Blue Light   | 470nm           | Calms and tightens skin, helps with acne.                                             |
| Yellow Light | 590nm           | Smooths skin, reduces redness.                                                        |
| Purple Light | 390nm           | Relaxing effect, improves lymph metabolism.                                           |

| Light Color   | Wavelength (nm) | Primary Function                                                   |
|---------------|-----------------|--------------------------------------------------------------------|
| Cyanine Light | 415nm-520nm     | Soothing, can help with allergies.                                 |
| White Light   | Mixed           | Accelerates tissue metabolism, tightens skin, improves fine lines. |



Figure 3: Overview of the 7 LED light functions and their benefits.

## 6. MAINTENANCE

- **Cleaning:** After each use, gently wipe the massage head with a soft, damp cloth. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the device in a cool, dry place, away from direct sunlight and extreme temperatures.
- **Battery Care:** To prolong battery life, avoid fully discharging the device frequently. Charge it when the low battery

indicator appears.

## 7. TROUBLESHOOTING

| Problem                          | Possible Cause              | Solution                                                                                                                                            |
|----------------------------------|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Device does not turn on.         | Low battery.                | Charge the device using the provided Type-C cable. The green light will flash during charging, and the red light will be steady when fully charged. |
| Device automatically shuts down. | No operation for 5 minutes. | This is a power-saving feature. The green light will flash for 5 seconds before shutdown. Simply turn the device back on to resume use.             |
| Device automatically shuts down. | Battery is too low.         | The red light will flash for 5 seconds before shutdown. Charge the device immediately.                                                              |

## 8. SPECIFICATIONS

|                    |                                                      |
|--------------------|------------------------------------------------------|
| Model Number       | LS-001, LS-002                                       |
| Product Dimensions | 8.11 x 3.38 x 1.29 inches                            |
| Item Weight        | 4.66 ounces                                          |
| Brand              | Fxtiaa                                               |
| Manufacturer       | Shenzhen Jingxin Electronics Co., Ltd                |
| Power Source       | Rechargeable (Type-C)                                |
| Light Colors       | 7 (Red, Green, Blue, Yellow, Purple, Cyanine, White) |

## 9. WARRANTY AND SUPPORT

For warranty information and customer support, please refer to the documentation included with your purchase or contact the manufacturer directly. Keep your proof of purchase for any warranty claims.

Manufacturer: Shenzhen Jingxin Electronics Co., Ltd

For further assistance, please visit the official Fxtiaa store on Amazon:[Fxtiaa Store](#)

### ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

**Email**

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

**Question**

Example: How do I reset this device or fix this error?

**Details**

Model number, symptoms, what you tried, or the section of the manual you are using.

**Submit question**