

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [Teglu](#) /
- › [Teglu Paraffin Wax Machine Model FHC-2222 User Manual](#)

Teglu FHC-2222

Teglu Paraffin Wax Machine Model FHC-2222 User Manual

Brand: **Teglu** | Model: **FHC-2222**

1. INTRODUCTION

This manual provides detailed instructions for the safe and effective use of your Teglu Paraffin Wax Machine Model FHC-2222. This device is designed for paraffin wax treatments for hands, feet, and elbows, offering a capacity of 4 liters and efficient heating capabilities. Please read this manual thoroughly before first use and retain it for future reference.



Image 1.1: The Teglü Paraffin Wax Machine FHC-2222 shown with included accessories such as paraffin wax blocks, thermal gloves, and a silicone brush.

2. SAFETY INSTRUCTIONS

WARNING: Failure to follow these safety instructions may result in electric shock, fire, or serious injury.

- Always unplug the machine from the power outlet when not in use, before cleaning, or before adding/removing parts.
- Do not operate the machine with a damaged cord or plug, or if the machine malfunctions or has been damaged in any manner.
- Keep the machine away from water or other liquids. Do not immerse the machine in water.
- Do not touch hot surfaces. Use handles or knobs. The wax and internal components will be hot during and after operation.
- Ensure the machine is placed on a stable, heat-resistant surface during operation.
- Test the wax temperature on a small area of your skin (e.g., wrist) before immersing hands or feet to prevent burns.

- This appliance is not intended for use by persons (including children) with reduced physical, sensory, or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Use only paraffin wax specifically designed for paraffin therapy.
- Do not add water or other substances to the wax.
- Do not use if you have open wounds, skin infections, or severe circulatory problems. Consult a physician if unsure.

3. PACKAGE CONTENTS

Verify that all items are present in the package:

- 1 x Teglu Paraffin Wax Machine (Model FHC-2222)
- 3 x 450g Paraffin Waxes (Lavender scent)
- 2 x Thermal Hand Gloves
- 100 x Plastic Gloves for Hands
- 100 x Plastic Gloves for Feet
- 1 x Silicone Brush
- 1 x Thermal Insulation Panel

4. SETUP

1. Unpack all components from the box.
2. Place the paraffin wax machine on a flat, stable, and heat-resistant surface. Ensure there is adequate ventilation around the unit.
3. Ensure the power cord is not damaged and is within reach of a suitable electrical outlet.

5. OPERATING INSTRUCTIONS

5.1. Melting the Paraffin Wax

1. Open the lid of the machine.
2. Place the desired amount of paraffin wax blocks into the aluminum liner. The machine has a 4L capacity, suitable for immersing hands, elbows, or feet (up to EU shoe size 47).
3. Close the lid securely.
4. Plug the machine into a power outlet. The display will illuminate.
5. Press the power button to turn on the machine.
6. Set the desired melting temperature and time using the control panel. The machine can adjust temperature from 30°C to 70°C (86°F to 158°F). For 1350g (3 blocks) of paraffin, the melting time is approximately 45 minutes due to the 200W fast heating element and aluminum liner.
7. The machine will begin heating. The display will show a countdown of the heating time.
8. Once the wax is melted, the intelligent automatic lid opening feature will activate, opening the lid to a 90-degree angle to accelerate cooling to a safe application temperature (typically around 53°C). A beep will sound

when the wax reaches the optimal temperature for use.

200W Fast & Gentle Heat

The ALUMINUM LINER reduces the melting time of paraffin to **30 mins(2 packs paraffin)-45mins(3 packs paraffin), 4 TIMES FASTER** than traditional plastic liner (2-3 hours)



ONE USE
0,1-0,15 KW/H



Image 5.1: Illustration of the 200W fast heating capability, showing 45 minutes for 3 packs of paraffin wax to melt.

Intelligent Automatic Lid Opening

Once the wax is melted,
the lid will be automatically opened to accelerate cooling



Wax melts



Automatic opening



Cool to fine temperature

Image 5.2: The intelligent automatic lid opening feature, which opens the lid to 90 degrees to facilitate cooling of the melted wax to the optimal temperature.



Image 5.3: A diagram illustrating the function description, including plugging in, pressing the switch, setting time, automatic lid opening, and cooling to 53°C.

5.2. Applying the Paraffin Wax

1. Once the wax has cooled to the appropriate temperature (around 53°C), place the thermal insulation panel into the melted wax. This panel helps prevent direct contact with the hot bottom of the machine.
2. Carefully dip your hand, foot, or elbow into the melted wax. Ensure the wax covers the desired area.

3. Remove the limb from the wax and allow the wax to solidify for a few seconds.
4. Repeat the dipping process 3 to 5 times to build up several layers of wax.
5. Immediately after the final dip, cover the waxed area with a plastic glove (provided).
6. For enhanced heat retention, place a thermal glove (provided) over the plastic glove.
7. Allow the wax to remain on the skin for approximately 15 minutes.
8. After 15 minutes, remove the thermal glove, then peel off the plastic glove and the solidified wax. The wax should come off easily.
9. Discard the used wax. Do not reuse wax that has been removed from the skin.

Smooth & Soften

Revitalizing rough hands and feet in minutes,
leaving skin hydrated and silky smooth

Hand



Sole



Feet



Elbow



Image 5.4: Step-by-step guide for using the paraffin wax machine, including placing wax, inserting the grid, dipping hands, and wearing gloves.

6. MAINTENANCE

6.1. Cleaning

1. Ensure the machine is unplugged and completely cooled before cleaning.
2. Once the remaining wax has solidified, gently remove it from the aluminum liner. You can often lift it out as a single block.
3. Wipe the interior and exterior surfaces of the machine with a soft, damp cloth. Do not use abrasive cleaners or immerse the unit in water.
4. Ensure all parts are dry before storing or next use.

6.2. Storage

Store the machine in a cool, dry place, away from direct sunlight and out of reach of children.

7. TROUBLESHOOTING

- **Machine does not turn on:** Check if the power cord is securely plugged into both the machine and the electrical outlet. Ensure the outlet is functional.
- **Wax not melting or melting slowly:** Ensure the lid is closed properly during heating. Verify that the temperature setting is appropriate for melting (e.g., 70°C).
- **Wax is too hot for application:** The automatic lid opening feature is designed to cool the wax to a safe temperature. If it still feels too hot, allow additional cooling time before use. Always test a small amount on your wrist first.
- **Wax is too cold for application:** Increase the temperature setting slightly and allow the machine to reheat.

8. SPECIFICATIONS

Feature	Specification
Product Dimensions	14.76 x 10.83 x 5.51 inches (37.5 x 27.5 x 14 cm)
Item Model Number	FHC-2222
Manufacturer	Teglu
Capacity	4000ml (4 Liters)
Voltage	US 110-120V; EU 220-240V (50/60Hz)
Power	200W
Temperature Range	30-70°C (86-158°F)
Temperature Control	NTC automatic temperature control

9. WARRANTY AND SUPPORT

No specific warranty information is provided in the product details. For warranty inquiries or technical support, please contact the manufacturer, Teglu, or your retailer directly. Retain your proof of purchase for any warranty claims.



