

## YA-MAN YJFM18N

# YA-MAN Photo Plus Shiny Neo Lift YJFM18N Beauty Device User Manual

Model: YJFM18N

## 1. INTRODUCTION

This manual provides essential information for the safe and effective use of your YA-MAN Photo Plus Shiny Neo Lift YJFM18N beauty device. Please read this manual thoroughly before first use and keep it for future reference.

The YA-MAN Photo Plus Shiny Neo Lift YJFM18N is a multi-functional RF beauty device designed to support skin care routines. It incorporates 3MHz RF (radio frequency) and EMS technology for enhanced daily skin care. The device also features a new "Triple RF" function, which utilizes YA-MAN's unique DYHP (Dynamic Hyper Pulse) technology for improved penetration and lifting effects. The LED brightness has been enhanced by 3 times compared to previous models.

## 2. SAFETY PRECAUTIONS

**WARNING: Failure to follow these safety instructions may result in injury or damage to the device.**

### 2.1 Do Not Use If You Are:

- Using medical electrical equipment (e.g., implantable pacemakers, life support systems like artificial heart-lung machines, or wearable electrocardiographs).
- Suffering from heart disease, disability, or suspected conditions.
- An infant or unable to express your will.
- Unable to perceive temperature or stimulation.
- Suffering from autoimmune diseases.

### 2.2 Consult a Doctor Before Use If You Are:

- Suffering from acute diseases, infectious diseases, malignant tumors, or febrile diseases.
- Suffering from high blood pressure or blood disorders.
- Experiencing fever.
- Suffering from sensory impairment due to severe peripheral circulatory dysfunction (e.g., diabetes).

- Suffering from neurological disorders, skin diseases, or atopic dermatitis.
- Undergoing medical treatment and taking medication.
- Having an allergic constitution, contact dermatitis, or other pre-existing skin conditions, or concerned about your skin during use.

### 2.3 Do Not Use On The Following Areas:

- Areas with skin injuries, rashes, inflammation, or other abnormalities.
- While intoxicated.
- Head, inside the mouth, private parts, or mucous membranes.
- Eyelids (orbital area) or eyeballs.
- Throat.
- Swollen moles or viral warts.
- Areas with bone fractures.
- Areas where metal has been implanted in the body due to surgical procedures (e.g., plastic, orthopedic, cosmetic surgery).
- Areas with internal bleeding, swelling, or scabs after cosmetic medical procedures.

*Note: Refrain from use for 2 weeks after such procedures. Consult a doctor before resuming use. (Examples: hyaluronic acid injection, silicone injection, Botox, thread lift, double eyelid surgery, HIFU, etc.)*

If you feel unwell or notice any skin abnormalities during use, immediately stop using the device and consult a doctor if the condition does not improve.

*Note: Some safety restrictions regarding users, usage areas, and usage conditions were relaxed from December 1, 2021, to allow more customers to use the product. There may be discrepancies between the instruction manuals included with some products and the data on the website. We appreciate your understanding.*

## 3. PACKAGE CONTENTS

---

Please check that all items are present in the package:

- Main Unit (YA-MAN Photo Plus Shiny Neo Lift YJFM18N)
- Charging Stand
- USB Charging Cable (USB Type-C)
- Cotton Stopper
- Instruction Manual (with Warranty)

## 4. PRODUCT OVERVIEW

---

The YA-MAN Photo Plus Shiny Neo Lift YJFM18N is a compact and powerful beauty device designed for comprehensive facial care. It integrates multiple technologies to support skin health and appearance.



Image: The YA-MAN Photo Plus Shiny Neo Lift YJFM18N beauty device shown with its charging stand and USB-C cable.  
The device has a sleek, ergonomic design with a treatment head and control buttons.

### Key Features:

- **RF (Radio Frequency) Technology:** Utilizes 3MHz RF for deep warming and skin conditioning.
- **EMS (Electrical Muscle Stimulation):** Delivers gentle electrical pulses to stimulate facial muscles.
- **Triple RF Function:** An advanced feature combining various RF frequencies for enhanced effects.
- **DYHP (Dynamic Hyper Pulse) Technology:** YA-MAN's proprietary technology for improved cosmetic ingredient penetration into the stratum corneum and lifting sensation.
- **Enhanced LED:** Features LED lights with 3 times the brightness for skin benefits.
- **Cordless Operation:** Designed for convenient use without being tethered to a power outlet.

## 5. SETUP

### 5.1 Charging the Device

1. Connect the USB Type-C charging cable to the charging stand.
2. Connect the other end of the USB cable to a compatible USB power adapter (not included) or a USB port.

3. Place the main unit onto the charging stand. Ensure it is seated correctly.
4. The charging indicator light on the device will illuminate to show it is charging.
5. A full charge takes approximately 3 hours. The indicator light will change or turn off once charging is complete.

### 5.2 Initial Patch Test

Before using the device on your face for the first time, always perform a patch test on an inconspicuous area of your skin (e.g., inner arm) to check for any adverse reactions. Follow the operating instructions for a short duration on the test area.

## 6. OPERATING INSTRUCTIONS

The device offers various modes for different skin care needs. Refer to the detailed instructions for each mode in the full manual. A typical usage sequence is as follows:

1. **Preparation:** Remove all makeup and cleanse your face thoroughly.
2. **Apply Cotton (if necessary):** If a mode requires it, attach a cotton pad to the electrode head using the cotton stopper. Apply an appropriate amount of lotion or serum to the cotton pad or directly to your skin.
3. **Power On:** Press and hold the power button to turn on the device.
4. **Select Mode and Level:** Press the mode button to cycle through available modes. Use the level adjustment buttons to select your desired intensity level. Start with the lowest level and gradually increase as comfortable.
5. **Treatment:** Gently glide the device over your skin according to the specific instructions for each mode. Avoid prolonged use on a single spot.
6. **Power Off:** After completing your treatment (approximately 30 minutes of operating time per charge), press and hold the power button to turn off the device.

## 7. MAINTENANCE

### 7.1 Cleaning the Device

- Always turn off the device and disconnect it from the charging stand before cleaning.
- Wipe the electrode head and main unit with a soft, dry cloth.
- For stubborn stains, lightly dampen the cloth with water or a mild, non-abrasive cleaner. Do not use alcohol, thinner, or corrosive detergents.
- Ensure no liquid enters the device's internal components.
- Allow the device to dry completely before storing or recharging.

### 7.2 Storage

- Store the device in a cool, dry place, away from direct sunlight and high humidity.
- Keep out of reach of children.

## 8. TROUBLESHOOTING

If you encounter any issues with your device, please refer to the following common problems and solutions:

Problem	Possible Cause	Solution
---------	----------------	----------

Problem	Possible Cause	Solution
Device does not turn on.	Battery is low or depleted.	Charge the device fully.
No sensation during use.	Intensity level is too low; insufficient contact with skin; skin is too dry.	Increase the intensity level; ensure full contact with the skin; apply more lotion or serum to keep skin moist.
Skin feels uncomfortable or irritated.	Intensity level is too high; prolonged use on one area; allergic reaction.	Reduce the intensity level; move the device continuously; stop use immediately and consult a doctor if irritation persists.
Charging indicator not working.	Cable or charging stand not properly connected; device not correctly placed on stand.	Check all connections; ensure the device is correctly seated on the charging stand.

If the problem persists after trying these solutions, please contact customer support.

## 9. SPECIFICATIONS

Feature	Specification
Model Number	YJFM18N
Brand	YA-MAN
Color	Gold
Dimensions (with cotton stopper)	Approx. W47 × D44 × H176 mm
Weight (with cotton stopper)	Approx. 210g
Materials	ABS resin, PET resin, Silicone rubber, Stainless steel
Power Source	Battery-powered (Lithium battery built-in)
Rated Input (Charging Stand)	DC5V 1.0A
Rated Output (Charging Stand)	DC9V 0.55A
Power Consumption (during charging)	Approx. 4.5W
Charging Time	Approx. 3 hours
Operating Time	Approx. 30 minutes
Country of Manufacture	Japan

## 10. WARRANTY AND SUPPORT

Your YA-MAN Photo Plus Shiny Neo Lift YJFM18N comes with a manufacturer's warranty. Please refer to

the warranty card included in your package for details on coverage and duration.

For technical support, troubleshooting assistance beyond this manual, or warranty claims, please contact YA-MAN customer service. Contact information can typically be found on the official YA-MAN website or on the warranty card.

Please retain your proof of purchase for warranty purposes.