

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

› [CYSSJF](#) /

› [CYSSJF Wireless Caregiver Call Button Watch SJ-300PLUS+SW 1+2 Instruction Manual](#)

CYSSJF SJ-300PLUS+SW

CYSSJF Wireless Caregiver Call Button Watch SJ-300PLUS+SW 1+2 Instruction Manual

Model: SJ-300PLUS+SW 1+2

1. INTRODUCTION

The CYSSJF Wireless Caregiver Call Button Watch is designed to provide reliable and convenient assistance for seniors, patients, and loved ones. This system includes emergency call buttons that can be worn on the wrist or as a pendant, paired with a portable pager for caregivers. It ensures timely alerts, offering peace of mind for both users and their caregivers.

SMART WIRELESS WATCH CAREGIVER PAGER



A good helper for caregiver.
Caregivers don't need to be around at all times.



Image 1: Overview of the Wireless Caregiver Pager System.

2. WHAT'S IN THE BOX

- 2 x Emergency Watch Call Buttons
- 1 x Portable Watch Pager (Receiver)
- 1 x Charger

3. PRODUCT FEATURES

- **Wireless Caregiver Pager with Watch:** The system includes two emergency watch call buttons and one portable watch pager.
- **Panic Button Watch:** Easy-to-use SOS watch call button for immediate assistance.
- **Watch Call Buzzer:** Features 4 ringtone prompt modes (shake & beep, beep, shake, mute/screen only) with a working range of up to 328 feet in open areas.
- **Fully Expandable System:** Multiple call bell watches can be added, and one emergency call button can be paired with multiple pagers.
- **Multi-scenario Application:** Ideal for patients, seniors, and loved ones in various settings, allowing caregivers freedom while ensuring timely alerts.

- **Waterproof Design:** Call buttons are waterproof for use in bathrooms or during activities where they might get wet.

PRODUCT SIZE OF RECEIVER

Small and portable



Image 2: Dimensions of the portable receiver unit.

4. SETUP AND INSTALLATION

1. **Charge the Portable Pager:** Ensure the portable watch pager is fully charged using the provided charger.
2. **Power On the Pager:** Turn on the portable watch pager.
3. **Pairing Call Buttons:** The call buttons are typically pre-paired with the pager. If re-pairing is needed, refer to the detailed instructions in the included manual.
4. **Placement of Call Buttons:** The emergency call buttons can be worn on the wrist using the included wrist strap or as a pendant with a lanyard. They can also be affixed to a surface using double-sided adhesive (not included in this specific package but often available for similar products).

DETAILS OF THE RECEIVER



Silicone material,
comfortable to wear



Sliding bolts,
easy installation and removal



Clear,
eye-protecting screen

Image 3: Close-up view of the receiver's features, including silicone material, sliding bolts, and eye-protecting screen.

5. OPERATING INSTRUCTIONS

1. **Initiating a Call:** To request assistance, simply press the SOS button on the watch call button.
2. **Receiving an Alert:** The portable watch pager will receive the signal and alert the caregiver according to its set mode (shake & beep, beep, shake, or mute/screen only). The screen will display the calling signal.
3. **Responding to a Call:** Caregivers should respond promptly upon receiving an alert to provide assistance.

CALL BUTTON FOR ELDERLY



Sturdy button

Get help with one click, easy to use and more comfortable.

Easy pairing

Follow the instructions to quickly pair and use.

Long standby time

Waterproof call button with a standby time of 6-12 months.



Image 4: Features of the call button, highlighting its sturdy design, ease of pairing, and extended standby time.

HOW TO USE



Press the watch button

↓
Calling

Press the watch button when the elderly/patient needs help.



Receive Signal
↓
Come for Help

Caregivers see the calling signal and provide timely assistance



Image 5: Visual guide on how to use the system, from pressing the watch button to receiving the signal on the pager.

6. MODES AND SOUNDS

The portable watch pager offers 4 distinct alert modes: shake & beep, beep, shake, and mute (screen only). It also features 58 different sounds, allowing for customization to suit various preferences and environments. The volume can be adjusted on the pager unit.

Official Product Videos

Video 1: Product Demonstration Video by CYSSJF. This video showcases the basic functions and appearance of the caregiver pager system.

Video 2: Product demonstration and operation video by CYSSJF. This video provides a detailed look at how to operate the device and its various features.

Video 3: Daytech Caregiver Pager Call Button. This video demonstrates the functionality of a similar caregiver pager system, including sound and light alerts.

Video 4: FullHouse Caregiver Pager Call Buttons for Elderly at Home. This video illustrates the use of call buttons in a home setting.

7. APPLICATION SCENARIOS

The CYSSJF Wireless Caregiver Pager is versatile and can be used in various environments where timely assistance is crucial:

- **Home:** For elderly individuals or those needing care at home, providing immediate contact with family members or caregivers.
- **Hospitals:** Can be used by patients to call nurses or staff for assistance.
- **Nursing Homes:** Facilitates communication between residents and care staff.
- **Rehabilitation Centers:** Supports patients in calling for help during recovery.

APPLICATION SCENARIO



Image 6: Examples of where the caregiver pager system can be effectively used, including hospitals, nursing homes, rehabilitation centers, and private residences.



Image 7: Illustration of the long-range transmission capability, showing effective coverage across multiple rooms and floors.

8. MAINTENANCE

- **Battery Replacement:** The call buttons require 1 CR123A battery (included). Replace batteries when the low battery indicator appears.
- **Cleaning:** Wipe the devices with a soft, dry cloth. Avoid using harsh chemicals or abrasive cleaners.
- **Waterproof Care:** While the call buttons are waterproof, avoid prolonged submersion. Ensure the battery compartment is securely closed after replacement to maintain water resistance.

9. TROUBLESHOOTING

- **No Alert from Pager:** Check if the pager is charged/powering on. Ensure call buttons are within range. Repair devices if necessary.
- **Weak Signal:** Ensure the pager is not obstructed by thick walls or large metal objects. Try repositioning the pager or call button.
- **Call Button Not Responding:** Replace the battery in the call button. Ensure it is not damaged.
- **Contact Support:** For persistent issues, refer to the full user manual or contact CYSSJF customer support.

10. SPECIFICATIONS

- **Voltage:** 12 Volts
- **Control Method:** Touch
- **Sensor Technology:** Wireless
- **Product Dimensions:** 1.97 x 2.01 x 0.63 inches
- **Item Weight:** 7 ounces
- **Item Model Number:** SJ-300PLUS+SW 1+2
- **Batteries:** 1 CR123A batteries required (included)
- **Power Source:** Battery Powered
- **Battery Description:** Lithium

11. WARRANTY AND SUPPORT

CYSSJF provides 18 months of free maintenance services for products and lifetime maintenance services. Within 18 months from the date of purchase, damages caused by poor material or design defects of the product can be

replaced or repaired free of charge. Repairs and replacements performed after 18 months will be charged the product cost. For any product issue, please contact the seller for warranty information.

Manufacturer: Quanzhou Shuojia Electronic Commerce Co., Ltd.