

## COMFIER CF-5003P

# COMFIER Shiatsu Feet Massager Machine

MODEL: CF-5003P

### Introduction

The COMFIER Shiatsu Feet Massager Machine, model CF-5003P, is designed to provide a comprehensive foot massage experience. It incorporates deep kneading, rolling, and compression massage techniques, complemented by a heating function, to help alleviate foot discomfort and promote relaxation. This manual provides essential information for the safe and effective use of your device.

Key features include innovative shiatsu rollers targeting acupuncture points, professional 360° compression massage with 3 adjustable levels for instep, heel, and ankles, and rapid heating up to 140°F. The device is easy to operate with four buttons and features a 15-minute auto-off setting for safety. Its portable design and removable, washable foot sleeves ensure convenience and hygiene.

### Setup

- 1. Unpack the Device:** Carefully remove the foot massager and all accessories from its packaging. Ensure all components, including the power adapter and remote control, are present.
- 2. Position the Massager:** Place the massager on a flat, stable surface. The bottom of the unit features non-slip material to prevent movement during use.
- 3. Connect Power:** Plug the power adapter into the designated port on the massager, then connect it to a suitable electrical outlet (110-120V for US models). The device is corded electric.
- 4. Prepare for Use:** Ensure the foot sleeves are properly inserted. For hygiene, it is recommended to wear socks during use.

# Deep Kneading & Rolling



Image: COMFIER Shiatsu Feet Massager Machine with remote and power cable. The image shows the foot massager unit, its remote control, and the power adapter, laid out on a surface.

# Professional 360° Full Wrap Massage



Image: Close-up of the COMFIER Shiatsu Feet Massager Machine's control panel, showing buttons for power, time, knead, compress, and heat.

Video: An overview of the Foot Massager with Heat, Kneading Rolling Compression functions. This video demonstrates the various features and how they operate.

## Operating Instructions

1. **Power On:** Press the Power button on the control panel or remote control to turn on the massager. The device will typically start in a default mode combining kneading, compression, and heat.

2. **Adjust Settings:** Use the dedicated buttons on the control panel or remote to customize your massage experience:

- **Knead:** Cycle through Low, Medium, High levels of kneading intensity.
- **Compression:** Cycle through Low, Medium, High levels of air compression intensity.
- **Heat:** Cycle through Low, High, or Off settings for the heating function. The heating function may take some time to reach its maximum temperature of 140°F.
- **Timer:** Select 15, 20, or 30-minute auto-off settings. The device will automatically shut off after the



selected time.

3. **Enjoy Massage:** Place your feet into the foot wells and adjust settings to your preference. The deep kneading and rolling functions target acupuncture points, while compression airbags provide a soothing squeeze to the instep, heel, and ankles.

**Warm Enough Max.140°F**

Fast Heating Up 1 Min

enhance blood circulation and improve sleep quality



Image: A person comfortably seated, using the COMFIER Shiatsu Feet Massager Machine. This illustrates typical usage in a home or office setting.



## Easy to Operate

Friendly for the Elderly

Image: A visual representation of the heat function within the COMFIER Shiatsu Feet Massager Machine, showing heat radiating around the feet.

Video: A user demonstrating the operation of the foot massager, highlighting the various settings and their effects.

## Maintenance

- 1. Cleaning:** Before cleaning, ensure the power supply is turned off and the plug has been pulled out from the socket. Use a soft, dry cloth to wipe the surface of the product and control panel. Do not use corrosive cleaning agents like benzene or thinner.
- 2. Washable Foot Sleeves:** The internal foot sleeves are removable and washable. Detach the sleeves via zipper and wash them by hand or machine. Allow to air dry completely before reattaching.
- 3. Storage:** Store the massager in a cool, dry place out of reach of children. Avoid contact with sharp edges or pointed objects which might cut or puncture the cover of the massager.





Image: A user demonstrating the removal of the washable foot sleeves from the COMFIER Shiatsu Feet Massager Machine for cleaning.

## Troubleshooting

Fault	Solutions
There is a clattering noise in the movement.	This is normal, caused by the working of the motor while in use.
The massager switches off suddenly.	The device will switch off when the timer runs out. Press the 'Power' button to continue using it after resting the unit for 15 minutes. Check if the device has been used for more than 30 minutes. If so, let the device rest for 15 minutes before continuing. The device may require longer time to cool in warmer environments.

Fault	Solutions
The massage device fails to start.	Check if adapter is plugged into socket. Check if adapter output is connected to device power inlet.

## Specifications

- **Product Dimensions:** 13.39 x 14.76 x 7.28 inches
- **Item Weight:** 7.32 Pounds
- **Item Model Number:** CF-5003P-US1
- **Manufacturer:** Comfier
- **Power Source:** Corded Electric
- **Material:** Plastic
- **Color:** Black
- **Max. Temperature:** 140°F
- **Foot Size Compatibility:** Fits up to Men's Size 12

## Safety Information

General safety precautions include using the massager on dry feet, avoiding sensitive areas or open wounds, checking for loose cords, and not using the massager near water. The product features an auto-off setting after 15 minutes for added safety. Always use on dry feet and avoid sensitive areas or open wounds. Store in a dry place. Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

## Warranty and Support

The foot massager comes with a 1-year warranty, which is extendable to 3 years. For any questions or support, the customer service team is available and will respond within 24 hours.