

Manuals.plus /

- › Oral-B /
- › Oral-B iO Series 3 Electric Toothbrush User Manual

Oral-B 8006540731536

# Oral-B iO Series 3 Electric Toothbrush User Manual

Brand: Oral-B | Model: 8006540731536

## 1. PRODUCT OVERVIEW

The Oral-B iO Series 3 electric toothbrush offers an advanced cleaning experience with its unique magnetic drive system. It delivers energy precisely to the bristle tips, which glide from tooth to tooth with gentle micro-vibrations for effective and comfortable brushing. This model features intelligent pressure control, multiple brushing modes, and a visual timer to ensure optimal oral hygiene.



Image: Oral-B iO Series 3 Electric Toothbrush, highlighting its packaging and 'Test Winner' status.

## 2. WHAT'S IN THE BOX

Your Oral-B iO Series 3 package typically includes:

- Oral-B iO Series 3 Electric Toothbrush Handle
- 1 x Oral-B iO Ultimate Clean Brush Head
- Charging Station
- Travel Case



Image: Contents of the Oral-B iO Series 3 package, including the toothbrush, charger, and travel case.

### 3. SETUP

1. **Unpack the Toothbrush:** Carefully remove all components from the packaging.
2. **Charge the Battery:** Place the toothbrush handle on the charging station. Ensure the station is connected to a power outlet. The battery indicator light will show charging progress. A full charge is recommended before first use.



Image: Diagram illustrating key components including the charging station.

## 4. OPERATING INSTRUCTIONS

3. **Apply Toothpaste:** Moisten the brush head and apply a small amount of toothpaste.
4. **Power On:** Press the power button to turn on the toothbrush.
5. **Use Pressure Control:** The intelligent pressure control warns with a red light if you apply too much pressure and illuminates green when optimal pressure is used. This protects your gums and ensures effective cleaning.
6. **Select Brushing Program:** Choose from 3 brushing programs for different oral care needs: Daily Clean, Sensitive, and Whitening. Press the mode button to cycle through the programs.
7. **Brush Your Teeth:** Guide the brush head slowly from tooth to tooth. Brush for the recommended 2 minutes. The visual timer will indicate when the recommended brushing time is reached.



Image: Intelligent Pressure Control indicating optimal pressure with a green light.

## BUNTER LEUCHTRING

Der 360° Ring blinkt in bunten Farben, wenn die von Zahnärzten empfohlene Putzzeit von 2 Minuten erreicht wurde.



Image: The 360° light ring indicating the 2-minute brushing timer.

## 5. KEY FEATURES

---

- **Advanced Cleaning Technology:** Removes 100% more plaque for healthier gums in just one week compared to a manual toothbrush.
- **Gum Protection:** The iO pressure control indicates optimal pressure with a green light, ensuring gentle yet effective cleaning.
- **Personalized Brushing:** Features 3 brushing programs (Daily Clean, Sensitive, Whitening) to suit individual needs.
- **Visual Timer:** The toothbrush lights up when the dentist-recommended 2-minute brushing time is achieved, and also reminds you when to change the brush head.
- **Round Brush Heads:** Oral-B's round brush heads reach areas that elongated manual toothbrushes cannot.
- **Long-lasting Battery:** Equipped with a reliable, long-lasting lithium-ion battery for extended use.

# Entfernt schonend 100% mehr Plaque\* in 1 Woche

mit der fortschrittlichsten Technologie von Oral-B



Image: Illustration of the brush head effectively removing plaque.

## 6. MAINTENANCE

- **Brush Head Replacement:** Replace the brush head as needed, typically every three months, or when the indicator bristles fade. The toothbrush will provide a reminder for brush head change.
- **Cleaning the Handle:** After each use, rinse the brush head thoroughly under running water. Wipe the handle with a damp cloth to remove any toothpaste residue.
- **Cleaning the Charging Station:** Periodically wipe the charging station with a damp cloth. Ensure it is unplugged before cleaning.

## 7. TROUBLESHOOTING

If you encounter any issues with your Oral-B iO Series 3 toothbrush, please refer to the following common solutions:

- **Toothbrush does not turn on:** Ensure the battery is charged. Place the handle on the charging station for at least 30 minutes.
- **Brush head does not vibrate/rotate:** Ensure the brush head is securely attached to the handle. If the issue persists, try a different brush head.

- **Pressure control light is always red:** You might be applying too much pressure. Try to brush more gently. The light should turn green when optimal pressure is applied.
- **Short battery life:** Ensure the toothbrush is fully charged. Battery performance can decrease over time with regular use.

For further assistance, please contact Oral-B customer support.

## 8. WARNINGS AND SAFETY INFORMATION

---

- Use with gentle pressure to avoid injury to gums and teeth.
- Do not use for purposes other than oral hygiene.
- Keep out of reach of children.
- Do not immerse the charging station in water.
- If the product is dropped, the brush head should be replaced before the next use even if no visible damage is apparent.

## 9. SPECIFICATIONS

---

Feature	Detail
Brand	Oral-B
Model Number	8006540731536
Power Source	Battery Powered (Lithium-Ion)
Number of Brushing Modes	3 (Daily Clean, Sensitive, Whitening)
Special Features	Brushing Timer, Intelligent Pressure Control, Micro-vibrations
Brush Head Shape	Round
Country of Origin	Germany
Dimensions (Packaging)	24.4 x 17 x 9.4 cm
Weight (Packaging)	580 Grams

## 10. WARRANTY AND SUPPORT

---

This product is new from an authorized dealer. Invoice including VAT is available. Shipping to packing stations is possible.

For support or warranty claims, please contact the manufacturer directly:

Sulzbacher Str. 40,  
65823 Schwalbach a. Ts.,  
Germany

For more information, please visit the official Oral-B website or refer to the legal notices provided with your purchase.

