

Medcursor MD-KP02

Medcursor MD-KP02 Heated Knee Massager User Manual

Model: MD-KP02 | Brand: Medcursor

1. INTRODUCTION

The Medcursor MD-KP02 Heated Knee Massager is a versatile, cordless device designed to provide targeted heat and vibration massage. This 3-in-1 unit can be applied to the knee, elbow, and shoulder areas to enhance comfort and support. This manual provides essential information for safe and effective use of your device.

2. IMPORTANT SAFETY INFORMATION

Please read all safety warnings and instructions carefully before using this product. Failure to follow these instructions may result in electric shock, fire, or serious injury.

- Ensure to use the massager for a maximum of 30 minutes to prevent overheating or excessive massage.
- Always keep the device away from water to avoid electrical hazards.
- Check the massager for any damage or defects before each use.
- Do not use the massager if you have any skin injuries, wounds, or infections in the targeted area.
- People with certain medical conditions, such as deep vein thrombosis, should consult their healthcare provider before using the massager.
- Always follow the manufacturer's instructions for proper use, handling, and maintenance.

3. PACKAGE CONTENTS

Verify that all items are present in the package:

- 1 x Medcursor Heated Knee Massager Unit
- 1 x 23-inch Extension Strap
- 1 x USB-A to USB-C Charging Cable (Note: Charging head not included)
- 1 x User Manual

Ideal Gift For Family and Friends



PACKAGE CONTENTS

- 1 x Knee Massager with Heat
- 1 x USB-A to USB-C Cable
- 1 x Extension Strap
- 1 x User Manual



Image: The Medcursor Heated Knee Massager package contents, including the massager unit, USB-C cable, extension strap, and user manual.

4. PRODUCT OVERVIEW

The Medcursor MD-KP02 features an ergonomic design with adjustable Velcro straps for a secure and comfortable fit. It integrates heating elements and vibration motors for comprehensive relief.

Key Features:

- **5 Heating Levels:** Adjustable temperatures ranging from 40°C (104°F) to 60°C (140°F).
- **3 Vibration Modes:** Features 2 vibration motors providing deep massage.
- **Cordless Design:** Built-in large-capacity rechargeable battery for portability.
- **3-in-1 Application:** Suitable for knees, elbows, and shoulders.
- **Easy Operation:** 4 intuitive buttons for power, temperature, vibration, and timer settings.
- **30-Minute Timer:** Automatic shut-off function for safety.

5. SETUP

5.1 Charging the Device

Before first use, fully charge the device. Connect the provided USB-A to USB-C charging cable to the charging port on the control panel and a compatible USB power adapter (not included).

The display will indicate charging status. A full charge ensures optimal performance and portability.

5.2 Attaching the Massager

The massager is designed for flexible application. Use the adjustable Velcro straps to secure it firmly but comfortably to the desired body part.

For larger areas or shoulder application, utilize the 23-inch extension strap.



Image: Medcursor Heated Knee Massager with charging cable and extension strap. The main unit is black with a digital display, accompanied by a USB-C cable and a rolled-up extension strap.

6. OPERATING INSTRUCTIONS

6.1 Basic Operation

1. **Power ON/OFF:** Long press the power button (usually the central button) for 2 seconds to turn the

device ON or OFF.

2. **Adjust Temperature:** Press the flame icon button to cycle through the 5 heating levels (40°C to 60°C).
3. **Select Vibration Mode:** Press the wave icon button to cycle through the 3 vibration modes.
4. **Set Timer:** Press the clock icon button to set the automatic shut-off timer (15, 30, or 45 minutes). The device will automatically turn off after the set time.
5. **Lock/Unlock Controls:** Short press the power button to lock or unlock the control panel. This prevents accidental changes to settings during use.

Video: A demonstration of the Medcursor Heated Knee Massager, highlighting its cordless design and ease of use for various activities.

6.2 Application Areas

The versatile design allows for use on multiple body parts:



Image: The Medcursor massager shown applied to a knee, a shoulder, and an elbow, demonstrating its 3-in-1 versatility.

- **Knee:** Position the massager over the knee joint and secure with straps.
- **Elbow:** Wrap the massager around the elbow, using the extension strap if needed for a snug fit.

- **Shoulder:** Use the extension strap to secure the massager comfortably around the shoulder area.

Video: A demonstration of the HOMRELEXA knee heating pad for knee pain, showing how to apply and operate the device on the knee. The video highlights the ease of use and the heating function.

7. MAINTENANCE

7.1 Cleaning

Wipe the massager's surface with a soft, damp cloth. Do not immerse the device in water or use harsh cleaning agents.

7.2 Storage

Store the device in a cool, dry place away from direct sunlight and extreme temperatures. Keep out of reach of children.

8. TROUBLESHOOTING

If you encounter any issues with your Medcursor MD-KP02, please refer to the following common problems and solutions:

Problem	Possible Cause	Solution
Device does not turn on	Low battery; Device is locked	Charge the device fully; Short press the power button to unlock.
No heat/vibration	Heat/vibration mode not activated; Device malfunction	Ensure respective mode buttons are pressed; Contact customer support if issue persists.
Device turns off unexpectedly	Timer expired; Low battery	Check timer settings; Recharge the device.







9. SPECIFICATIONS

- **Brand:** Medcursor
- **Model Number:** MD-KP02
- **Product Dimensions:** 28 x 14.57 x 1 inches
- **Item Weight:** 1 Pounds
- **Material:** Fabric
- **Power Source:** Battery Powered (1 Lithium Ion battery included)
- **Heating Levels:** 5 (40°C/104°F to 60°C/140°F)
- **Vibration Modes:** 3
- **Timer:** 15, 30, 45 minutes auto-off

10. WARRANTY AND SUPPORT

MEDCURSOR provides customer support for its products. For any inquiries, technical assistance, or warranty claims, please refer to the contact information provided in your product packaging or visit the

Related Documents - MD-KP02

	<p>Medcursor MD-KP02 Knee Massager with Heat User Manual - Operation, Care, Troubleshooting</p> <p>Comprehensive user manual for the Medcursor MD-KP02 Knee Massager with Heat. Learn about operation, safety precautions, care, maintenance, troubleshooting, and warranty information for this cordless massager.</p>
	<p>Medcursor MD-HM02 Hand Massager User Manual</p> <p>Comprehensive user manual for the Medcursor MD-HM02 Hand Massager, covering important safeguards, operation instructions, care, maintenance, troubleshooting, and warranty information.</p>
	<p>Medcursor MD-LM01 Air Compression Leg Massager with Heat - User Manual</p> <p>Comprehensive user manual for the Medcursor MD-LM01 Air Compression Leg Massager with Heat. Learn about product features, safety instructions, operation, troubleshooting, and care.</p>
	<p>Medcursor MD-MG06 Deep Tissue Massage Gun User Manual</p> <p>User manual for the Medcursor MD-MG06 Deep Tissue Massage Gun, providing instructions on safe operation, features, maintenance, troubleshooting, and warranty information.</p>
	<p>Medcursor Kneading Massager User Manual and Safety Guide</p> <p>Comprehensive guide for the Medcursor Kneading Massager, covering safe operation, features, specifications, maintenance, and warranty. Learn how to get deep kneading relief for hard-to-reach muscles.</p>
	<p>Medcursor Foot Massager User Manual - Model MD-52230</p> <p>Comprehensive user manual for the Medcursor Foot Massager (Model MD-52230). Learn how to operate, maintain, and troubleshoot your device. Includes safety warnings and warranty information.</p>